

OPER 0.4 MRI system is a permanent magnet MRI. The 0.4 Tesla has all functions of the middle field MRI and low operating costs.

www.shimadzu.eu medical@shimadzu.eu

Shimadzu Europa GmbH Albert-Hahn-Str. 6-10 47269 Duisburg, Germany

Phone +49 - 203 - 7687-0 Fax +49 - 203 - 7687680

Brilliant

The unique and brand-new BRANSIST salire digital x-ray angiography system family provides outstanding high-quality clinical images to support safer and quicker intervention procedures. This is achieved particularly with the latest in flat panel detector (FPD) technology called "Direct-conversion."

Through flexible system configurations, BRANSIST safire meets perfectly the needs of today's clinical demands: it is available with a 9 inch or 17 inch FPD, as ceiling or floor mounted version, or even as a bi-plane system.

Functions and features such as 3D-imaging, No Mask DSA, SYNCHRONAVI, SUREengine Image Enhancement and CT-like imaging support the full range of intervention applications. BRANSIST safire provides flexible system operation, fast digital image processing and patient comfort:

- Outstanding image quality achieved particularly by "Direct-conversion" FPD
- High-speed response and various applications ensured by latest digital technology
- Fast, easy and precise C-Arm positioning up to 60 degrees per second
- Advanced dose management ensuring safety of patients and operators.




BRANSIST safire – outstanding image quality

www.shimadzu.eu









Dear Reader,

There is no modality in radiology that in the past few years has seen such stunning quantum leaps in technological development as computed tomography. Studies, however, show that with each additional row radiation exposure for the patient increases.

With the introduction of the new premium models that offer innovative scanning technologies doses - particularly in cardiac examinations - started to decrease. Two to three years ago the standard radiation dose of a heart scan was approximately 10 mSv. Today, it is 1 to 2 mSv. An enormous success!

Maybe the volume scanners will end the slice race once and for all and research will focus again on patient benefit. While in projection radiography there are clear guidelines, computed tomography does not apply a quantifiable image quality standard which is directly related to patient dose.

As far as radiation exposure is concerned, what might in fact be more important than official guidelines is professional handling of the systems. Modern CTs combined with sophisticated post-process-

> ing workstations offer a host of functionalities that ensure high quality images at low doses. The different technological approaches the manufacturers chose to master the dose challenge are described in the chapter on CT in this issue.

> The RadBook – as usual – presents an independent overview of all systems and accessories, accompanied by many insightful articles. And best of all: as of now RadBook is available online.

> > **Imprint**

See you at www.radbook.eu

Daniela Zimmermann Guido Gebhardt

Published by:

EUROPEAN HOSPITAL Verlags GmbH Theodor-Althoff-Straße 39, 45133 Essen, Germany, Phone +49-201-87126-851, www.european-hospital.com

and XT Media Guido Gebhardt, Adalbert-Stifter-Weg 2 b, 85661 Forstinning, Germany, Phone +49-8121-907630, www.radiologieforum.de

Editor-in-Chief: Guido Gebhardt, Daniela Zimmermann, Executive Director: Daniela Zimmermann, **Advertising**: Denise Hennig (D, A, CH), Gavin Hua (China), Simon Kramer (Europe), Hanna Politis (U.S.A, Canada), Hannah Wizer (Israel), **Subscription:** Janka Hoppe, Subscription rate: € 19.- plus postage

© 2010 by EUROPEAN HOSPITAL Verlags GmbH and XT Media, Guido Gebhardt. All rights reserved



Editorial	3
Computed tomography	5
ASIR – FORGET COMPROMISE: HIGH IMAGE QUALITY AT LOW DOSE	6
PHILIPS DESIGNS FOR MAXIMUM DIAGNOSTIC RESULTS WHILE LOWERING RISKS	8
DIGITAL VOLUME TOMOGRAPHY	IU
MISSION POSSIBLE: REDUCING RADIATION DOSE IN CT.	٦١
Low Dose and Robustness of Aquilion ONE for Cardiac CT	Iŏ
MR systems	33
THE POLESTAR SURGICAL MRI SYSTEM	34
Injectors	46
Interventional systems	52
Mammography	63
TEACHING HOSPITAL EVALUATES 3D IMAGING FOR BREAST CANCER SCREENING	
Breast Tomosynthesis	55
PIONEERS BY TRADITION.	1211
Mammography at its Best	68
R/F systems	79
SHIMADZU – PIONEERS OF VISION	11111
Computed & digital radiology	92
WORLD'S FIRST PURPOSE-BUILT DIRECT DIGITAL RSA X-RAY SYSTEM	93
Molecular imaging	113
Displays /printers	120
Ultrasound systems	130
COMPANIES & SUPPLIERS	ין / ר





is a powerful tool for making sound diagnoses and following disease evolution. Yet for years, its imaging power has been constrained by the need to limit the radiation dose delivered to patients.

Now, at last, clinicians can have it both ways. GE Healthcare introduces Adaptive Statistical Iterative Reconstruction (ASiR), an image reconstruction technique that ASiR enables dose reduction up to 50% for DiscoveryTM CT750 and up to 40% for LightSpeedTM VCT and BrightSpeedTM Elite across the whole body.

ASiR combined with SnapShot Pulse allows dose reduction up to 83% for cardiac exams and boosts image quality by improving Low Contrast Detectabilty.

ASiR aligns with the clinical consensus around ALARA (As Low As Reasonably Achievable) practice that aims to limit patient's X-ray exposure as much as possible.

ASiR opens new perspectives in CT imaging. Unlike conventional dose reduction techniques like collimators and mA modulation, its benefits are universal, independent of patient size or anatomic region.

ASiR™ - Compact

ASiR is a new CT reconstruction engine that helps reducing dose while improving image quality. ASIR models the noise root causes and subtracts it throughout multiple iterations. ASiR is not a simple filter but a true raw data iterative reconstruction.


- Integrated Vectorial Reconstruction Board
- Up to 35 fps reconstruction
- Raw data space reconstruction
- Patient specific dose optimization
- Improved Low Contrast Detectability
- Available on DiscoveryTM CT750 HD, LightSpeedTM VCT Series, BrightSpeedTM Elite
- Lowers dose across all anatomy: Neuro, Cardiac & Colonography, Chest, Abdo/Pelvis

Experience in routine

- More than 400.000 patients examined to date
- Dose and image quality results confirmed in more than 20 international clinical studies

ASiR technology is especially beneficial in higher risk populations, such as children and young women, in whom the greater sensitivity of growing tissues makes strict dose limitation a real necessity. Based on wide experience in patients scanned with ASiR, GE customers can attest to consistently high image quality with low dose levels.

Conventional CT image reconstruction techniques are simple and fast, but have limitations, as they are sensitive to noise and artifacts.



ASiR extracts noise by modelling its root causes for each patient and application type.


Runoff vascular exam

3.5 mSv* effective dose (DLP: 300 mGy.cm)
ASiR allows lower kV on many applications such as vascular examinations, resulting in a significant dose reduction.
Courtesy: Dr JL. Sablayrolles, Centre Cardiologique du Nord, Saint Denis, France.
* Obtained by EUR-16262 EN, using an adult pelvis factor of 0.019*DLP and leg factor of 0.008*DLP.



Thoracic aorta angio and heart exam

1.08 mSv* effective dose (DLP: 27.6 mGy.cm)
7-day-old girl.
ASiR allows up to 60% dose reduction in these congenital heart disease exams.
Courtesy: Dr F. Stålhammar
Queen Silvia Children's Hospital,
Gothenburg, Sweden
* Dose conversion obtained
from ICRP 0-year-old Chest factor of 0.039*DLP.



Heart exam


0.78 mSv* effective dose
(DLP: 46 mGy.cm)
Prospective Gating mode, 100 kV, 87,5
mAs, ASiR Reconstruction.
Courtesy: JK. Min, MD, Cornell
University Medical Center, New York
Presbyterian Hospital
* Obtained by EUR-16262 EN, using

an adult chest factor of 0.017*DLP.




CT colonography exam


1 mSv* efffective dose (DLP: 52 mGy.cm, total dose for prone and supine series) Courtesy: Dr V. Barrau, Centre Cardiologique du Nord, Saint Denis, France.
* Obtained by EUR-16262 EN, using an adult pelvis factor of 0.019*DLP.



Thoracic Aorta and heart exam on obese patient (BMI=34)

2.2 mSv* effective dose (DLP=133 mGy.cm)
Whole thoracic aorta acquired in prospective gating mode, 100 kV, 175 mAs, ASiR reconstruction.
Courtesy: Dr. Pontone,
Centro Cardiologico Milano, Italy
* Obtained by EUR-16262 EN, using an adult chest factor of 0.017*DLP





RAD-BOOK 2010

Philips Designs for Maximum Diagnostic Results While Lowering Risks

hile advances in medical imaging are continuously opening new doors in both early diagnosis and disease management, regulators, clinicians and patients all have concerns about the potential risks of radiation dose. In the US, the Food and Drug Administration has just announced initiatives to address reports of acute overdoses in some American hospitals and to reduce lifetime exposure to radiation. Recent studies have shown that the average American's total radiation exposure has nearly doubled in the last 30 years and that medical radiation now accounts for more than half the American population's total radiation exposure whereas it used to account for just one-sixth.

These concerns about radiation dose risk and reduction make it increasingly important for radiologists to have access to imaging systems which limit exposure while maintaining image quality. Researchers at Philips have spent decades designing innovative products which produce maximum diagnostic results while lowering risks. This research has resulted in Philips' DoseWise radiation management which is a set of techniques, programmes and practices based on the ALARA (As Low As Reasonably Achievable) principle.

Philips' Most Advanced Solution: iDose for CT

The iDose iterative reconstruction technique is the latest addition to the DoseWise portfolio of radiation management tools. iDose uses a number of Philips technologies to lower X-ray dose by up to 80% while maintaining high image quality and fast reconstruction times.

iDose works in concert with other DoseWise tools including filters to optimize beam distribution by individual patient; detectors to provide ultra-high signal-to-noise ratio; collimators to eliminate unnecessary dose at the beginning and end of the examination and to reduce scatter and Automatic Current Selection (ACS), which adapts tube current to patient size. Dose modulation tools automatically adjust the dose delivered to the patients, compensating for their individual physiology and optimising the dose by anatomic region.


Together these tools enable the iDose to provide an equivalent diagnostic image quality at a fraction of the dose.

Original FBP reconstruction of an abdomen



iDose reconstruction of an abdomen









iDose reconstruction of a pelvis (VC)

In a blinded study, 20 radiologists at 14 global sites compared datasets from both full-dose scans and low-dose scans using iDose technology. The study showed that there was no significant difference found in the diagnostic image quality between the conventional full-dose images and the reduced dose images.

As Dr Rajiv Shah, Associate Professor of Radiology at the CWRU Metro Health Medical Center in Cleveland, Ohio commented: "It's a no-brainer. The processed images at the decreased dose are as good as the full-dose images." While other users, such as Dr Shawn Teague from the Indiana University School of Medicine, Indianapolis and Dr Shigeru Suzuki from the Teikyo University School of Medicine, Tokyo felt that the low dose images were better than full-dose images in some cases. Dr Suzuki, in particular, preferred the sharpness of lung parenchyma on the iDose image.

iDose in Clinical Routine

In order to streamline workflow for radiologists, iDose was designed to provide the same look and feel as conventional, higher-dose radiation images, therefore no additional training is needed. Radiologists simply need to plan their work as normal and then enable the iDose, using its dedicated low-dose protocols. The dose-lowering features of the iDose come into play automatically during the scan.

The iDose reconstruction technique is made possible by the RapidView IR console that contains hardware advances designed specifically to provide the performance requirements and processing power needed for low-dose scanning. It automatically generates statistical noise and anatomic structure models to improve operational efficiency. A comparative dual model (noise and anatomic) technique removes noise while preserving morphologic information.

The RapidView IR balances image quality, dose reduction, full-dose appearance and workflow. And it is fast – very fast. In fact, the RapidView IR provides radiologists with reconstruction times that are 20 times faster than those provided by conventional hardware.

The iDose is available both for clinics acquiring new Brilliance iCT and 64-channel scanners and for clinics with these already installed. Installation simply involves bringing in a new computer central processing unit and software, and deliveries can be expected in the second half of 2010.

DoseAware for Interventional Radiologists

The radiation exposure risk of hospital staff is a very real concern and DoseAware, another new addition to the DoseWise group of solutions, has been designed to help interventional radiologists track their X-ray exposure in real

time. This personal dose measuring system, which is the first of its kind, allows healthcare professionals performing potentially lengthy interventional procedures to manage their radiation exposure on the spot.

Using the DoseAware technology, all staff working in an X-ray environment wear a personal dose meter which measures scatter radiation and transmits this information to a base station mounted in the examination room, where it is displayed on an LCD screen. The staff can see from the screen if their X-ray exposure is in the green, yellow or red zone and can take immediate action to limit their exposure. The data can also be viewed on a computer screen and it can be retrieved from the base station for archiving, analysing and reporting

DoseAware does not replace the thermoluminescent dose meter (TLD) as the legal dose meter which reports accumulated exposure for a period of time but it does manage radiation exposure in realtime, putting staff in control of the way they work and allowing them to immediately change their behaviour and work patterns.

Research and innovation from Philips demonstrate the company's unique dedication to supporting its customers with systems that help reduce the radiation dose to patients and hospital staff, while maintaining exceptional image quality to aid in the diagnosis and treatment of patients.



COMPUTED TOMOGRAPHY

RAD-BOOK 2010

Digital Volume Tomography

tate-of-the-art radiological diagnostics require state-of-the-art equipment. Cone beam technology well-nigh revolutionized dental and medical x-ray diagnostics. Meeting the continuously increasing demand by orthodontists and ENT specialists for high quality and high resolution x-ray images, digital volume tomography (DVT) has proven to be a valuable modality.

DVT enjoys an excellent reputation with orthodontists and clinicians alike, in particular with regard to image quality in 3D skull base diagnostics. Superior to conventional computed tomography (CT), it offers significantly better resolution and requires less dosage. Its major advantage, however, is its impressive efficiency and safety record.

"This innovative technology provides radiologists with an opportunity to expand the service portfolio for referring physicians. The economies of scale that can be achieved with DVT make sure that the owners' investment enjoys long-term protection," explains Markus Hoppe, managing director of SCS GmbH).

Morita's 3D Accuitomo 170 is a DVT scanner especially designed for the specific requirements of dento-maxillo-facial diagnostics. Leading dentists and orthodontists have been very successfully using this scanner with different image sizes 40x40 mm, 60x60 mm and 80x80 mm). The

High precision + low dose = superior efficiency



even larger formats of 100x100 mm and 170x120 mm with voxel sizes of 80 μ m are ideally suited for non-dental applications. 360° rotation combined with upright sitting position of the patient ensure the high quality of the image data.


Low dose - high resolution

Images in landscape format are well suited for examinations of the entire head and neck area, particularly in implant therapy, restorative dentistry and surgery. Hard and soft tissue can be visualized exactly due to increased image dynamics. The 3D Accuitomo 170 is equipped with a special Flat Panel Detector (FDP) and offers 14 bit grayscale capabilities.


Furthermore, the systems features zoom reconstruction. Diagnostic data are quickly available – independent of data

volume. In any 80×80 mm image with a voxel size of $160 \mu m$ and 170×120 mm image with a voxel size of $250 \mu m$ details can be zoomed to a voxel size of $80 \mu m$. Thus follow-up detail imaging, which puts an additional burden on the patient's health and on the physician's time and financial resources, can be avoided.

"Combining extremely highly spatial resolution with extremely low radiation dose, digital volume tomography has the potential to replace conventional CT for examinations of the skull," Axel Hoppe predicts. While well accepted in dentistry, not even ten percent of the ENT specialists have adopted DVT technology. Nevertheless: High resolution and low dose are a convincing proposition. Paired with diagnostic competence – be it delivered by radiologists, dentists or any other specialised physician – DTV can be a most powerful tool.







The need for low-dose high-contrast images leads to an increasing demand for DVTs that can be met only by radiologists.





3D ACCUITOMO H17



CONSULTING SOLUTION EDUCATION SERVICE DIAGNOSTIC INFRASTRUCTURE



Standard equipment

- High-end capture PC
- High-end server PC
- Diagnostic display 30° TFT
- □ i-Dixel software
- 3-D image reconstruction
- 3-D real-time imaging
- Resticing
- Volume rendering
- One data viewer EX

DICOM RIS/PACS-Integration



Nine image formats


- o 40 mm × 40 mm, o 60 mm × 60 mm, o 80 mm × 80 mm, o 100 mm × 100 mm, o 170 mm × 120 mm
- Consistant image quality in large field of view
- High dynamic and high contrast imaging in hard and soft tissue
- Highly detailed 3-D digital images within 18 seconds

Patient chair

- Remote controlled XYZ-motion
- Dimensions: 1,62 m x 1,34 m x 2,18 m
- Weight: ca. 400 kg

COMPUTED TOMOGRAPHY

RAD-BOOK 2010



Mission Possible: Reducing Radiation Dose in CT

Over the past decade, Siemens has been a pioneer in creating a host of innovative technical features that significantly reduce radiation exposure in CT scans.

SOMATOM Sessions, the customer magazine from Siemens computed tomography, recently requested feedback on some of the most important of these innovations from physicians in Germany and the U.S.A. who have had experience with them.

By Catherine Carrington

ots of people talk about radiation dose and CT. But for more than a decade, Siemens Healthcare has made dose reduction a mission. The result: an impressive portfolio of innovations in scanner hardware, software, and imaging protocols that together have cut patient radiation exposure to a fraction of what it once was.

"Reducing radiation dose has always been a concern for Siemens," says Thomas Flohr, PhD, Director of CT Physics and Applications for Siemens Healthcare in Forchheim, Germany. "CT is the imaging modality of choice in many situations, and it would be used even more if not for the concern about radiation dose."

Siemens' focus was intensified in the late 1990s, when the company began to systematically search for new ways to reduce radiation dose. A timeline shows not only how relentless Siemens has been in pursuing this goal over the years, but also how creative Siemens Research & Development was. Key milestones include:

1994: Introduction of DOM, later extended to CARE Dose4D, a fully automated, real-time, anatomical dosemodulation technology that reduces radiation dose, depending on the area of the body, by 20% to 68% – without degrading image quality.

1997: Introduction of an ultra-fast ceramic (UFC) detector designed with a new gadolinium-oxy-sulfite scintillator. The UFC detector – still a key component

of multidetector and Dual Source CT systems – cut radiation dose by 30% when compared to previous generations of CT detectors.


1999: Introduction of ECG-pulsing, a technique that synchronizes tube current to the electrocardiogram. Used during spiral cardiac CT scans, ECG-pulsing maintains nominal tube current only during targeted phases of the cardiac cycle, markedly reducing tube current during phases that will not be used for image reconstruction. Dose savings: 30% to 50%.

2005: Introduction of the SOMATOM® Definition Dual Source CT scanner, which offers further dose efficiencies in cardiac CT through faster scanning, Adaptive ECG-Pulsing, and automated adaptation of table speed to heart rate. Dose savings: up to 50%, compared to single source CT.

2007: Introduction of the Adaptive Cardio Sequence, a prospective ECG-triggered "step and shoot" technique that reduces the average dose for CT coronary angiography to about 2.5 mSv.

2007: Introduction of the Adaptive Dose Shield, a technique of asymmetric collimator control that eliminates overscanning at the beginning and end of the CT spiral. Depending on the length of the scan, it reduces dose by 5% to 25%.

2008: Introduction of the SOMATOM Definition Flash CT scanner. With dual detectors and a table speed of



up to 45 cm/s, the Flash cuts radiation dose for coronary CT angiography to less than 1 mSv in many patients.

2008: Introduction of X-CARE, organ-based dose modulation that reduces output of the X-ray tube when it is directly in front of the breast and other dose-sensitive organs, such as the thyroid gland and eye lens. Reduces radiation dose to the breast by 30% to 40%.

2009: Introduction of Iterative Reconstruction in Image Space (IRIS). By "cleaning up" image noise, iterative reconstruction makes it possible to reduce radiation dose by up to 60% and still produce high-quality images.

Several dose-reduction strategies deserve special attention, including CARE Dose4D, the Adaptive Dose Shield, the SOMATOM Definition Flash CT scanner, and IRIS. Each of these is an example not only of Siemens' commitment to minimizing radiation exposure but also its track record of innovation.

CARE Dose4D

Determining the right tube current and, therefore, the right radiation dose, has always been crucial, says Marilyn J. Siegel, MD, Professor of Radiology and Pediatrics at the Mallinckrodt Institute of Radiology, Washington University School of Medicine, St. Louis, Missouri, USA.

But achieving that goal was much more difficult before CARE Dose4D, because adjustments in tube current had to be made empirically.

"CARE Dose4D has really been a great advantage for a number of reasons," Siegel says. "We get great image quality, reduced dose, and increased patient comfort. And it's automated, so it's easier for the technologist."

CARE Dose4D automatically adapts radiation dose to the size and shape of the patient, achieving optimal tube current modulation in two ways. First, tube current is varied on the basis of a topogram, by comparing the actual patient to a "standard-sized" patient. As might be expected, tube current is increased for larger patients and reduced for smaller patients. Differences in attenuation in distinct body regions are taken into account. For example, in an adult patient, 140 mAs might be needed in the shoulder region, whereas 55 mAs would be sufficient in the thorax, 110 mAs in the abdomen, and 130 mAs in the pelvis.

In addition, real-time angular dose modulation measures the actual attenuation in the patient during the scan and adjusts tube current accordingly – not only for different body regions, but also for different angles during rotation. This is particularly important in efficiently reducing dose in the shoulder and pelvic region, where the lateral attenuation is much

higher than the anterior-posterior attenuation.

Siemens has further refined this process with CARE Dose4D. Clinical experience has shown that the relationship between optimal tube current and patient size is not linear. Larger patients clearly need a higher dose than average-sized patients, but they also have more body fat, which increases tissue contrast. Smaller patients need a lower dose than averagesized patients, but they have less fat and less tissue contrast, which would result in noisy images if the dose were too low. Therefore, during real-time dose modulation, CARE Dose4D reduces radiation dose less than might be expected for smaller patients, while increasing the dose less than might be expected for larger patients. This maintains excellent diagnostic image quality while achieving an optimal radiation dose.

"CARE Dose4D is different from dose modulation approaches used by other vendors," says Flohr. "It uses measured attenuation data in real time, not just information from topograms; it makes use of a wide mA-range; and it can fine-tune dose on the fly."

Nowhere is CARE Dose4D more important than in pediatric imaging, where the risk associated with radiation exposure is many times higher in children than in adults. A child's smaller body bsorbs more of the radiation dose than does a

14

COMPUTED TOMOGRAPHY

RAD-BOOK 2010



"With children, you want as low a dose as possible but also excellent spatial resolution. CARE Dose4D allows us to reduce radiation exposure in all three planes without impairing diagnostic image quality."

Marilyn J. Siegel, MD, Professor of Radiology and Pediatrics at the Mallinckrodt Institute of Radiology, Washington University School of Medicine, St. Louis, Missouri, USA

larger body. In addition, the cancer induction risk is higher in children, because they have a longer lifespan ahead of them. At the same time, because children's anatomy is smaller it can be more difficult to visualize.

"With children, you want as low a dose as possible but also excellent spatial resolution," explains Siegel. "CARE Dose4D allows us to reduce radiation exposure in all three planes without impairing diagnostic image quality, and that's key."

Publications in scientific journals have shown that in adults, CARE Dose4D reduces radiation dose by 68% in the cervical spine, 37% in the lumbar spine, 30% in the thorax, and 38% in the abdomen and pelvis. In pediatric scans of the heart, a 58% dose reduction has been reported for CARE Dose4D.¹

Adaptive Dose Shield

In spiral CT, it is routine to do an extra half-rotation of the gantry before and after each scan, fully irradiating the detector throughout, even though only part of the acquired data is necessary.

As a result, the wide-cone beam exposes tissue that will never be part of reconstructed images. Until recently, no one gave much thought to this needless radiation exposure to patients. Such "overscanning" beyond the targeted scan range was simply accepted as an inevitable part of spiral CT.

Siemens took a fresh look at the problem and, in 2007, introduced the Adaptive Dose Shield, a technology based on precise, fast, and independent movement of both collimator blades. Instead of exposing patients to unnecessary radiation, the Adaptive Dose Shield asymmetrically opens and closes collimators at the beginning and end of each scan, tempo-rarily blocking those parts of the X-ray beam that are not used for image reconstruc-



"There's more and more awareness about the amount of radiation used for CT scanning. Siemens has thoroughly looked into this and is one of the first vendors to implement the tools we need to improve our scanning."

Christoph Becker, MD, Professor of Radiology and Section Chief of CT and PET/CT at Munich University Hospital in Germany:

tion. As a result, only the targeted tissue is irradiated. Like many other dose-saving innovations, it is a feature pioneered by Siemens.

"There's more and more awareness about the amount of radiation used for CT scanning," says Christoph Becker, MD, Professor of Radiology and Section Chief of CT and PET/CT at Munich University Hospital in Germany. "Siemens



"This ultra-low dose wasnever possible before, but with SOMATOM Definition Flash – with its high temporal resolution and improvements in the X-ray tube and detector – it is now possible." Jörg Hausleiter, MD, Cardiologist, Associate Professor of Medicine, German Heart Center, Munich, Germany

has thoroughly looked into this and is one of the first vendors to implement the tools we need to improve our scanning."


At Munich University Hospital, Becker has two Siemens scanners equipped with an Adaptive Dose Shield, the SOMATOM Definition AS+ and the SOMATOM Definition Flash. Although the Adaptive Dose Shield reduces the radiation dose in every study, the savings are especially notable over shorter scan ranges. Dose savings can reach 25% or more in cardiac imaging, for example.


The Adaptive Dose Shield is especially well suited to pediatric imaging. "In any circumstance in which children have to be investigated, I would always prefer to use a scanner with the Adaptive Dose Shield," Becker says. "It's always on, and it always reduces the radiation dose."

Flash

At the German Heart Center, Jörg Hausleiter, MD, has been using a SOMATOM Definition Flash CT scanner since April. With this revolutionary scanner, he can image the heart in a quarter of a single heart beat. Equally impressive, he has been able to achieve a radiation dose of 1mSv or less in a large proportion of patients undergoing CT coronary angiography. "That's unbeatable compared to other CT scanners," says Hausleiter, an Associate Professor of Medicine at the Munich-based hospital.

The SOMATOM Definition Flash gets its name from its flash-fast speed. Equipped with two detectors, two X-ray sources, and a gantry that rotates in 0.28





1: Instead of just taking into account the patient's external dimensions and apparent size, CARE Dose4D analyzes the cross-sectional anatomy in real-time and adjusts the emitted X-ray dose accordingly – providing excellent image quality with minimized exposure.

seconds, the scanner boasts a temporal resolution of just 75 ms. Moreover, thanks to an innovation unique to the SOMATOM Definition Flash, the patient table no longer slowly inches forward during scanning. Instead, in low-dose Flash Spiral mode, the table can glide along at 45 cm/s while the scanner integrates data from both detectors, achieving a gap-free scan even though each spiral is wide open.


Still, according to Hausleiter, the key question is whether excellent image quality can be achieved at such a high scan speed and low dose. With the SOMATOM Definition Flash, the answer is clearly yes. "This ultra-low dose was never possible before, but with this scanner – with its high temporal resolution and improvements in the X-ray tube and detector – it is now possible," he says.


Of the first 100 coronary CT scans performed on the Definition Flash at the German Heart Center, more than 70% could be done in Flash mode. As a result, the average radiation dose for all coronary CT scans – including longer scans needed for presurgical evaluation and triple rule-out studies – dropped from a median of 5 to 7 mSy down to 1.8 mSy.

Of the 70% of patients scanned in Flash mode, approximately half could be scanned at 100 kV. (In general, a tube voltage of 100 kV is suitable for patients with a body mass index of less than 30 or a body weight of less than 90 kg). In these patients, Hausleiter found that the median radiation dose was just 1 mSv. The other half of the patients were scanned at 120 kV, and received a radiation dose of 1.6 to 1.8 mSv, still far lower than the typical radiation dose for coronary CT angiography.

The PROTECTION I study highlights how much progress has been made. In 2007, Hausleiter and an international group of researchers from 50 medical centers set out to determine the typical radiation dose for patients undergoing coro-







2: Pediatric imaging: no breath hold and no anesthesia was necessary for the scan with 0.37s scan-time by using only 1 mSv (Fig. 2A); Split-second thorax scan by using only 1.65 mSv (Fig. 2B and 2C).


nary CT angiography, using CT scanners manufactured by a variety of vendors. Published in the February 4, 2009, issue of JAMA, the study showed that the median dose was 12 mSv.

"It's important to realize the large steps we've taken," says Hausleiter. "The dose we can achieve today is one-tenth of what it was in the PROTECTION I study. That's a major improvement."

Such a low radiation dose could expand CT's horizons in the evaluation of heart disease. For example, for patients with high heart rates and irregular heart rhythms, the "step and shoot" Adaptive

COMPUTED TOMOGRAPHY

RAD-BOOK 2010



3: Image cardio sequence: Fully flexible X-ray pulsing in combination with 75ms temporal resolution results in low dose cardio scan (0.36 mSv dose).

Cardio Sequence, with prospective ECG-triggering and arrhythmia detection, is ideal and keeps radiation dose to about 2.5 mSv. For patients with reasonably low and stable heart rates, the Flash Spiral is the method of choice. But even for patients with mild arrhythmia, Hausleiter thinks the Flash mode, which captures all necessary data in a single heart beat, may be fast enough to do the job, and at a radiation dose of 1mSv. And, if that one heart beat happens to be an extra unwanted beat generated by the arrhythmia, the Flash's low radiation dose means there is little risk in repeating the study.

A radiation dose of below 1 mSv also raises the possibility of using CT for screening patients at risk for heart disease. "We need to start thinking about that question," Hausleiter says. "With coronary CTA, we would gain information on calcification, the location of plaques, and the presence of noncalcified plaques – the type we really worry about. In the end, screening could reduce the number of heart attacks."

Image data reconstruction of an abdominal scan with Standard FBP at full dose (Fig. 4a) and scanned at 60% lower dose while reconstructed with Iterative Reconstruction in Image Space (Fig. 4b). Despite the fact that Fig. 4B was acquired at significantly lower dose it shows the same low noise compared to the standard FBP at full dose.

Iterative Reconstruction

Iterative reconstruction, which Siemens is slated to debut at the 2009 RSNA meeting in Chicago, is the latest success story in the company's mission to reduce radiation dose. Essentially, iterative reconstruction introduces a correction loop in the image generation process that cleans up artifacts and noise in low-dose images.

Other vendors are working on iterative reconstruction, but Siemens has developed a unique method. A typical approach to iterative reconstruction is to measure data in the reconstructed image and compare it to the original data, using differences to identify ways to improve the image. This approach is time-consuming because, with each iteration, new measurement data must be calculated.

Siemens instead takes the original data and reconstructs a super-high-resolution image. The image is very noisy, because the filtering that ordinarily reduces image noise is not used, in order to avoid any loss of information. Then prior knowledge of the scanned object is used to smooth the image and reduce noise within homogeneous regions, while contrast edges are preserved. This process is repeated over several steps, or iterations.

"Why is Siemens' approach better?" Because we start with a super-high-resolution image and clean it up," says Thomas Flohr. "We can fine-tune the



"With Siemens Iterative Reconstruction I can save up to 60% dose for wide range of routine applications while maintaining excellent image quality" says U. J. Schoepf, MD, Professor of Radiology and Cardiology and Director of CT Research and Development at the Medical University of South Carolina.


process, so we don't lose object information. We maintain image texture that is familiar to readers, so the resulting image looks like a standard CT image and doesn't have the plastic-like look that is often the drawback of other iterative reconstruction approaches. And the process is very fast and efficient." Most important, Siemens' iterative reconstruction technique can reduce radiation dose by up to 60%, depending on the body region and the original scan dose.


Future Directions


The next automated tool for dose reduction is likely to be automatic kV adaptation to the patient's size and the exami-











5: Single Source CT requires slower table feeds to prevent gaps in the acquired volume (top, center). Dual Source CT combines the data from 2 detectors for faster table feeds above a pitch of 3 (bottom).

nation type. Researchers are beginning to understand and further evaluate its effect on image quality and dose. In the PROTECTION II study, for example, Hausleiter and his colleagues randomly assigned 400 patients to undergo coronary CT angio-graphy with either a 100 kV protocol or the more conventional 120 kV protocol. Reported at the 2009 American College of Cardiology Annual Scientific Session, the study showed that the use of 100 kV reduced radiation dose by 31%, while image quality scores were virtually identical. "This proves you can use 100 kV very liberally when looking at the coronary arteries," says Hausleiter. At the Mallinckrodt Institute, Siegel has also been evaluating the radiation savings possible through use of a lower tube voltage. Her work with Lucite phantoms that simulate the size of various body regions in children has shown that at a tube voltage of 80 kV the radiation dose is reduced when compared to a tube voltage of 140 kV, even when the tube current is increased to ensure good image quality.


At the St. Louis Children's Hospital, Siegel has been using the SOMATOM Definition AS 64-slice CT scanner to scan pediatric patients. She will continue her research in pediatric phantoms and in patients with this newer generation scanner to determine the impact on radiation dose and image quality of modulating kV. She anticipates that with this newer-generation scanner, the quality of CT studies will improve even further as radiation dose is decreased. "There is an old saying, 'Beautiful pictures come at the cost of higher radiation dose,' "Siegel says.

"We've already disproved that, and we intend to further disprove it."

6: To accelerate the convergence of the reconstruction IRIS applies the raw data reconstruction only once. During this newly developed initial raw data reconstruction a so called master image is generated that contains the full amount of raw data information. The following iterative corrections known from true iterative reconstruction are consecutively performed in the image space. In addition, the noise texture of the images is comparable to standard well-established convolution kernels. The new technique results in artifact and noise reduction, increased image sharpness and dose savings up to 60% for a wide range of clinical applications.

Medical writer Catherine Carrington holds a master's degree in journalism from the University of California Berkeley and is based in Vallejo, California.

1 Mulkens et al.: Use of an Automatic Exposure Control Mechanism for Dose Optimization in Multi-Detector Row CT Examinations: Clinical Evaluation, Medical Physics.



RAD-BOOK 2010

Low dose and Robustness of Aquilion ONE for Cardiac CT

fter ten years of research and development, Toshiba released the Aquilion ONE™ CT scanner at the RSNA 2007. This scanner was able to cover 16cm in the 'z' direction using an array of 320 x 0.5mm detectors. This was, and still is the only CT scanner in the world capable of capturing images of the entire heart and coronary arteries in a single rotation without any need for patient motion during the scan.

This paper reports huge robustness of the system. Furthermore, it is now possible to provide diagnostic images of the heart and coronary arteries at low radiation dose in all patients. Figure 1 shows the principal of generating slice data from a wide volume acquisition. Double number of slices (640) can be generated out of volume data without increasing radiation dose as this is done in the reconstruction¹. This is as opposed to the z-sampling technique which increases the exposure dose due to larger penumbras.

Retrospective Spiral Technique

This is the original method used in cardiac CT which involves a spiral acquisition using a very low pitch in order to acquire an entire image dataset of the heart and coronary arteries at all phases of the cardiac cycle. A major drawback of this technique is radiation dose which is commonly above 10 mSv, even when using ECG dose modulation techniques (which reduces tube current during systole and increases it during diastole usually the phase at which the heart is stillest and the highest quality images are obtained).

Prospective-gating ('step and shoot')

Prospectively-gated techniques have been introduced on all commercially available 64-MDCT in an attempt to address these radiation issues. This involves acquiring an axial block of data at end-diastole using the full width of the 64 debeats with the scanner moving between each acquisition. These separate blocks of data are acquired with some overlap allowing the images to be 'stitched' together afterwards in an attempt to create a single dataset of the entire heart and coronary arteries. Although this technique typically reduces radiation doses to the 5-10mSv range, this comes at the price of robustness as the data has been acquired at multiple separate timepoints. Conventional prospective gating techniques require a low resting heart rate (ideally <60bpm) and are very sensitive to poor breath-holding and to any variation in heart rate during the study. They are completely unsuitable for patients with multiple ectopic beats or those in AF.

Aquilion ONE, in contrast, uses a unique prospective technique in which the entire dataset is acquired in one heartbeat. This provides huge advantages compared with conventional techniques both in terms of radiation dose and robustness.


$HR \leq 65bpm$

A single-beat prospective radiation pulse is chosen that enables data to be acquired between 70% and 80% of the R-R interval. This technique is used without patient selection and provides reliably excellent images of the heart and coronary arteries at radiation doses which are typically below 3 mSv as reported by Dr. Bull².

"Using the conversion factor k = 0.014mSv/mGy/cm (ICRP102) the estimated patient radiation exposure was E = 0.887mSv averaged for all patients. An increase with patient's BMI was observed with a range between 0.22 mSv and 1.57 mSv", reported by Lembcke et al³.

65bpm < HR < 75 bpm

This technique is used for patients in which beta blockers are either contra-indicated or not completely effective. The radiation pulse is widened in this scenario to enable data to be acquired in a single beat between 30% and 80% of the R-R interval. This option is unique to the Aquilion ONE as it enables end-systolic images to be reconstructed in addition to conventional end-diastolic images still using prospective ECG gating in a single heart beat. End-systolic isovolumic relaxation is



Robustness

All conventional spiral CT systems 'build-up' an image of the heart and coronary arteries over multiple heart-beats. This can be achieved in two ways²:

tectors (3.2-4cm coverage). As this coverage is insufficient to cover the entire heart, the image is again built up using multiple (typically 4-5) similar acquisitions performed over multiple heart

Fig. 2:
Arrhythmia
rejection: the
system skips
the irregular
RR-intervals
until it finds
the best one
(green line)
within clinical justified
time level





Fig. 3:
AIDR was
introduced at
the RSNA
2009 to save
the radiation
dose up to
60% for
a wide range
of clinical
applications.



Adaptive Iterative Dose Reduction

Although low radiation dose as low as 0.22 mSv has been reported by Charite Cardiac group, Toshiba keeps working on further dose reduction. At the RSNA 2009 Adaptive Iterative Dose Reduction (AIDR) has been introduced which result in dose saving up to 60% (Fig. 3).

Most importantly is this iterative dose reduction algorithm does not blur the edges nor change the texture but only suppresses the noise.

References

- 1 J. Blobel, H de Vries, R. Irwan, J Mews, Y Ogawa, "640 Multislice Reconstruction with Dynamic Volume CT", VISIONS 13:2009
- 2 R. Bull, "Robustness of Aquilion ONE for cardiac CT", VISIONS 15:2010
- 3 A. Lembcke, P. Rogalla, J. Mews, J. Blobel, "Dose Optimized Single Beat Coronary CT Angiography", VISIONS 15:2010

relatively short but unlike end-diastolic relaxation this is 'fixed' and does not shorten with increasing heart rate. Radiation doses are typically below 5 mSv using this technique with excellent image quality^{2,3}.


$HR \ge 75bpm$

A feasible option is to use a single-beat prospective technique acquiring data throughout a single entire cardiac cycle. This enables data to be reconstructed at any phase of the cardiac cycle. The reconstruction options (0-100% of R-R interval) are similar to conventional spiral 'retrospective' techniques but at greatly reduced radiation dose as this is still a prospective technique with no 'oversampling' required (typical radiation doses are 4-6 mSv with good image quality)².


Arrhythmias rejection

The Aquilion ONE's unique ability to acquire an entire cardiac dataset in a single beat with no table motion has particular advantages in arrhythmias. As long as ventricular rate is reasonably well controlled (approximately < 75 bpm) reliably excellent images are acquired in all patients with no motion artefact and no step artefacts. In effect the Aquilion ONE is able to 'wait' for a sufficiently long R-R interval in which to deliver its prospective radiation pulse.

Clinical Cases



Case 1 55 year old woman with atypical chest pain. Negative exercise test but continued clinical concern. Single rotation, 70%-80% acquisition. Images show critical LAD stenosis subsequently treated with PCI.



Case 2
80 year old man. Persistent fever post aortic valve replacement.
Two step protocol covering whole chest, HR 80bpm. Each step single rotation, reconstructed at 40% of R-R interval. Images show large false aneurysm secondary to aortic root abscess.

Source: Dr Russel Bull, Royal Bournemouth Hospital UK

LightSpeed Channels Power in kW Coverage/ rotation in mm VCT XT 100 40 mm isotropic per 0.4 s (Cardiac: 0.35 s), 0.35 mm resolution incl. SnapShot Pulse &

VCT Pro 100 40 mm isotropic per 0.4 s, (Cardiac: 0.35 s) 0.35 mm resolution

VCT 64 85 40 mm isotropic per 0.4 s, 0.35 mm resolution VCT XTe 85 40 mm isotropic per 0.4 s, 0.35 mm resolution incl. SnapShot Pulse, Volume Helical Shuttle & ASiR

VCT Select 32 85 40 mm, 20 mm isotropic per 0.4 s, 0.35 mm resolution



▶ Highlights

▶ GE Healthcare LightSpeed Series

- VCT XTe: Up to 35 frames reconstruction / sec VCT XTe: Up to 35 frames reconstruction / sec VCT XTe: Adaptive Statistical Iterative Reconstruction (ASiR) VCT XTe: Volume HelicalShuttle up to 500 effective slices
- VCT XT: 5-beat cardiac acquisition, up to 70% dose reduction
- VCT XT: from one to 6 mSv for cardiac CTA
 VCT XT: expanded coverage to 80 mm for CT perfusion
- VCT: 5-beat cardiac imaging: resolution, coverage, fast scanning VCT: a true 64-channel DAS

▶ GE Healthcare BrightSpeed Series BS 8 **BS 16** BS 16 Elite Channels 16 8 16 16 Power in kW 42 53.2 42 42 2 x 0.625, Coverage/ 2 x 0.625, 16 x 0.625 16 x 0.625 rotation in mm 4 x 1.25, 8 x 1.25 & & 16 x 1.25 & 16 x 1.25 4 x 2.5 & 8 x 2.5 mm mm per 0.8 s mm per 0.5 s 4 x 5 mm per 0.8 s per 0.8 s ▶ Highlights - LightSpeed VCT Technology inside - BS 4: perform long

- coverage and high grade CT-A
- -BS 8: any organ in a breath-hold
- -BS 16: Sub-mm microvoxels for incredible detail
- -BS 16 Elite: faster routine scanning with variable speed rotations



▶ GE Healthcare HiSpeed CT Dual

Channels Power in kW Coverage/ rotation in mm

2 x 10 mm

24

HiSpeed CT/e Dual



Highlights

- Xtream productivity
- As easy as 1-2-3 with GE's SmartGantry
- Patient concious designMulti-slice scanning for everyone and everywhere
- Increased speed and greater clinical flexibility

▶ GE Healthcare LightSpeed RT

Channels Power in kW Coverage/ rotation in mm LightSpeed RT 16

100 16 x 1.25 mm

per 0.8 s

LightSpeed RT 16 Pro

16 x 0.625 mm per 0.5 s



- Multi-slice wide bore radiation oncology CT scanner
- Large 80 cm opening
- -65 cm field of view - Full RT connectivity
- Complete radiotherapy simulation solution



▶ GE Healthcare Discovery CT-750 HD


Views / s > 7.000 Power in kW 150 **Resolution** | 230 Microns

See more, know more, less dose New Gemstone Detector plus complete new imaging chain

– 230 Micron Resolution over 2 m coverage

– Dual Energy with 1 Tube and 500 cm FOV

- Volume HelicalShuttle up to 500 effective slices
 Adaptive Statistical Iterative Reconstruction (ASIR)
- Dose reduction up to 50% over the intire body up to 83% for cardiac examination
- -Up to 35 frames reconstruction / sec



- Sub-second, real-time image reconstruction
- Minimum scan time 0.8 sec and maximum field of view 500 mm
- Preventive examination supported by fatPointer or riskPointer
- -Straight forward patient registration and easy system handling

Morita 3D Accuitomo 170

Voxelsize 80µm Scan time 18 s Scan-Volumina

Ø40 x H40 mm, Ø60 x H60 mm, Ø80 x H80 mm. Ø100 x H100 mm,

Ø170 x H120 mm



▶ Highlights

- Offers high definition 3D-CT images with low patient dose
- Displays both hard and soft tissue
- A wide dynamic range and precise grayscale differentiation capability
- Enables comprehensive examination for diagnosing temporal bone, paranasal sinuses, mandible, skull base, etc.
- Compact Floor space: 1,620 mm X 1,200 mm

Brilliance iCT

▶ NeuroLogica CereTom® Portable CT-Scanner

Slices Power 1.4 kW Adjustable slice

thickness

1.25, 2.5, 5, 10 mm



- Rotation speed 0,5 sec KV Range 80-140 kV at 7,5 mA
- Patient dose CTDI 41mGy
- Field of view 25 cm
 Images compatible with surgical navigation systems
- -Wireless connectivity to PACS and DICOM 3 compliant with modality worklist
- Intraoperative scanning capabilities with DORO® CERETOM® Intraoperative Cranial Stabilisation system
- Immediate 2D, 3D and MPR images with custom pre-set protocols
- Advanced visualisation solutions powered by Barco's Voxar 3TTM

▶ Philips Brilliance iCT scanners

Channels Power in kW Coverage/ rotation in mm

256 120 80

Brilliance iCT SP 128 100



▶ Highlights

- Enhanced performance for routine and emerging applications
- Patient-specific acquisition protocols to balance image quality and dose utility
- Revolutionary AirGlide Gantry for whisper-quiet performance at 220 rpms
- Exclusive dose-saving features like the Eclipse Collimator, Step & Shoot Cardiac and Dedicated Pediatric Protocols
- Life-cycle benefits through a scalable hardware and software platform

▶ Philips Brilliance CT 64

Channels Power in kW Coverage/ rotation in mm Brilliance 64 with Essence technology



Highlights

- Myocardial perfusion, CTA and whole brain perfusion
- -CT stroke assessment, 80 mm coverage through Jog Scan
- Brilliance CT workspace user environment improves productivity by working the way you do
- DoseWise design delivers optimal dose efficiency without compromising image quality
- Submillimeter isotropic accuracy
- Essence technology with the latest x-ray tube, detector and reconstruction advancements improving image quality

▶ Philips Brilliance CT 16

Channels Power in kW Coverage/ rotation in mm

48 24





▶ Highlights

- Advanced performance systems for routine imaging needs with fast acquisition and high quality image results

 Diagnoses of small lesions with submillimeter slices
- Brilliance CT workspace user environment improves productivity by working the way you do
- DoseWise design delivers optimal dose efficiency, without compromising image quality
- Scalable platform for growth and future applications, making it a secure, long-term investment







One beat cardiac imaging

Aquilion ONE: the world's first dynamic volume CT



Aguilion ONE

Toshiba Medical System's Aquillon ONE is a quantum leap in CT imaging that can capture a 3D image of the heart in just one beat.

The wide coverage provided by the Aquillion ONE's 16cm detector, which has 320 detector rows, can scan the brain or heart in less than a second. So you can see an entire organ in 3D with perfect continuity along the z-axis. Or see it in 4D, moving as time passes. Or see it extremely fast, with a lower contrast medium dose and exposure dose.

The Aquillon QNE will bring you dynamic views of the body you could not see before.

See it for yourself at www.toshiba-medical.eu.


Toshiba: Made for Patients, Made for You, Made for Life!



Philips MX CT scanners MX 4000 Dual MX 6000 Dual Channels Power in kW 28 42 Coverage/ 20 20 rotation in mm ▶ Highlights - Multislice technology doubles performance for hard-to-image anatomy Advanced applications right at console: 3D, MIP, MPR, volume rendering, CTA, virtual endoscopy, etc. – 0.8* sec rotation; 0.9* sec reconstruction; 0.8* mm slices

- People-centric design improves user productivity and patient

comfort





▶ Siemens Somatom Definition Flash Channels 2 x 128 slices 200 kW Power in kW 48 cm (with A4DS) Coverage/ rotation in mm 0.24 mm resolution ▶ Highlights Flash speed. Lowest Dose. Split-second thorax imaging without the need for breath hold

for beta-blockers - Full cardiac detail at half the dose required for a

conventional single

source CT

- Widest detector coverage in 16-slice class
 One of the industry's smallest site requirements at 18 square meters
- Fully compatible with Brilliance Workspace, Extended Brilliance Workspace and the Brilliance Workspace Portal
- Sub-mSv heart scanning to cover the entire heart in only 250 ms Single dose Dual Energy for a 2nd contrast in daily routine
- Iterative Reconstruction in Image Space delivers up to 60% additional dose reduction or significantly improved image quality

▶ Siemens Somatom Definition AS+

▶ Siemens Somatom Definition Channels 2 x 64 slices Power in kW 160 kW Coverage/ 20 cm (with A4DS) rotation in mm 0.24 mm resolution ▶ Highlights The world's first Dual Source CT - Faster than every beating heart without the need

One-stop diagnosis in acute care for acute chest pain,

abdominal pain, trauma, and neuro imaging
Beyond visualization with Dual Energy by characterizing,

highlighting and quantifying material for the first time

Power in kW rotation in mm

Channels

Coverage/

128 slices

27 cm (with A4DS)

0.24 mm resolution

100 kW

▶ Highlights -The world's first adaptive scanner Adapts to any patient, e.g. cardiac, pediatric, obese, trauma, intervention with 78 cm bore, up to 300 kg patient table load and up to 2 m scan range

– Eliminates over-radiation in every spiral scan, additionally reducing the dose by up to 25% - Iterative Reconstruction in Image Space delivers up to 60% additional dose reduction or significantly improved image quality

-Allows whole organ perfusion studies and 4D dynamic studies with up to 27 cm covarage



Then look at the latest edition of the most popular CT in the world – the SOMATOM Emotion Excel Edition

Siemens' SOMATOM Emotion Excel Edition incorporates all of Siemens' renowned creativity and technical competence to combine highest 16-slice CT performance in an affordable package. Our research competence is your opportunity to own leading technology with a substantially lower total cost of ownership.

www.siemens.com/healthcare +49 69 797 6420

Answers for life.



▶ Siemens Somatom Definition AS

Channels Power in kW Coverage/ rotation in mm

64 slices 40 slices 100kW 80 kW 67 mm (with A4DS) 0.24 mm resolution



▶ Highlights

- The world's first adaptive scanner
- Adapts to any patient, e.g. cardiac, pediatric, obese, trauma, intervention with 78 cm bore, up to 300 kg patient table load
- and up to 2 m scan range

 Eliminates over-radiation in every spiral scan, additionally reducing the dose by up to 13%
- Enables 3D guidance in minimally invasive procedures with Siemens' unique Adaptive 3D Intervention Suite
- Adapts to your space with on-site upgradeability from 40/64 to 128 slice configuration

▶ Siemens Somatom Definition AS Excel Edition

Channels Power in kW Coverage/ rotation in mm



▶ Highlights

- Affordable Performance
- Adapts to your space with on-site upgradeability from Excel Edition to AS 64 or 128 slice configuration.

▶ Siemens Somatom Sensation

Channels Power in kW Coverage/ rotation in mm

64 slices 80 kW 28.8 mm 0.24 mm resolution 40 slices 70 kW 28.8 mm

0.33 mm resolution



▶ Highlights

- -Trusted performance. Without compromise.
- Speed, resolution, coverage and low dose without compromise
- Industry's highest isotropic resolution
- -0.33 second rotation time

▶ Siemens Somatom Sensation Open

Channels Power in kW Coverage/

40 slices 50 kW 28.8 mm

rotation in mm 0.33 mm resolution



▶ Highlights

- -82 cm large bore with 82 cm FOV
- 280 kg high capacity patient table
- Designed for RT planning
- Available as sliding gantry for trauma and inter-operative installations.

▶ Siemens Somatom Definition AS 20

Channels Power in kW Coverage/

20 slices 80kW

rotation in mm 67 mm (with A4DS)



▶ Highlights

- The world's first adaptive scanner
- Adapts to any patient, e.g. cardiac, pediatric, obese, trauma, intervention with 78 cm bore, up to 300 kg patient table load and up to 2 m scan range
- Eliminates over-radiation in every spiral scan, additionally reducing the dose
- Enables 3D guidance in minimally invasive procedures with Siemens' unique Adaptive 3D Intervention Suite
- Adapts to your space with on-site upgradeability from 20 to 40,64 or 128 slice configuration

▶ Siemens Somatom Emotion 6/16

Channels Power in kW Coverage/ rotation in mm

16 slices 6 slices 50 kW 42 kW



▶ Highlights

- The Most Popular CT in the World
- Perfection in image detail with the smallest tube focal spot size and up to 68% dose reduction with CAREDose 4D
- Clinical efficiency simplified with the breathing indicator, CT storage box in gantry, up to 16 images per second recon and remote access to acquisiton workplace Savings in every scan with the smallest space required for
- installation and lowest power and air-conditioning requirements

What inspired our iCT innovation? The healthcare needs of over 6.7 billion people.

Disease doesn't discriminate so the Philips Brilliance iCT doesn't either. Capture crisp images from all over the body,



from underweight infants to overweight adults. Be confident that you are caring for patients with the quality they deserve - a low dose without sacrificing image quality. Visit us at www.philips.com/CT or call 1-800-229-6417.

*Because our innovations are inspired by you.





▶ Siemens Somatom Emotion Excel Edition

Channels Power in kW Coverage/ rotation in mm 16 slices 50 kW



▶ Highlights

- -The Most Popular CT in the World now even more affordable
- Perfection in image detail with the smallest tube focal spot size and up to 68% dose reduction with CAREDose 4D
- Clinical efficiency simplified with the breathing indicator and CT storage box in gantry and up to 8 images per second reconstruction
- -Savings in every scan with the smallest space required for installation and lowest power and air-conditioning requirements

▶ Siemens Somatom Spirit Channels 2 slices Power in kW 40 kW Coverage/ rotation in mm 10 mm

▶ Highlights

- Join the World of CT
- New technology, more performance, less cost
 Easy user interface
- New level of cost-effectiveness
- New dimension in patient-friendliness



- Lateral table movement*
 Isophasic ONE-beat cardiac imaging in ONE rotation
- Isophasic 16 cm dynamic volume imaging with 20 volumes/s
- Whole organ perfusion without table movement, e.g. brain,
- liver, pancreas, kidney, etc.*

 Dual Energy with 50 cm Field of View*
- -8 cm ultra helical with outstanding image quality
 -0,5 mm detector technology with best low contrast resolution 2 mm @ 3HU
- Morphology and moving joints



- Cardiac volume imaging
 8 cm ultra helical with outstanding image quality

- Wide area organ and heart perfusion*
 Dual Energy with 50 cm Field of View*
 Isophasic 8 cm dynamic volume imaging with 20 volumes/s
- 4D-Volume shuttle technology for areas larger 8 cm
- Lateral table movement*
- -0.5 mm detector technology
- with best low contrast resolution 2 mm @ 3HU



- SureCardio-Prospective for helical cardiac imaging with lowest dose (option)

 - Up to 35 ms temporal resolution (option)
- Patient specific automatic optimization of cardiac scan parameter
- -Variable helical pitch to combine two scans in one run (e.g. ECG-gated and Run-off)
 -0,5 mm detector technology with best
- low contrast resolution 2 mm @ 3HU



- SureCardio-Prospective for helical cardiac imaging with lowest dose (option)

 up to 35 ms temporal resolution (option)
- Patient specific automatic optimization of cardiac scan parameter
- Variable helical pitch to combine two scans in one run
- (e.g. ECG-gated and Run-off) (option)
 -0,5 mm detector technology with best low contrast resolution 2 mm @ 3HU


Balance image quality and dose like never before.

Today, with breakthrough ASiR** technology from GE, you can lower patient dose in CT imaging without sacrificing image quality. So while leading-edge clarity lets you diagnose with confidence, a reduced dose enables better patient care. That's our kind of balance—and thanks to ASiR, only GE delivers it.

Learn more at gehealthcare.com/lowdoseCT







- CT DSA with SureSubtraction (option)
- Helical cardiac imaging, incl. cardiac CTA
- up to 40 ms temporal resolution (option) - Patient specific automatic optimization of cardiac scan parameter -Advanced 3D, automatic bone removal, etc. - Ultra low dose scanning

−0,5 mm detector technology with best low contrast resolution 2 mm @ 3HU

▶ Toshiba Aquilion Large Bore Slices 16 Coverage/rotation 3,2 cm Rotation speed 0,5 s ▶ Highlights 4D CT - for respiratory triggered radio therapy

- 90 cm gantry bore - 85 cm extended field of view (option) - 70 cm standard field of view


- Superior homogeneity for high precision

radiation therapy planning

- 0,35 mm isotropic spatial resolution

- 0,5 mm detector technology with best low contrast resolution 2 mm @ 3HU


Real time multislice open bore fluoroscopy



- 0,35 mm isotropic spatial resolution - Easy "Ready-Set-Go" user concept - Ultra low dose scanning

- Real time multislice fluoroscopy

-0,5 mm detector technology with best low contrast resolution 2 mm @ 3HU



- Powerful 3D software with auto bone removal - Easy "Ready-Set-Go" user concept

- Ultra low dose scanning

- Real time multislice fluoroscopy

-0,5 mm detector technology with best low contrast resolution 2 mm




Alliance Medical – flexible diagnostic imaging services



Highlights

- Static diagnostic imaging centers MRI, CT, PET, PET/CT Interim services for bridging downtimes
- Regular "routing" services

▶ Hologic ImageChecker D



Highlights

- Industry GOLD standard for CAD
- Highest sensitivity at comparable false marker rates
 Seamless integration of ImageChecker and review workstation
 Powered by R2 technology,
- the most trusted and fieldproven CAD technology
- Customized CAD results

▶ IBA Dosimetry Dosimax plus A HV

Dosimeter for measuring simultaneously dose, dose rate, exposure time and dose length product



▶ Highlights

- Designed according to IEC 61674
- For use with solid state detectors or ionization chambers
- For CDTI determination in combination with head and body phantom

▶ PTW CT Dosimetry

Quality control equipment for CT dose measurements



▶ Highlights

- Combined head & body phantom available
- CT chamber for precise dose length product measurements and CTDI determination

RTI Electronics CT Slice Probe



The CT Slice Probe is designed to make CTDI measurement more exact and has also the ability to further analyze the result. Following parameters are achived from a single exposure: CTDI100, CTDIvol, DCTIw, CT dose profile, DLP, Point Dose, Performance of the AEC, FWHM and Scatter Index.




- ► **Highlights** All in One Shot
- Quick and Simple Set up
- Accurate and Sensitive
- -No limitations due to the beam width

▶ TeraRecon Aquarius

Product Technology Resolution

Aquarius workstation 3D diagnostic workstation

AquariusNET server Client-server 3D architecture



▶ Highlights

- -Thin client-server solution enterprise-wide
- VolumePro: handling many 3D sessions at once
 Rendering on central server, results streamed to PC
 Fast and efficient in the reading workflow
- AquariusNET runs on almost any standard PC

CT ACCESSORIES

Tomovation - Modular building solutions



▶ Highlights

-Engineering, rental, sale of modular buildings MRI, CT, PET, PET/CT including or excluding diagnostic equipment

▶ ulrich medical – CO₂ Insufflator for virtual coloscopy

Pressure Insufflation rate Setting 0-30 mmHG, infinitely variable, preselectable

1-4 l/min, arbitrary

supported by voice confirmation system



▶ Highlights

Automatic insufflation of CO2 into the colon for virtual coloscopy examinations in CT

- Significant improvement of diagnostic results compared to manual room air insufflation

- Increase of patient comfort due to automatic adjustment of over pressure and faster resorption
- Easy setting of gas volume and pressure
- Display of gas consumption
- Four adjustable flow rates

▶ Ziosoft Ziostation

Full Suite of Clinical Applications

3D, 4D, MIP, MPR, CPR, Multi-Modality Viewer

▶ Highlights

- Multi-volume fusion
- Vessel analysis
- Coronary analysis
- Colon analysis
- PET/CT viewer
- CT Cardiac Function Analysis
- CT Brain Perfusion

- Calcium Scoring
 Reporting functions
 Customizable button palettes and keyboard/mouse shortcuts to configure the system individually
- Automatic segmentation algorithms employ anatomical analysis, targeting specific objects to maximize accuracy

 – Full Web-Based Capabilities

 – Action macros perform multiple processes with one or zero
- button clicks

RAD-BOOK 2010

The smart solution

for QA in Radiodiagnostics



Multimeter MagicMax + Test device Primus

Your best choice for constancy and acceptance tests!

MagicMax Multimeter is a Dosimeter and kV-meter in one unit.

Primus Image Quality Test Device at digital / conventional fluoroscopic and radiographic X-ray units.



www.iba-dosimetry.com info@iba-dosimetry.com



Visit us at the ECR in Vienna 4 to 8 March 2010, Booth 305, Expo 10

MAGLIFE Serenity

PERFECTION IN MRI MONITORING -SAFE, USER-FRIENDLY, MOBILE

The MAGLIFE Serenity provides the guarantee for highest safety:

For years, SCHILLER has been one of the world's leading suppliers of MRI-compatible patient monitoring systems. The new MAGLIFE Serenity guarantees highest ECG quality during magnetic resonance imaging (MRI) scanning even under strongest gradient influence. It monitors all vital parameters during anaesthesia in an MRI environment.









The Art of Diagnostics

RAD-BOOK 2010



RAD-BOOK 2010

Navigating on intraoperative images

The PoleStar® Surgical MRI System

agnetic Resonance images have proven to be indispensable for brain disease diagnosis since MRI was invented in the 1970's. Neurosurgeons use MR images to understand the location of pathological tissue in relation to healthy parenchyma and to functional areas, and to plan surgical access to lesions. With the introduction of intraoperative image-guided navigation systems in the early 90's, MR images registered to the patient enable more precise incision planning, minimally invasive surgeries, and precise access to small, deep-seated lesions. Today, most neurosurgical operating rooms are equipped with a navigation system and their use has become a de-facto standard of care.

However, intraoperatively, changes of the brain, which occur due to tissue resection and brain shift, make preoperative images less relevant. To provide up-todate information to the surgeon, the images have to be acquired during surgery.

Satisfying the Neurosurgeon's **Needs for Intraoperative Imaging**


PoleStar® Surgical MRI Systems deliver the power of MR imaging to the neu-

rosurgeon in real time during surgery. Surgeons can assess the degree of achievement of their surgical goals and the absence of any complications before they even close the patient. Whether the goal is complete tumor resection, or just a reduction in tumor size, MR images acquired during surgery provide the information to decide whether additional resection is necessary. Furthermore, the newly acquired images are displayed by the integrated

navigation system to provide accurate navigation throughout the surgery, even in presence of the inevitable brain shift. The navigation system can simultaneously display the preoperative diagnostic scans as well as any intraoperative scans so the surgeon can always review the progress of surgery.

PoleStar® users report that images acquired by the system enable them to achieve higher tumor resection rates. Senft et al¹ report:

"In our study, GTR was the preoperative surgical goal in 30 patients harboring contrast-enhancing tumors. In 8 (26.7%) of these patients, we were able to detect residual contrast-enhancing tumor by iMRI at a point when the surgeon believed GTR had been achieved. These intraoperative findings led to further tumor resection, resulting in GTR in all 30 patients (100%) as confirmed by postoperative 1.5-T MRI.



Intraoperatively acquired images were highly reliable. After intraoperative scanning, tumor resection was continued in an additional 4 patients not amenable to complete tumor removal, thus at least leaving less tumor tissue behind. The surgical goal of GTR or STR was achieved in all 42 patients with contrast-enhancing, histologically confirmed, high-grade gliomas." These results are similar to those reported for other intraoperative MRI systems.


Baumann et al² in their series of six giant pituitary adenomas resected with the benefit of PoleStar® imaging report "a much higher total resection rate of 66% (between 27.9% and 40% according

to literature), accompanied by at least as favorable endocrino-


logical outcome and symptom control as with conventional pituitary surgery" and conclude: "We believe these favorable outcomes are supported by the use of intraoperative imaging, which serves as an intraoperative quality control to resection and to preservation of surrounding structures beyond that of the operating microscope or



integrated with the StealthStation® S7™ System



Resection control in pituitary surgery. Left column – preoperative images, middle column – intraoperative images revealing residual tumor, right column – postoperative images demonstrating complete resection



Preoperative T1 weighted images with contrast agent

endoscope. This intraoperative information translates into more radical resections that may lead to fewer severe complications, fewer additional procedures (such as transcranial complementary operations, hormonal replacement therapy, radiotherapy), and to reduction of overall costs of therapy of this difficult group of patients."

The PoleStar® system enables surgeons to improve patient outcomes and positions the hospital as a regional leader in neurosurgery. This may attract more patients and leading surgeons to the hospital.

The PoleStar® Surgical MRI System

Starting with the introduction of the PoleStar® N-10 system in 1999, through the PoleStar® N20 announced in 2003, and on to the new PoleStar® N30 system, over 6000 cases have been performed at the more than 50 sites worldwide.

The PoleStar® system consists of a moveable MRI scanner and an integrated StealthStation® navigation system. Prior to surgery the scanner is moved to the OR table and positioned under the head of the table. With the push of a button the scanner's magnets move into position around the patient's head for imaging. Since the magnet moves to the patient, the patient can be positioned in a variety of standard surgical positions to minimize changes in surgical technique. Most of the hospital's standard instruments can also still be used when operating with the system.

The system's simple interface, which is typically operated by the surgical staff or surgeon, provides T1 and T2 weighted images, in addition to FLAIR and SSFP protocols. After the typical 10 to 15 minute imaging session the magnet is retracted to its home position under the OR table and

the surgery continues. Since the magnet is tracked by the navigation system, no patient registration is required and the acquired images are immediately ready for navigation

PoleStar® systems are installable in any existing operating room with only minimal room modifications. When the system is not in use, the scanner is stored in a special storage cabinet, making the room available for any other surgery there is no need to dedicate an OR to iM-RI use. RF shielding is achieved with either conventional whole-room shielding or with a novel StarShieldTM local RF shield. It is located, retracted, at the feet of the patient for most of the time, then is closed around the patient and magnet just during the scan sessions. This enables computers and other electronic equipment to be left turned on all the time, again minimizing changes to conventional surgical technique.

The New PoleStar® N30 Surgical MRI

The PoleStar® N30 is the latest addition to the PoleStar® system line. It provides better imaging while further reducing the magnet's size and weight and enhancing its integration into the OR. Due to a greatly improved gradient subsystem, the T2 and FLAIR images acquired by PoleStar® N30 are about 30% better than the images acquired by its predecessor systems. At the same time the magnet weighs 80kg less than

PoleStar® N20, providing easier handling of the scanner in the OR. The system's dimensions are smaller and the scanner can be moved laterally to simplify its positioning near the OR table.

The PoleStar® N30 system is released with a host of new features and accessories, including a new more ergonomic patient head holder. The PoleStar® N30 scanner can be used with hospital's existing Treon®plus or S7™ navigation system, eliminating the need to purchase additional navigation hardware.

Prof Moshe Hadani from Sheba Medical Center in Tel Aviv says:

"The Polestar N30 is the latest version of the system. It is smaller in weight and dimensions and is easily manipulated in the operating room.

We operated with this new generation system on 20 patients, for brain and pituitary tumor removal, for aspiration of brain abscesses and for accurate navigation to small lesions. We found that the quality of the images of T2 and FLAIR sequences improved dramatically, rendering them useful in the resection of low grade gliomas. When intraoperative MRI is used for the right indication, I can not overestimate its important contribution to a successful operation."

¹ Usefulness of Intraoperative Ultra Low-Field Magnetic Resonance Imaging in Glioma Surgery, Christian Senft, Volker Seifert, Elvis Hermann, Kea Franz and Thomas Gasser, Neurosurgery 63 (ONS Supplement 2), ONS 257–ONS267, 2008

² Intraoperative magnetic resonance imaging-guided transsphenoidal surgery for giant pituitary adenomas, Fabian Baumann ,Christoph Schmid, and René-Ludwig Bernays, Neurosurgical Review, 33:83-90, 2009



-Only 9 m² room size needed







▶ GE Healthcare Discovery MR750 3.0T & MR450 1.5T

MR450

50 mT/m 200 T/m/s 200 T/m/s

1.5T



Powerfully simple - Express preparation exam - »Can't miss« applications and HD coils

Highlights

Field

Gradient

Slewrate

MR750

50 mT/m

3.0T

- **Simply powerful** Shorter TE/TR & faster acquisitions with unique gradients architecture Faster reconstruction
- 27 % more SNR with optical RF technology



▶ Highlights

- 70 cm bore
- excellent images and a user-friendly experience access to full 50 cm FOV imaging and
- application performance
- Next generation of clinical applications for reduced exam time, and improved patient comfort
- Ability to handle more patients with fewer scans helps boost efficiency and makes scheduling more predictable.

▶GE Healthcare Optima™ MR360 1.5T Field 1.5T 33 mT/m Gradient Slewrate 100 T/m/s ▶ Highlights - Remarkable flexibility and efficiency to match a wide range of imaging needs high image quality and lower total cost of ownership

Technologists benefit from ease of use and confidence
 Radiologists benefit from expanded diagnostic capabilities

Administrators benefit from more satisfied patients, efficient throughput, and opportunities for growth





generation systems, GE »ecomagination« certified








- Broad choice of multi-channel phased array coils

Latest advanced spine applications










MR-Touch extends the principles of physical palpation with a precise and non-invasive visual evaluation of tissue stiffness of the liver

- Third party certified healthymagination product



owners having something different from others

▶ Hitachi Medical Systems Oasis



▶ Hitachi Medical Systems Aperto Eterna 0.4 T Field 22 mT/m Gradient 55 mT/m/ms Slewrate

▶ Highlights

- Wide, 320 degrees open permanent MRI system: the first class permanent full body MRI system allowing a different solution to HF applications
- It features the highest field strength 0.4T amongst the permanent MRI systems presently on the market.
- Newly developed built-in technologies keep Aperto Eterna delivering image quality comparable with entry level HF MRI
- Fast processing chain allows increasing patient throughput
- Reduced running costs allowing fast return of investment



- Comfort class permanent open MRI system, which keeps enhanced capabilities meeting sophisticated open design
- Allows newly developed technologies available at an excellent

Low cryogen boil-off technology

- Allows newly developed technologies a small constraint of the cost of ownership
 High magnetic field homogeneity
 Environment friendly: extremely low power consumption and reduced installation requirements
- Reduced running costs allowing fast return of investment

▶ Medtronic Polestar® Surgical MRI System

Magnetic Field Gradient Amplitude Gradient Slew Rate

23.5 mT/m 80 T/m/S Compact and mobile


0.15 T




Highlights

- Designed for integration in most OR's using mobile RF shielding
- Compatible with most existing surgical equipment.
- Perfect match to neurosurgical workflow.
- Fits under OR table and moves up to patient for imaging.
- Standard patient positioning, no patient movement needed during procedure.
- -Integrated StealthStation® Image Guided Surgery System to maintain navigational accuracy throughout the surgery.
- Auto-registration of images to patient
- Confirmation of completeness of resection and absence of complications prior to closing.






- The only one FDA-cleared, truly dedicated MRI system designed specifically for breast imaging
 Fully integrated, user friendly MRI guided interventional
- system
- Řevolutionary acquisition technology, SpiralRODEO™, provides far superior image quality Integrated Aurora-CAD™ software, the new standard in
- breast MRI software





- Easy to install, extremely reliable, quick financial break-even











▶ Philips Achieva 1.5T SE

Field Gradient Slewrate

1.5 T 33 mT/m 122 mT/m/ms



▶ Highlights

A true value-for-money 1.5T

system with comprehensive imaging capabilities
-Smarter economics with PowerSave (reduces energy bill by up to 50%) and compact siting (only 27m²)

- Built on proven Achieva platform offering wide choice of easy and economical upgrade paths

▶ Philips Intera 1.5T

1.5 T Field 33 mT/m Gradient 80 mT/m/ms Slewrate



Highlights

- SmartExam 1 click for consistent and reproducible MR scans
- NetForum community access with ExamCards for all studies
- Investment value FreeWave platform based ready for new applications
- SENSE for fast imaging in all applications

▶ Philips Panorama HFO

Field Gradient Slewrate

1.0 T 28 mT/m 120 mT/m/ms



▶ Highlights

- -Patient friendly: three times larger patient aperture than conventional MR to handle stressed and claustrophobic patients, children, elderly and large patients
- High-field performance comparable to 1.5 T in a truly open configuration
- Increased productivity with SmartExam, one-click planning, scanning and processing
- Enables unique applications not possible with cylindrical systems

▶ Siemens Magnetom Skyra, A Tim+Dot System

Field Gradient

3T 45 mT/m 200 T/m/s Slewrate



▶ Highlights

- Boost productivity with Tim[®] (Total imaging matrix) and Dot™ (Day optimizing throughput)

 - Ultra-light and short 3T system

 - Greater patient access and comfort with 70 cm Open Bore
- design
- TrueForm® Design for increased homogeneity
 DirectRF™ for higher signal purity and improved stability
 New fixed and Tim Dockable Table options

▶ Siemens Magnetom Verio, A Tim System

Field Gradient Slewrate

3T45~mT/m200 mT/m/s



Highlights

- Unique combination of 3T, 70 cm Open Bore and Tim (Total imaging matrix) Over 300 Verio systems installed
- Very flexible imaging from high throughput clinical to
- advanced research
 TrueForm Design offers enhanced image quality at 3T by optimizing the homogeneity

▶ Siemens Magnetom Trio, A Tim System

Field Gradient Slewrate

3T45 mT/m 200 T/m/s



▶ Highlights

- Excellent 3T magnet with unmatched homogeneity and strong gradients with AudioComfort
- Tim with up to 32 RF channels for outstanding image quality, speed and flexibility
- New trendsetting applications make the extraordinary routine

▶ Siemens Magnetom Aera, A Tim+Dot System

Field 1.5 T 33 mT/m Gradient 125 T/m/s Slewrate



Highlights

- Boost productivity with Tim (Total imaging matrix) and Dot (Day optimizing throughput)

Greater patient access and comfort with 70 cm Open Bore

- Ultra-light and short 1.5T system
- TrueForm Design optimizes the 50x50x45cm FoV

DirectRF - for higher signal purity and improved stability

New fixed and Tim Dockable Table options

▶ Siemens Magnetom Espree, A Tim System

1.5 T Field 33 mT/m Gradient 170 T/m/s Slewrate



Highlights

-First MR with 70 cm Open Bore proven with more than 900 installations

- Shortest system length of 125 cm only
- Especially appreciated by obese or claustrophobic patients
- Tim (Total imaging matrix) coils can be flexibly combined

-Field of view up to 205 cm with syngo TimCT

▶ Siemens Magnetom Avanto, A Tim System



- Leading applications with Tim (Total imaging matrix)
 8 out of the top ten U.S. hospitals work with the MAGNETOM Avanto
- -500 mm field of view, zero eddy-currents
- AudioComfort: ear protection not mandatory
- Parallel imaging from head to toe

▶ Siemens Magnetom Symphony, A Tim System



Highlights

Imginglits
Improved image quality with Tim
Complete range of applications powered by Tim
Advanced applications, like non-contrast enhanced applications syngo ASL (Arterial Spin Labeling) and syngo NATIVE
Whole-body capabilities with up to 200 cm Field of View

- Higher throughput with Tim's revolutionary coil concept

▶ Siemens Magnetom Essenza, A Tim System

Field 1.5 T Gradient 30 mT/m Slewrate 100 mT/m/s



▶ Highlights

- The most affordable*, all-new 1.5 T MRI low initial investment, low set-up and running cost, reduced siting costs
- Helps to enhance business opportunities by providing a broader range of applications and higher throughput with Tim

 – Enhance patient comfort with a ultra-short (145 cm)
- 1.5 T system and ultra-light weight coils
- Optimize image quality with the IsoCenter Matrix
- * Results may vary. Data on file

▶ Siemens Magnetom Espree - Pink

Field 1.5 T Gradient 33 mT/m Slewrate 170 T/m/s



▶ Highlights

- First dedicated 70 cm Open Bore breast scanner with shortest system length of only 125 cm
 Pink Comfort with Open Coil design and variable coil
- geometry (VCG) for both imaging and biopsy Sentinelle Vanguard for Siemens
- Pink Applications offering a complete portfolio for all needs
 Pink Workflow including dedicated workplace for reading and biopsy planning: syngo BreVis and syngo BreVis Biopsy

▶ Siemens Magnetom C!

Field Gradient Slewrate 0.35 T 24 mT/m 55 mT/m/ms



▶ Highlights

- Patient friendly exams due to side loading
- 270° accessibility assistance in reach of the patient
 Can be sited in less than 30 m²/323 sqft
 Offers a comprehensive set of clinical applications

- Excellent return-on investments: decreased costs - optimized profitability



- Connectivity of up to 128 coil elements with 16 or 32 channel-readout
- Image reconstruction rate of up to 4.000 images/sec
- -FBI MR angiography without contrast medium
- -55 x 55 x 205 cm scanning region







Image reconstruction rate of up to 4000 img/sec

optical data transfer



▶ Xingaoyi Ningbo (XGY) Oper - 0.5T

0,5 T Field 24 mT/m Gradient Slewrate 70 mT/m/ms




▶ Highlights

- The first mid-field Permanent Magnet MRI system used in clinical application in the world
- Richer performed scanning sequences, More complete function, Clearer images
- Much more patients with shorter scanning time

► Xingaoyi Ningbo (XGY) Oper - 0.4T

Field 0,4 T 20 mT/m Gradient Slewrate 66 mT/m/ms



Highlights

- Higher SNR & Larger Imaging range with Multi-RF Channels
 Excellent Image with fully performed scanning sequences
- Low power consumption, low failure rate, high operating ratio
- Small installation site

MRI ACCESSORIES

▶ Alliance Medical – flexible diagnostic imaging services



▶ Highlights

- Static diagnostic imaging centers MRI, CT, PET, PET/CT
- Interim services for bridging downtimes
- Regular "routing" services

▶ Xingaoyi Ningbo (XGY) Oper - 0.3T


Field 0,3 T 15 mT/m Gradient 48 mT/m/ms Slewrate



▶ Highlights

- Complete function, Excellent images, Fully performed scanning sequences
- Quick scan image is much clearer with high slew rate
- Extremely low power consumption and very low failure rate

▶ Invivo DynaCAD



▶ Highlights

- Comprehensive tool for automatic image processing

- Comprehensive tool for automatic image processing
 Improves workflow and diagnostic confidence
 Automatic calculation of subtractions, MPRs and MIPs
 Volume analysis feature for therapy response monitoring
- Fast, accurate interventional planning

▶ Invivo Precess





- Precess brings you the future of MRI patient monitoring Wireless ECG and SPO2
- Wireless remote display5.000 Gauss Compliant
- -3.0 T compatible



▶ Invivo Eloquence



▶ Highlights

- Eloquence is a complete, fully integrated fMRI solution
 Critical goals may be achieved with ease
 Integrated in-magnet compatible components

- Operator room control of both experiment and analysis
- High performance tools for postprocessing fMRI data

▶ Invivo Foot and Ankle Array Coil



▶ Highlights

- Produces exquisite MR images of the foot, ankle and toes
- »Ski Boot« design for fast and easy patient set-up Three 5° cradle tilt settings for patient comfort
- Patient comfort pads to reduce motion artifact
- Optimized for parallel imaging

▶ Schiller Maglife Serenity

MRI compatible up to 3 Tesla Mains and Battery driven (1,5 and 6 hours) 12,1" colour Display



- ▶ Highlights Optical core and skin temperature
- Configuration for Anaesthesia, car-Configuration for Anaestresia, cardiac und Intensive care applications
 Patented artefact inhibition
 6 optical Gating outputs
 Optimized for Adult children and

- neonates

▶ Schiller Maglife light

MRI compatible up to 3 Tesla Parameter: SPO2 and/or NIBP Mains and Battery driven (1,5 hours)



- HighlightsOptimized for day to day application
- No installation necessary
- HTML printing function
 Optimized for Adult children and neonates

▶ Tomovation - Modular building solutions



Engineering, rental, sale of modular buildings MRI, CT, PET, PET/CT including or excluding diagnostic equipment

▶ Ziosoft Ziostation


Full Suite of Clinical Applications

3D, 4D, MIP, MPR, CPR, Multi-Modality Viewer

▶ Highlights

- Multi-volume fusion
- Vessel analysis
- Coronary analysis
- Colon analysis
- PET/CT viewer
- CT Cardiac Function Analysis

- CT Cardiac Function Analysis
 CT Brain Perfusion
 Calcium Scoring
 Reporting functions
 Customizable button palettes and keyboard/mouse shortcuts to configure the system individually
- Automatic segmentation algorithms employ anatomical analysis, targeting specific objects to maximize accuracy
 Full Web-Based Capabilities
 Action macros perform multiple processes with one or zero
- button clicks



RAD-BOOK 2010

ACIST MEDICAL SYSTEMS	COVIDIE Surrously Types Head		
	Medicor	Performance. For life."	
MED (TRON AG			
	MNemo	to	
	TOMOVATION'		<i>lrich</i>

Acist EmpowerMR

Syringe Volume 1 – 100 ml in user-specified increments of 1 ml Pressure 40 – 300 psi in user-specified increments of 1 psi Flow rate 0.1 – 10 ml/sec in increments of 0.1 ml/sec

▶ Highlights

7 T tested, no minimum distance requirement from magnet

- Hydraulic technology – no motor, no battery in the MR

No interference with magnetic field, no magnetic attractive force

- Double-barrel injector, floor stand mounted

Ease-of-operation through intuitive operators interface distributed by TOMOVATION GmbH



Syringe Volume 1 - 200 ml in

user-specified increments of 1 ml

Pressure 40 - 300 psi in

user-specified increments of 1 psi

Flow rate 0.1 - 10 ml/sec in increments of 0.1 ml/sec



Highlights

Double-barrel injector, floor stand or ceiling mount

- Ease-of-operation through intuitive operators interface

- 10 ml/sec maximum flow rate with »change on the fly« control

-Extravasation detection (EDA) stops injecting if an extravasation is detected

 Networkable through IrisCT and CANopen interfaces distributed by TOMOVATION GmbH

Covidien CT 9000™ ADV CT injector

Capacity 30, 50 ml hand held; 50, 75, 100, 125 ml

high pressure prefilled; 200 ml empty

25 – 300 psi (1.7 – 20.7 bar) Pressure limit

0.1 - 9.9 ml/sFlow rate



▶ Highlights

Contrast delivery injector for computed tomography

New easy to read color touchscreen
Requires only 3 simple key strokes to activate injection

Accommodates front loading 200 ml disposable syringes and prefilled syringes

- Highly visible LED display

▶ Covidien Angiomat™ Illumena™ Angiography injector

Capacity 50, 75, 100, 125 ml high

pressure prefilled;

150 ml, 200 ml empty 75 – 1200 psi (5.17 – 82.74 bar) Pressure limit

in angio-cardiac and peripheral

modes, 75 – 300 psi (5.2 – 20.7 bar) in CT mode

Flow rate

0.1 - 40.0 ml/s angio-cardiac and peripheral modes; 0.1 - 10.0 ml/s CT mode

▶ Highlights

Injector for angiography, cardio-logy and CT contrast delivery

Digital powerhead display

- Automatically »flips« as powerhead is rotated

- Fill control bar allows easy, one finger operation

 Latex free and transparent syringes provide crystal clear view of the contrast medium

- Sensitive touchscreen display for direct and easy setup

▶ Covidien Optivantage™ DH Dual Head CT injector

Capacity 50, 75, 100, 125 ml high

pressure prefilled; 200 ml empty Pressure limit 50 – 325 psi (3.5 – 22.4 bar)

0.1 - 10.0 ml/sFlow rate



▶ Highlights

Contrast delivery injector for dual head injector protocols

Fully programmable powerhead: color coded display
 Patency check feature:

saline flush prior to injection

Timing bolus feature: to determine ideal scan

- Auto-fill feature: automatically filling of syringes

- Drip mode: changing drip rate, volume and duration

▶ Covidien Optistar™ Elite MRI injector

Capacity 10, 15, 20, 30 ml high pressure prefilled; 60 ml empty

Pressure limit 20 – 150 psi (1.4 – 10.3 bar)

for 60 ml syringes

Flow rate 0.1 - 0.8 ml/s

▶ Highlights

– Dual syringe injector for magnetic resonance tomography

Battery free operation

Single bolus and dual phase injections

Full color touchscreen for a greater visibility Switch from injection to drip mode at any time

- Drop in syringe loading reduces preparation time

► Medrad Intego™ PET Infusion System

Flowrate Radiation profile Max radioactivity

Max FDG volume Max FDG activity concentration Dose range FDG dose accuracy

1 mL/s

< 0.02 mSv/hr at 30.5 cm 27.7 GBq, shielding effectiveness

25.9 GBq, dose preparation 30 mL 1.48 GBq/mL

37 - 925 MBq +/- 2% of the

measured dose

▶ Highlights

- Automatically on-demand measurement of patient-specific FDG doses

- Safety features to reduce radiation exposure

of dose preparation > 20%

- Consistently and accurately administers FDG within +/-2% of the measured dose – providing precise, repeatable control of the FDG dose delivered.

► Medrad Dual syringe CT injector Stellant® DX

Capacity Pressure limit Flow range

A: 200 ml – B: 200 ml 325 psi (22.1 bar)

A: 0.1 - 10 ml/s in 0.1 ml/s increments B: 0.1 - 10 ml/s in 0.1 ml/s increments

Highlights

- Saline Flush Capability for contrast efficiency, improved image quality and reduced artifacts

Automated loading, filling and priming

Stores and recalls up to 32 protocols

- Precisely times contrast delivery with real-time display of injection pressure

MEDŘAD's P3T Personalized Patient Protocol software enables to program the contrast amount

and the delivery rate for each individual patient, enhancing image quality while maintaining efficient workflow.



► Medrad Spectris Solaris® EP MR Injection System

Capacity

65 ml contrast medium,

115 saline

Prefilled 5, 10, 15 20 ml

Pressure limit

325 psi (22.1 bar) 0.01 - 10mL/s in 0.01 mL/s increments Flow range

between .01 and 3.1mL/s; 0.1 mL/s increments between 3.1 - 10 mL/s



-Six user-programmable phases for added flexibility

 Optimized color touch screen support at-a-glance identification of protocol and quick check of fluid levels.

-Flexible power management allowing to switch to either battery or wall power in seconds

- 3T compatible



Medtron Accutron CT-D

200 ml (NaCl) Easy Loading Syringe (ELS)

200 ml (CM),

21 bar (304 psi) For both injection units:

Capacity Pressure limit

150 ml Selectable pressure

increasement

Flowrate

300 to 1200 psi / bar Variable 1 to 10 ml/s



Contrast and saline flush

cardiovascular power injector Precise fluid delivery, enhanced air management, fluid level sensing and gross air detection.

- Accurate injection pressure

control with user adjustable pressure limits

Bolus sharpness delivering exact variable and fixed contrast via a responsive syringe

Color graphical user interface with on screen tutorial for simplified setup



▶ Medrad Avidia

Capacity Pressure limit Flowrate 150 ml 1200/82 psi / bar 0.1 to 50 ml/s



▶ Highlights

 Fully automatic microprocessor controlled contrast medium injector for angiography applications with volume, flow and time control

-Power cable free, battery operated, fully mobile on a pedestal, with interface

Up to 20 single phase or multi phase injection programs possible
 Display of injector parameters on the injector head
 User control console with display and start button

▶ Highlights

Delivery Pressure

Capacity

Flow range

Absolutely wireless injector unit with rechargeable batteries

- Integrated heated syringe holder for Easy Loading Syringe (ELS)

- Wireless touchscreen remote control

-Use of prefilled syringes (as an option)

– Up to 6 phases


Secured injection position (built-in sensor)

Alternatively, display of injection parameters or pressure graph

Aluminium housing

Wall or ceiling syspension system (as an option)

CANopen Interface (as an option)



▶ Medtron Accutron MR

Capacity

65 ml or 200 ml (CM),

Delivery Pressure

65 ml or 200 ml (NaCl) Easy Loading Syringe (ELS)

Flow range

21 bar (304 psi) For both injection units: 0.1 – 10 ml/s, programmable in steps of 0.1 ml/s

▶ Highlights

- Absolutely wireless injector unit with rechargeable batteries

-Touchscreen control panel with different languages

- Wireless touchscreen remote control

Up to 6 phasesSecured injection position (built-in sensor)

– Use of prefilled syringes (as an option)

Alternatively, input of flow rate or phase duration

Pressure graph

Aluminium housing



Medtron injectron 82 HP

Capacity Delivery Pressure 200 ml (NaCl) Angio mode: 83 bar (1203 psi),

Flow range

CT mode: 21 bar (304 psi) Angio mode: 0.1 – 30 ml/s, CT mode: 0.1 – 10 ml/s

▶ Highlights

- Integrated heated syringe holder with Easy Loading Syringe (ELS) 200 ml - Fully digital, user programmable injector

- Remote control (as an option)

 Pressure jacket for prefilled syringes (as an option)

-Wall or ceiling suspension system (as an option) – Interface on request (as an option)

- 100 injection protocols can be defined and stored by the user (50 protocols Angio-Mode/ 50 protocols CT-Mode)

Aluminium housing



C O N T R A S T

Contrast agent injectors for:

- Computed Tomography
- Angiography
- Magnetic Resonance Imaging

Disposables

Meet us at Expo E, Booth 564



at the right time - at the right place



Hauptstr, 255 · D · 66128 Saarbrücken · Phone: +49 681 97017 · 0 · Fax: +49 681 97017 · 20 · info@medtron.com · www.medtron.com

▶ Medtron Accutron HP-D

Capacity

200 ml (CM), 200 ml (NaCl) Easy Loading Syringe (ELS)

Delivery Pressure Flow range

83 bar (1203 psi) 0.1 - 30 ml/s,programmable in steps of 0.1 ml/s

Highlights

 Absolutely wireless injector unit with rechargeable batteries

 Multiphase program controlled injection of CM and NaCl Single or multi injection mode

- Integrated heated syringe holder for Easy Loading Syringe (ELS)

-Touchscreen control panel with different languages

- Wireless touchscreen remote control

-Up to 3 phases

- Pressure graph

- Secured injection position (built-in sensor)
- 60 injection protocols can be definded and stored by the user

– Interface (as an option)

Aluminium housing

Medtron Accutron CT

Capacity **Delivery Pressure** Flow range 200 ml Easy Loading Syringe (ELS)

21 bar (304 psi)

0.1 – 10 ml/s, programmable in steps of 0.1 ml/s

▶ Highlights

- Absolutely wireless injector unit, rechargeable batteries

- Integrated heated syringe holder with Easy Loading Syringe (ELS) 200 ml
- Touchscreen control panel with

different languages

Wireless touchscreen remote control

- Secured injection position (built-in sensor)

Up to 6 phases

- Use of prefilled syringes (as an option)

 Alternatively, input of flow rate or phase duration

Display of injection parameter or pressure graph at the remote control

Interface capability (as an option)

Aluminium housing



Syringes

A: Contrastmedia A: 200 ml, 100 ml with adapter

B: Saline 50 ml

Pressure Throughput A: 300 psi, B: 300 psi A: 1-100 ml/200 ml in 1 ml steps



▶ Highlights

Needle positioning testProgammable autofill function

Program memory on CF

memory card

Throughput

- Creation of an optimized program by input of injection parameters

▶ Nemoto Dual Syringe-CT-injector DualShot Alpha

Syringes

A: Contrastmedia A: 200 ml, 100 ml with adapter/125

ml with prefilled syringe adapter

B: Saline 200 ml, 100 ml with adapter A: 300 psi, B: 300 psi

Pressure Throughput

A: 1 - 100/125/200 ml in 1 ml steps

B: 1 - 100/200 ml in 1 ml steps



Highlights

Needle positioning test

- Programmable autofill function - Program memory on CF

memory card


Advanced Programming Functions

- Timing Bolus option

– Auto prime function

Nemoto CT-injector A 60

200 ml, 100 ml with adapter Syringes Pressure 300 psi



Pressure Throughput

▶ Nemoto MR-injector Sonic Shot GX

A: MR Contrastmedia A: 60 ml standard; 10, 15, 20 Syringes

ml with adapter (for prefilled syringe)

B: Saline

A: 200 psi, B: 200 psi

A: 0.1 – 10 ml/s in steps of 0.1 ml/s B: 0.1 - 10 ml/s in steps of 0.1 ml/s



Highlights

- LCD-display

- Real time monitoring of the injection parameters

Economical entrance model

Highlights

Needle positioning test

- Progammable autofill function

Injector with MRI-compatible ceiling suspension available
 To be used for prefilled syringe

Nemoto Angio-injector RemPress

Syringes Pressure

150 ml 50-1200 psi

Throughput 0.1-25 ml/s in 0.1-ml/s steps



ulrich medical CT/MRI injector mississippi (XD 2000)

CA max. 2 x 1.000 ml (for CT), Media containers

CA max. 2 x 100 ml (for MRI) NaCl max. 1 x 2.000 ml max. 400 ml/patient

Injection volume 16 bar Pressure Flow rate

0.2 – 8.0 ml/s, by increments

of 0.1 ml/s



- Roll pump injector for CT and MRI

 Several injections consecutively out of one media container (multi dosing)

- Battery operated

- Proven hygienic safety

- Different software options available



• ulrich medical CT injector missouri (XD 2001)

Media containers

▶ Highlights

- Test Shot Mode

- Easy handling

- Infusion Mode

- Various installation styles

-21 cm touch screen monitor

CA max. 2 x 1.000 ml, NaCl max. 1 x 2.000 ml max. 400 ml/patient

Injection volume Pressure Flow rate

16 bar 0.2 - 8.0 ml/s,

by increments of 0.1 ml/s



▶ Highlights

- Roll pump injector for CT
- Several injections consecutively out of one media container (multi dosing)

Economic consumption of disposables

- Proven hygienic safety

- Different software options available

ulrich medical CT injector ohio tandem (XD 2002)

CA max. 2 x 1.000 ml,


Media containers

NaCl max. 1 x 2.000 ml max. 400 ml/patient

Injection volume Pressure Flow rate

16 bar 0.2 - 8.0 ml/s,

by increments of 0.1 ml/s



▶ Highlights

- Roll pump injector for CT

- Several injections consecutively out of one media container (multi dosing)

-Tandem function for different contrast agents without previous change of media containers

- Proven hygienic safety

- Different software options available

▶ ulrich medical CT/MRI injector ohio M with tandem function

Media containers

CA max. 2 x 1.000 ml (for CT), CA max. 2 x 100 ml (for MRI) NaCl max. 1 x 2.000 ml

Injection volume Pressure

max. 400 ml/patient 16 har

Flow rate 0.2 - 8.0 ml/s, by increments of 0.1 ml/s



▶ Highlights

Roll pump injector for CT and MRI

- Several injections consecutively out of one media container (multi dosing)

Battery operated

- Tandem function for different contrast agents without previous change of media containers

Proven hygienic safety

ulrich medical MRI injector tennessee (XD 2003)

CA max. 2 x 1.000 ml (for CT). Media containers

CA max. 2 x 100 ml (for MRI)

NaCl max. 1 x 2.000 ml Injection volume max. 400 ml/patient

Pressure 16 bar

Flow rate

0.2 - 8.0 ml/s, by increments

of 0.1 ml/s



- Roll pump injector for MRI accumulator free

- Several injections consecutively out of one media container (multi dosing)

Ready for use anytime

- Smooth workflow without interruption of daily workflow

– Proven hygienic safety



RAD-BOOK 2010









Design Detector Resolution Size

Biplane Angio system


1.5k a-Si 2.7 lp/mm

31 x 31 cm frontal, 31 x 31 cm lateral



Highlights

- Optimal detector size for dedicated neuro applications
 Highest DQE for significant dose savings
- FP-CT techniques
- Powerful 3D-processing tools
- Total in-room control



▶ GE Healthcare Innova 2121 IQ

Design Detector Resolution

2121 IQ

Biplane Cardiac system

1k a-Si Size 2.7 lp/mm

20 x 20 cm frontal, 20 x 20 cm lateral



▶ Highlights

- Image quality goldstandard in cardiac imaging
 Smart gantry for optimal C-arm positioning

- Highest DQE for significant dose savings
 Complete integration of Intra-Vascular-Ultrasound
- Seamless data exchange for comprehensive workflow solution

▶ Philips Allura Xper FD 20/10

Design Biplane Detector a-Si Resolution 0.154 mm Size 38 x 30 cm

Total in-room control



- Live 3D guidance with XperGuide for needle planning and guidance and Dynamic 3D Roadmap for easy navigation
 Isoft tissue information available with XperCT, controllable at
- Xper Workspace integration for instant recall and viewing of all related multimodality images – from CT, MR and precious x-ray cases – at tableside during the intervention
 Flat detector with 2k digital imaging chain featuring; complete
- 2048 x 2048 pixel, digital imaging chain

▶ Philips Allura Xper FD 20

Design Detector Resolution Size

Ceiling-mounted a-Si 0.154 mm 38 x 30 cm



▶ Highlights

- Live 3D needle guidance, bringing back needle guided interventions into the angio suite with XperGuide

Complete 2048 x 2048 pixel, digital imaging chain

- Image area of 30 x 40 cm adjustable to a square image of 16 cm - Powerful set of diagnostic tools, e.g. Bolus Chase, Rotational Scan

Accessibility to innovations such as high-speed XperCT and 3D Roadmapping

Shimadzu Bransist safire VB/HB

Detector Resolution Size Direct-conversion flat-panel detector (a-Se)

3.3 lp/mm

9" x 9" (22 x 22 cm)



Highlights

Bi-plane Angio/Cardio systemFlat-panel detector technology for

outstanding image quality
Grid controlled X-ray tube for superior dose management

Exclusive »Cyber Chase« feature to keep the ROI automatically in bi-plane rotation

Excellent coverage without moving a patient

- Motion-artefact-free by patented mask-less DSA technology

▶ Shimadzu Bransist safire VC/HC

Detector Resolution Size

Direct-conversion flat-panel detector (a-Se)

3.3 lp/mm

9" x 9" (22 x 22 cm) or 17" x 17" (43 x 43 cm)



▶ Highlights

- Ceiling-mounted C-arm

- Flat-panel detector for outstanding image quality

Grid-controlled X-ray tube for superior dose management

- High speed C-arm up to 60°/sec. rotational DSA

- Excellent coverage without patient moving

▶ Shimadzu Bransist safire VF/HF

Detector Resolution Direct-conversion flat-panel detector (a-Se)

3.3 lp/mm

9" x 9" (22 x 22 cm) or 17" x 17" (43 x 43 cm) Size



Highlights

- Floor-mounted C-arm

- Flat-panel detector for outstanding image quality

- Grid controlled X-ray tube for superior dose management

- High speed C-arm up to 60°/sec. rotational DSA

- Excellent coverage without patient moving

Flexibility and reliability with triple pivots

▶ Siemens Artis zee floor-mounted

Design Detector Resolution Size

Universal floor-mounted flat detector angiography system

2k a-Si with CsI scintillator 1920 x 2480 pixel, 3.25 lp/mm

30 x 40



▶ Highlights

 Advanced 3D imaging at low dose
 Slim-line design and flexible positioning capabilities for easy patient access with full body coverage

New ergonomic system controls for smooth table-side operation - Complete 3D portfolio including cross-sectional imaging with syngo DynaCT

▶ Siemens Artis zee ceiling-mounted

Design Detector Resolution Size

Universal ceiling-mounted flat detector angiography system

2k a-Si with CsI scintillator 1920 x 2480 pixel, 3.25 lp/mm

30 x 40



▶ Highlights

 Advanced 3D imaging at low dose
 Slim-line design and flexible positioning capabilities for easy patient access with full body coverage

New ergonomic system controls for smooth table-side operation Complete 3D portfolio including cross-sectional imaging with syngô DynaCT

IS CARDIOVASCULAR

Siemens Artis zee biplane

Design Detector Resolution Size

Biplane flat detector angiography system

2k a-Si with CsI scintillator 1920 x 2480 pixel, 3.25 lp/mm

30 x 40 per plane



▶ Highlights

- Advanced 3D imaging at low dose
 Slim-line design and flexible positioning capabilities for easy patient access with full body coverage
- Largest biplane anatômical coverage available today Clinical flexibility – from neurovascular to spine and abdominal imaging
- New ergonomic system controls for smooth table-side operation
- Complete 3D portfolio including cross-sectional imaging with syngo DynaCT

▶ Siemens Artis zeego

Design Detector Resolution Size

Multi-axis flat detector angiography system

2k a-Si with CsI scintillator 1920 x 2480 pixel, 3.25 lp/mm



▶ Highlights

- -Cross-sectional imaging with Large Volume syngo DynaCT to visualize the whole liver or the whole lumbar spine
- Complete 3D portfolio Small footprint and multiple park positions
- Ideally suited for the OR environment
- Flexible working height reduces fatigue associated with long
- Pre- and post-operative high-end imaging directly in the OR

▶ Toshiba Infinix - CFi/BP

Biplane C-Arm + Omega-Arm Design Detector 20 x 20 cm; 30 x 30 cm Tube 3 MHU, 200 mA pulsed



▶ Highlights

- Single User OperationSequential Navigation for fast and easy throughput
- Hyper Handle for One Hand Operation

- FollowME Concept
 Full 3D capacity for Angiography
 Full range of Dose optimization techniques

SAFE CAREFUL Automatic insufflation of CO, for virtual colonoscopy Easy adjustment of gas volume and pressure Premium diagnostic results ulrich GmbH & Co. KG Buchbrunnenweg 12 89081 Ulm | Germany Phone: +49 (0)731 9654-0 Fax: +49 (0)731 9654-2805 insufflator@ulnchmedical.com www.ulnchmedical.com



Design Mono C-Arm floor mounted Detector 20 x 20 cm or 30 x 30 cm Tube 3 MHU, 200 mA pulsed



Highlights

- Five Axis System for maximum freedom and flexibility
- Single User OperationSequential Navigation for fast and easy throughput
- Hyper Handle for One Hand Operation
 FollowME Concept
 Full 3D capacity for Angiography

- Full range of Dose optimization techniques

▶ Toshiba Infinix - VCi/BP

Design Biplane C-Arm + Omega-Arm Detector 30 x 40 cm; 30 x 30 cm Tube 3 MHU, 200 mA pulsed



- FollowME Concept Low Contrast Imaging
- 3D rotational angiography
- Full range of Dose optimization techniques

▶ Toshiba Infinix - VFi/SP

Design Mono C-Arm floor mounted Detector 30 x 40 cm or 30 x 30 cm 3 MHU, 200 mA pulsed Tube



▶ Highlights

- Single User Operation
- Five Axis System for maximum freedom and flexibility
- Sequential Navigation for fast and easy throughput
 Hyper Handle for One Hand Operation
 FollowME Concept

- Low Contrast Imaging
- 3D rotational angiography
 Full range of Dose optimization techniques

▶ Toshiba Infinix – CCi

Design Mono C-Arm ceiling mounted Detector 20x20 cm or 30 x 30 cm Tube 3 MHU, 200 mA pulsed



- ► Highlights

 Single User Operation

 Sequential Navigation for fast and easy throughput

 Hyper Handle for One Hand Operation

- FollowME Concept Full 3D capacity for Angiography
- -Full range of Dose optimization techniques

▶ Toshiba Infinix - VCi

Mono C-Arm ceiling mounted Design Detector 30 x 40 cm Tube 3 MHU, 200 mA pulsed



▶ Highlights

- Single User Operation
 Sequential Navigation for fast and easy throughput
 Hyper Handle for One Hand Operation
- FollowME Concept
- Low Contrast Imaging 3D rotational angiography
- Full range of Dose optimization techniques

IS MOBILE C-ARM

▶ GE Healthcare OEC 9900 Elite 15 kW Power II-Format 11, 15, 23 and 31 cm ▶ Highlights - DRM (Dynamic Range Management) - Intuitive touchscreen interface - Comfortable viewing with flat screen monitors - Easy archiving: CD/DVD and DICOM - Fully motorized imaging system



- Innovative connectivity solution: CD/DVD, USB and DICOM - Available as Compact, Compact², Compact+ and Series












▶ Philips BV Pulsera with 3D-RX

15 kW Power II-Format 31/23/17 cm



▶ Highlights

- Highlights
 Based on the BV pulsera
 Unique combination of conventional 2D c-arm flexibility and top-quality 3D imaging in a single compact system
 12" image intensifier for largest 3D reconstruction, expanding clinical applications
 Motorized propeller movement of 200 degrees in only
- 30 seconds acquisition




- widest range of applications, including cardiac interventions
 SmartVision a fully digital imaging chain including powerful
- image processing functions

 High quality images at a low dose, time after time

 Pulsed acquisition 30 pulses/sec


- Rotating anode power











- Optional DICOM compatibility


- NaviLink interface for 3D navigation




-Full DICOM compatibility





































Data evaluation by means of the DiaControl expert software



▶ Radcal ACCU-PRO™

X-Ray Analyzer Simultaneous dose, rate, time, kVp. HVL, filtration, mA/mAs, and more



▶ Highlights

- Use for manufacturing, installation, QA, and service R/F, mammography, CT, dental, leakage Ion chamber based dosimetry, no corrections required

- Correctly measure AEC fluoro and filtered beams
- Remote control, waveforms, and archiving with XLPRO software
- Compact, easy to use

▶ Radcal PDC-DAP/KAP verification meter



Highlights

- Newly patented Patient Dose Calibrator
- -Use to calibrate DAP/KAP and rate
- Also measures dose and rate
- Optical and radiographic alignment markers
- Simple to use with optional computer control

▶ RTI Electronics Piranha

The Piranha is designed as a truly self-contained, all-in-one, X-ray multi-function meter that assures accurate results in one shot. kV, time,


▶ RTI Electronics Barracuada

The Barracuda X-ray multimeter has a cabinet that can house up to six different application modules, and can measure on all modalities; R/F, mammography, flouroscopy, pulsed flouroscopy, dental, panoramic dental and CT systems



- ► Highlights

 All in One, All at Once
- Auto-Compensation
- Enhanced Graphical PDA Display R&F, Mammo, Dental and CT
- Ionization Chambers
- Built-In Bluetooth for PC and PDA
- -mAs, and Light Probes
- Fits in the Palm of Your Hand





RAD-BOOK 2010



Large Teaching Hospital Evaluates 3D Imaging for Breast Cancer Screening

King's College, London, UK

ing's College Hospital in London, UK has begun a trial program utilizing breast tomosynthesis, a new 3D digital x-ray technology for breast cancer screening and diagnosis. The first national hospital in the UK to undertake this kind of trial, King's College anticipates that if the trial is successful, 3D breast imaging could offer new hope in the fight against breast cancer. King's College is one of the UK's largest and busiest teaching hospitals, with over 6,000 staff assisting approximately 700,000 patients a year.

Breast cancer is the leading cancer in the UK with 46,000 new cases diagnosed each year and over 12,000 dying from the disease. The hospital hopes the trial will prove that 3D technology can assist doctors in reducing the estimated 70,000 to 80,000 women every year who are mistakenly told that something unusual has been found, as well as the small number of women mistakenly given the all clear. The system used in the trial is a Hologic Selenia® Dimensions® breast tomosynthesis system, now commercially available in Europe and other areas of the world and in the FDA process in the U.S.



King's College Hospital hosted the first ever hands-on breast tomosynthesis users training meeting, attracting radiologists from tomosynthesis sites throughout Europe.

2D Imaging Shows Room for Improvement

At present, 2D mammography is the standard breast x-ray used in the UK. Although recognized as safe and reliable in detecting the early signs of cancer, the "anatomical noise" associated with 2D imaging can sometimes hide cancers, or produce shadows that falsely create the suspicion of cancer. Some cancers remain undetectable with 2D technology. The trial hopes to prove the efficacy of 3D in overcoming this deficiency.

Digital breast tomosynthesis is a 3D imaging technology that enables doctors to look at separate "slices" of the breast. The system acquires individual images of a stationary compressed breast at multiple angles during a short scan. Individual images are reconstructed into a series of thin high-resolution slices that can be displayed individually or in a dynamic ciné mode

Reconstructed tomosynthesis slices reduce or eliminate the problems caused by tissue overlap and structure noise in single slice 2D mammography imaging. Digital breast tomosynthesis also offers a number of exciting opportunities, including improved diagnostic and screening accuracy, fewer recalls, greater radiologist confidence and 3D lesion localization.

"This is exciting new technology that could improve the accuracy of breast cancer screening and help save more lives by detecting more breast cancers when they are small and at an early stage, when they can be more effectively treated," states Dr. Michael Michell, director of

breast screening at King's College. "In human terms, this technology could spare a lot of women a lot of heartache, and also save the NHS (National Health Service) valuable resources through people not having to attend follow-up appointments and undergo further tests."

Sarah Sellars, assistant director of NHS Cancer Screening Programs, adds, "The NHS Breast Screening Program welcomes research into the use of new technology for screening, and we await the outcome of this study with interest."

The trial use of digital breast tomosynthesis at King's College has been authorized by the hospital's Ethics and Research and Development committees. The clinical results of the trial at King's will be considered by the NHS Breast Screening Program.

Update

In November 2009, Dr. Michell reported on the first results of the trial - a reader study involving women that attended an assessment clinic for evaluating a mammographic abnormality found in a routine film-screen exam.* Ninetyone percent of the women that were eligible for the study agreed to participate in a follow-up bi-lateral 2D digital mammography exam and a breast tomosynthesis exam.

The radiologists reading the results of the first 300 cases found that their diagnostic confidence increased with breast tomosynthesis. These results, Dr. Michell said, supported the need to move forward with a larger, multicentre study.

^{*} M.J. Michell et. al. "Digital Breast Tomosynthesis: A Comparison of the Accuracy of Digital Breast Tomosynthesis, Two-Dimensional Digital Mammography and Two-Dimensional Screening Mammography (film-screen)", Royal College of Radiologists Breast Group Annual Scientific Meeting, Belfast, Northern Ireland, November 1-3, 2009.

Professor Patrice Taourel of the University Hospital Lapeyronie of Montpellier on

Breast Tomosynthesis

ne of the largest teaching hospitals in France, the University Hospital Lapeyronie of Montpellier prides itself on utilizing state of- the-art technology to provide the highest levels of patient care, teaching, and research. Six years ago, the hospital replaced its analog mammography system with a Hologic Selenia digital mammography.

The difference in the quality of images between analog and digital mammography is dramatic, and digital mammography has proven to be a superior

modality for detecting breast cancer. But 20 – 30 percent of cancers remain undetected even with digital mammography, and too many women receive a call back for a second exam. So, when Hologic's tomosynthesis three-dimensional technology became available, the University Hospital became one of the first sites in France to implement the technology.

2D and 3D Technology: Better Views of the Breast

Tomosynthesis technology provides major improvements over digital mammography. Breast tomosynthesis reduces the interference of overlapping tissue by providing three-dimensional views of the breast and making it easier for doctors to see cancer especially in higher-risk women with dense breasts.

"Tomosynthesis supplements mammography," states Professor Patrice Taourel, head of the Radiology Department. "We use it on a systematic basis for patients who need more detailed views. From the beginning, we established a protocol to use tomosynthesis for all women who fit the stereotype of:

- It is their first mammogram and they do not have comparison films;
- They have risk factors, such as family history;

Well known for his work in breast imaging MRI, Professor Taourel thinks MRI is not necessary or practical as a breast imaging tool for all women. On the other hand, tomosynthesis can be used easily for patients. It takes less time and can be used at the first point of contact with a patient.





- They have a history of breast cancer and are looking for other lesions, or
- They want a second opinion.

"Our workflow has not changed with tomosynthesis," continues Professor Taourel. "It only requires several more seconds. For the patient, the change is almost transparent. For our radiologists, they need only to familiarize themselves with the new images they can see with tomosynthesis. It was the same when we went from analog to digital technology. The images are different, you see a lot more, and radiologists need to be properly trained to become familiar with these new images."

Tomosynthesis and MRI as Breast Screening Tools

Well known for his work with breast imaging MRI, Professor Taourel thinks MRI is still the best imaging modality, but it is not necessary or practical as a breast

imaging tool for all women.

"Tomosynthesis can be used easily for patients. It takes less time and can be used at the first point of contact with a patient. It is very helpful for breast health."

Professor Taourel feels more confident in his diagnosis that the lesion is either suspicious or non-suspicious with breast tomosynthesis. And tomosyn-thesis reduces the number of women the University Hospital of Montpellier recalls. His dream is to use tomosynthesis to guide biopsies.

Looking much like a digital mammography system, the x-ray tube of the Selenia Dimensions digital breast tomosynthesis system rotates around the breast taking 15 images from numerous angles in an exam that takes less than 4 seconds.

RAD-BOOK 2010

Pioneers by Tradition

istory does not reveal the number of millisieverts used by the X-ray machine presented to Dr. Karl Rohrbacher of Vienna at the beginning of the 1940s. However, this machine laid the foundation for one of the most modern radiology practices in Austria after the turn of the century, the Diagnosezentrum Urania (DZU). And today, the center uses the mammography system with the lowest radiation dose available, Sectra MicroDose Mammography.



Dr. Wolfgang Drahanowsky, Technical Director at the Diagnose Zentrum

Looking at the history of the practice, it comes as no surprise that the installation of the Sectra system in 2004 marked Sectra MicroDose Mammography's arrival in the Austrian market. Ever since Dr. Rohrbacher became an established radiologist, the philosophy has been to make use of technological progress. This approach has intensified since Prof. Dr. Peter Riedl took over the practice in 1979 and developed it into today's center with the help of his partners. For example, Riedl was the first private CT operator in Vienna and the DZU began to operate Austria's first 3-Tesla MR as early as in 2004

When Swedish developers from Sectra created a new, low-dose form of mammography using photon-counting technology, it was only a matter of time before it also became available in Vienna. "CT be-

From examination to diagnostic findings in only 30 minutes.

came routine; 3-Tesla became routine. With Sectra MicroDose, we were once again at the forefront of future-proof technology, which we are sure will become even more widespread in the future," comments Dr. Wolfgang Drahanowsky, Radiologist and Technical director at the Diagnosezentrum Urania (DZU) and one of the leading figures at the facility, alongsideUniv. Doz. Dr. Peter Barton , Univ. Doz. Dr. Lucas Prayer and Prof. Dr. Peter Riedl .

Because less is more

In the past three years, Sectra Micro-Dose has won over the team, primarily because of its extremely high image quality and the distinct 50% reduction in radiation compared with other digital systems. "It's the best dose reduction available in the medical technical equipment market. In 2009, we participated in a comparative study conducted by medical physicists into the radiation impact of various mammography machines and this provided confirmation," says Dr. Drahanowsky, referring to the currently unrivaled low radiation doses with no loss of quality using the photon-counting technology devel-

oped by Sectra. In this, the special multislit detector counts all of the x-rays directly, without previously converting the radiation to light. The effect of the detector design is that scatter is reduced by 97% instead of the 70% reduction common with flat panel digital systems. At the same time, with no conversion steps, there is zero electronic noise.

Equipping for the national screening program

In screening programs, minimizing radiation is the top priority. In Austria, the process of so-called "gray" mammography is still practiced: Every woman can have a mammography at the institution of her choice on a physician's referral. However, this year, Austria will follow the current trend in adopting a national screening program for the early detection of breast cancer and the DZU, with its two Sectra installations regards itself as fully prepared. "The referring physicians particularly appreciate the low dose we offer their patients. We receive direct referrals as a result of our MicroDose equipment," says Dr. Drahanowsky.

"Since there is emphasis on close cooperation with local hospitals, we also always try to maintain the top standard at clinical level in the form of the latest technology." Accordingly, the mammography unit at the DZU has not only received the "Mammadiagnostic Quality Certificate" from the Austrian Radiology Society and





Routine in the Diagnose Zentrum Urania: Every picture is diagnosed by at least two radiologists

the Federal Group of Radiology Specialists of the Austrian Medical Chamber, but has also achieved excellent reliability and performance reviews: "We conduct up to 100 mammography examinations per day.

This means that we cannot accept any

From image to result in 30 minutes

technical delays."

At the DZU, the team provides continuous double reading by two specialists. The diagnoses are made immediately after the tests and directly on-screen by the examining physician and an independent

second opinion is also provided. Both physicians are aided in their work by a CAD System (Computer Aided Detection) that is integrated with the Sectra PACS. The CAD system acts as a third reader, further ensuring that nothing is missed. "With Sectra PACS as a tool, the workflow for double diagnosis works beautifully." says Dr. Drahanowsky. "This means that the first and second experts can examine the data independently and make a simultaneous independent evaluation of the image material." After only 30 minutes, the patient receives an envelope with her results. It couldn't be faster!

The Sectra PACS is fully integrated with the DZU's overall

IT system. Images and data from various imaging technologies or manufacturers, such as from the mammography unit's long-term archiving systems, can be accessed at any time without any problems. This means that comparative images of patients, who have already been examined at the DZU, can again be used in evaluation. Even the smallest changes in breast tissue over time can be found in this way

One of the most modern institutions of its kind in Vienna: The Diagnose Zentrum Urania



In the meantime, Dr. Drahanowsky's assumption regarding the establishment of Sectra technology has become reality: other radiological centers in Austria are already catching on and are also investing in MicroDose Mammography. The progressive approach of the institute has once again proven correct.

Benefits of Sectra MicroDose Mammography[™] at a glance

- Excellent high resolution, high-contrast images
- Reduction of at least 50% in the radiation dose compared to other available systems
- Efficient workflow and higher patient processing, with up to 18 exams per hour
- Excellent ergonomics
- A curved and heated patient support for increased comfort and ease of positioning

MicroDose conquers the Vienna market

The "Regional Structural Plan for Health in Vienna" envisages a dramatic reduction from 54 radiology practices and 14 diagnostic centers to an anticipated 16 diagnostic centers! Each of these 16 diagnostic centers would provide for an area

with about 70,000 women. To continue to offer the patients effective mammography screening, reliable and efficient equipment with a high rate of processing must be acquired. As a result, four radiology practices in the Austrian capital rapidly decided in favor of a Sectra MicroDose system. In addition to the DZU, which decided to acquire a second MicroDose, the Diagnosticum Gersthof, the Hietzing Diagnostic Center and Röntgen Liesing also opted for the Sectra mammography equipment. Accordingly, Sectra with its partner, Enomed, has managed to increase its market share for digital mammography machines to 25% in a very short time.

Mammography at its Best

Seamless Image Logistics, Optimized Reporting, and Structured Archiving in Mammography

n western countries, the carcinoma of the breast is the most common cancer in women. In order to improve the cure rate, it is of major importance to detect the disease at the earliest possible stage. As a result goverments have initiated nation-wide mammography screening programs. For this purpose, it is necessary to link a variety of IT systems, and to ensure smooth workflow across a multitude of locations. But curative mammography, too, entails particular requirements regarding image logistics and image display.

Optimized Reporting Based on Digital Mammography

Mammographies are captured, to an increasing extent, digitally through Computed Radiography (CR) or Direct Radiography (DR) systems. The complementary exams relevant for mammography are usually also carried out in digital format (MRTs or ultrasound exams). In combination with information systems linked via worklists, this provides the cornerstone for the realization of a smooth holistic digital workflow. Patient admission is done in the information system, as is the planning of, and request for, the exam. The worklists generated automatically based on this information ensure the necessary consistency of the data. This safeguards that it is possible to locate all the data pertaining to a case. This applies to the technologist in radiology as well as to the reporting radiologist.

During reporting, routine work steps are automated by highly specialized hanging protocols. These have been programmed on the basis of know-how from many expert reporting physicians. In order to make sure that the physician "feels

at ease" with his reporting environment, the hanging protocols as well as their graphic user interface can be adapted to the individual needs.

The main support for the reporting process is provided by automated prefetching of previous exams and their standardized display on screen. A multitude of image processing are there to provide additional support. In conclusion, this leads to increased certitude in reporting and therefore to optimized quality in reporting. Further support can be achieved by integration of computer-aided detection systems (CAD). For additional benefits regarding ergonomics, the application can be managed via keyboard as well as, in addition, a keypad.

Reporting documentation is realized in a structured manner and in accordance with international standards (the ACR scheme and BI-RADS classification). The categorization defines, on the one hand, mamma density (ACR) and, on the other hand, the risk for carcinomas is appraised. Subsequent care strategies for the patient can be deduced from this.

Flexible Archiving, Rapid Access

Archiving of images and reports is done in a flexible manner, adapted completely to the needs of the individual customer. For example, data can be stored on CD or DVD in a single media archive, or in networked structures in the Storage Area Network (SAN). The most recent option is Storage Service Providing (SSP). This service enables the customer to store his data in an external archive operated by a service provider. This type of data storage enables cost-optimized operation for the user. With full orientation of fees towards use, the customer is not required to make investments and tie up capital. Cost can be planned because there are no additional expenditures for migration, administration, or maintenance.


No matter which archiving technology the user selects, his access to all previous exams is ensured at any time. This includes information from a third-party PACS, too. In addition, X-ray films generated and scanned externally can be integrated seamlessly into the workflow.

Increased Requirements Regarding Image Logistics in Screening

Mammography screening entails extended requirements regarding integration and image logistics. Software applications for screening campaign purposes which organize the invitation procedure and manage the client data.

www.visus.com





▶ Fujifilm Amulett

Technology Resolution Size

Direct Optical Switching, a-Se

50 μm, pixel size 18 x 24 cm / 24 x 30 cm



▶ Highlights

- 50 micron pixel at extremely low radiation dose
- High DQE and high MTF
 Direct Optical Switching technology
- Excellent ergonomic design
 Optimized workflow and quick, comfortable examination
- -Stereo-Option available

▶ GE Healthcare Senographe Essential

Technology Resolution Size

a-Silizium 100 µm 24 x 31 cm



▶ Highlights

- High patient throughput
- Dual track tube Mo/Rh Automatic Optimization of Parameters (AOP)
- -Ergonomic paddles that shape to the breast
- Stereo-Option available

▶ GE Healthcare Senographe DS

Technology Resolution Size

a-Silizium 100 µm 19 x 23 cm



▶ Highlights

- High patient throughputDual track tube Mo/Rh
- Automatic Optimization of Parameters (AOP)
- Ergonomic paddles that shape to the breast
- Stereo-Option available

▶ GE Healthcare Senographe 2000D

Technology Resolution Size a-Silizium 100 µm 19 x 23 cm



▶ Highlights

- Automatic and fast image processing
- Dual track tube Mo/Rh
- Automatic Optimization of Parameters (AOP)
- DICOM connectivity for digital workflow
- Proven technology

▶ GE Healthcare Performa

Anode Filter kV Range

▶ Highlights

- Compact, small system

Мо Mo/Rh 20 - 35



▶ Highlights

Combines traditional mammography and multi slice 3D tomosynthesis imaging capabilities

 Reconstructed tomosynthesis slices reduce or eliminate

the problems caused by tissue overlap and structure noise in mammography imaging offering improved diagnostic confidence and enhanced patient care

Selenia Dimensions breast tomosynthesis is cleared for sale in the European Community and awaiting FDA clearance in the U.S.





RAD·BOOK 2010

▶ Highlights

▶ Hologic Selenia® digital mammography



- -Flexible and interactive tools for the mammographer to use to increase productivity and optimize workflow
- Seamless integration with hospital infrastructure, assuring the best possible care
- -Smart påddle system streamlines technologist workflow

▶ Hologic MammoSite® multi-lumen targeted radiation therapy



▶ Highlights

- 5-day radiation therapy after breast lumpectomy instead of the conventional 5-7 week whole breast radiation
- Multi-lumen device gives radiation oncologists the ability to shape the radiation does for typical cases
- Treats patients who are otherwise not appropriate candidates for traditional brachytherapy
 Targets area where cancer is most likely to recur
- Spares healthy tissue and organs from the side effects of radiation

▶ Hologic Digital StereoLoc® II stereotactic upright biopsy



- Precise needle guidance
- Provides easy transition from Selenia® digital mammography to stereotactic breast biopsy

▶ Hologic SecurView® diagnostic workstation



▶ Highlights

- Review digital images from mammograms, MRI, PET, and ultrasound
- Flexible, intuitive image review capabilities
- Unlimited hanging configurations Work interactively and intelligently
- through information-sharing
- Supports the display, manipulation and interpretation of DICOM images from multiple modalities including ultrasound, MRI, PET and nuclear medicine

▶ Hologic R2™ computer-aided detection (CAD)



▶ Highlights


- R2 CAD sophisticated pattern recognition software proven to help find breast cancers at an earlier stage
- Identifies features and brings them to the radiologist's attention in order to decrease false negative readings
 Often compared to a second pair of eyes, it serves as an interpretive aid during image review

▶ Hologic ATEC® automated tissue excision and collection system



Highlights

- -One system for use under all 3 imaging modalities
- (stereotactic, MRI and ultrasound)
 Consistently larger and more complete specimens for a confident diagnosis
- Fast tissue acquisition every 4.5 seconds
- -Safe, fully closed system limits fluid exposure
- Saline lavage irrigates and aspirates the cavity
 Pain medication can be continuously delivered direct to the biopsy site



Hologic has taken another significant leap forward in breast imaging with the introduction of Selenia[®] Dimensions[®] digital mammography system, the first practical tool to deliver on the extraordinary promise of breast tomosynthesis. Selenia Dimensions is the embodiment of state-of-the-art technology and flexibility, delivering:

- Revolutionary fusion imaging capabilities that combine digital mammography and tomosynthesis images in one exam, under one compression
- Exceptionally sharp images that allow you to visualize the finest details
- A sophisticated, ergonomic design for streamlined workflow and enhanced comfort for both the technologist and the patient

Learn More call +32.27.11.46.95 e-mail womensimaging@hologic.com visit www.hologic.com





▶ Hologic Eviva™ vacuum-assisted breast biopsy device



▶ Highlights

- A fast, safe and simple breast biopsy device that removes suspicious breast tissue for pathological evaluation
- The only vacuum-assisted device with a fully integrated end-deploy biopsy site marking solution

▶ Hologic MammoPad® radiolucent breast cushion



▶ Highlights

- Enables 5 million women to experience a softer, warmer mammogram
- Because patient is relaxed, technologists and radiologists often find better tissue acquisition

 – Does not create a need for increased dose
- Does not compromise high level of image quality needed for a routine mammogram

▶ Hologic MultiCare® Platinum prone biopsy table



▶ Highlights

- Minimally invasive stereotactic breast biopsy table
- Pinpoint accuracy– Enhanced patient comfort and optimized efficiency
- Superb image quality
 Adjustable breast tray compensates for breast thickness
- Cartesian coordinates system ensures accurate targeting
- True 360 degree access

▶ Hologic SecurView® technologist workstation




review capabilities Work interactively and intelligently through information-sharing

► Highlights - Powerful technologist workstation for the

mammography suite

Fast access to patient images
- Flexible, intuitive image

▶ Hologic Quantra™ volumetric assessment for breast density



- A break-through technology that estimates volumes in the breast and calculates the volumetric fraction of fibroglandular
- Aggregates volumetric measurements from each view in a study into a simple, concise assessment for each breast

▶ Konica Minolta Regius Pureview M

Anode Mo Mo/Rh Filter kV Range 20 - 35



– Revolutionary new x-ray mammography system based on phase contrast technology

-Sharpness and spatial resolution highly improved by the use of phase contrast technology

Reading at 43.75 µm thus equivalent to resolution of around 70 million pixels
Flex AEC 48 independent

detectors



▶ Konica Minolta Acies Mammo Workstation



▶ Highlights

- Easy to use diagnostic workstation for digital mammography
- Advanced configurable hanging protocol
 Multi-viewing modality
- DICOM compliance
- Software upgrades included

▶ medigration MammoView



▶ Highlights

- Extremely easy to use and manage
- Access to all images (including previous images) in seconds
- -Direct findings in the image
- Hanging protocols can be configured individually to automate your routine procedures
- Outstanding image quality (2048 greyscale)
- Excellent price/performance ratio

▶ Philips MammoDiagnost DR

Technology Resolution Size

Amorphous Selenium 85 µm 24 x 30 cm





- Comfortable and efficient workflow thanks to the intuitive Eleva User Interface and ergonomic design award winning system
- Excellent, UNIQUE-processed images to help with diagnosis
- Smooth procedures make your patients feel at ease

▶ Philips MammoDiagnost

Anode Mo/W (Dual Track) Filter Mo/Rh kV Range 23 - 35

▶ Highlights

- Efficient workflow for screening, diagnosis and interventions
- Perfect image quality to help with a diagnosis
- Smooth procedures make your patients feel at ease

▶ Philips PCR Eleva CosimaX

Technology Size

Resolution

Dual-side reading 18 x 24 cm / 24 x 30 cm HR-BD cassette/plate up to 50 lm pixel size



▶ Highlights

- Cost-efficient path to go digital
 Intuitive workflow with the Eleva User Interface
 Excellent image quality with UNIQUE image processing
 Multi-slot reader with up to 80 plates per hour
 Multi-purpose reader for mammography and other
 applications such as pediatrics and extremities

▶ Planmed Nuance Excel

Anode Filter

Mo or W (optional) Mo/Rh or Rh/Ag (optional)

kV Range 20 - 35



▶ Highlights

- Low dose FFDM Unit with 23.9 x 30.5 cm a-Se detector and fully automatic Flex-AEC with tissue type recognition Acquisition Workstation (AWS) with 3 megapixel
- TFT monitor and optional Nuance Acquire Station with motorized height adjustment

 - Integrated MaxView Breast Positioning System

 - Side Access for optimal patient positioning and ergonomics

- Optional stereotactics with Nuance Excel DigiGuide

▶ Planmed Nuance

Anode Mo or W (optional) Filter Mo/Rh or Rh/Ag (optional)

kV Range 20 - 35

▶ Highlights

- FFDM Unit with 17.1 x 23.9 cm a-Se detector and fully automatic Flex-AEC with tissue type recognition

- Acquisition Workstation (AWS) with 3 megapixel TFT monitor and optional Nuance Acquire Station with motorized height adjustment

- Integrated MaxView Breast Positioning System

- Side Access for optimal patient positioning and ergonomics Optional stereotactics with

Nuance DigiGuide



▶ Planmed Nuance Classic

Anode Mo Filter Mo/Rh kV Range 20 - 35

▶ Highlights

- High-end analog mammography unit with Flex-AEC

Field upgradeable to Full Field Digital Mammography

Side Access Patient

Positioning

Optional MaxView Breast Positioning System

Film or digital stereotactics system DigiGuide available as an add-on

Network ID camera and CR interface available



▶ Planmed Sophie Classic MVR

Anode Filter Mo/Rh kV Range 20 - 35



- Versatile mid-tier film unit with multiple options

Optional Flex-AEC with tissue

type recognition

 – Öptional MaxView or TwinComp compression systems

 – Öptional magnification and film or digital stereotactics

- Optional CR interface

▶ Planmed Sophie Classic Mobile

Anode Filter Mo/Rh kV Range 20 - 35



▶ Highlights

 Robust, mobile analog unit with integrated, telescopic radiation protection screen

- Optional battery backup
- Active brake system with lockable front wheels
- Versatility with optional Flex-AEC, magnification,
MaxView Breast Positioning System or TwinComp and film or digital stereotactics

- Network ID camera and CR interface available

▶ Planmed Sophie Classic S

Anode Мо Mo/Rh Filter



▶ Highlights

- Unique photon counting detector


-50% dose reduction compared with other FFDM systems

design enabling unsurpassed throughput

Complete screening solution with Sectra Breast Imaging PACS







MICRODOSE MAMMOGRAPHY. NO IMAGE DISTORTIONS TO HAUNT YOU

Ghosts belong in creaky old buildings. Not in X-rays. Our patented photon-counting technology effectively banishes them. The detector is fast enough to be ready when the next photon arrives, so there's nothing to interfere with image interpretation. Along with ghosts, you can banish

inefficiencies, too. Specially designed for the smoothest possible workflow, you'll find MicroDose is a superb allround screening solution.

For more on the highest image quality at the lowest dose, visit sectra.com/medical



▶ Sectra mammography workstation

Technology Resolution Size

Photon-counting 50 μm, 14 bit 24 x 26 cm



▶ Highlights

- A multimodality workstation with dedicated software tools to meet the special workflow and throughput requirements of mammography

- Fast image display Dedicated keypad

- Automatic display protocols

▶ Siemens Mammomat Inspiration*

W/Rh, a-Se Technology Resolution 85 µm Detector size 24 x 30 cm

▶ Highlights

-Screening, upgradable to stereotactic biopsy and 3D imaging with tomosynthesis

 Dual target anode W/Rh
 reduces dose up to 50% especially for dense breasts

Comprehensive system solution with syngo-based acquisition workstation Streamlined workflow: One-click-to-image

- Special MoodLight function

Broad range of accessories

*not available in the USA



▶ Siemens Mammomat 3000 Nova

Anode Mo/W Filter Mo/Rh kV Range | 23 - 35

▶ Highlights

- Pivoting buckys, easy switching between 18 x 24 and 24 x 30

- Prepared for stereo biopsy

- Opdose auto-selects best anode/filter combination (Mo/Mo, Mo/Rh, W/Rh)

 Opcomp – Siemens'
- exclusive optimized compression system
- Optional digital spot imaging with syngo Opdima



▶ Siemens Mammomat Inspiration with 3D Tomosynthesis*

Technology W/Rh, a-Se Resolution 85 µm 24 x 30 cm Detector size

▶ Highlights

- Platform for multiple mammography applications: Screening, diagnostics, stereotactic biopsy and Tomosynthesis in one system and one Acquisition Workstation

3D imaging via the acquisition of breast images taken at multiple angles $(+25^{\circ} \text{ to } -25^{\circ})$: Improved capability to diagnose especially very dense breasts

The only installed system which offers all applications and can be upgraded to 3D tomosynthesis

The largest angular range in industry increases depth resolution and contrast *not available in the USA



▶ Siemens Mammomat 1000

Anode Mo Mo/Rh Filter kV Range 23 - 35



▶ Highlights

 Cost-efficient system, ideally suited for screening

- Opcomp – Siemens' exclusive optimized compression system

-Soft speed - two speed compression

▶ Siemens Acuson S2000 Automated Breast Volume Scanner*

Technology Ultrasound

Size 15.4 cm x 16.8 cm Transducer



▶ Highlights

Acquisition of full-field volumes

Acquisition of rull-field volumes of the breast automatically, quickly and comfortably
 Reporting on ABVS Workplace - Intuitive volume image analysis and manipulation, comprehensive BI-Rads reporting
 Ideally suited to image patients with dense breast tissue and/or a history of breast disease
 Patient friendly - minimal compression

No radiation

▶ Siemens syngo MammoReport

CPU Displays User Interface Dual-core Xeon 5130 2.00 GHz or faster Dual high-contrast 5 MP monitors, 2.5 x 2 k Dedicated keypad customized to user's individual workflow



- Customization to personal workflow and image arrangement, controllable with a single key Supports CAD and a
- CAD-driven workflow
- Multi-modality and third-party viewing, including viewing of 3D tomosynthesis, 3D ultrasound, and MR breast images
- Advanced display layouts for different demands
- DICOM and ÎHE workflow compliance

▶ Agfa HealthCare iCAD SecondLook Digital



▶ Highlights

- SecondLook^R Digital is a second opinion tool for Agfa HealthCare's CR Mammography
- Adds value in daily routine
- Can be used in mammography screening as well as by diagnostic HealthCare providers
- Markers that indicate potential microcalcifications and/or masses are visible on Âgfa HealthCare's IMPAX Breast **Imaging Workstation**

▶ IBA Dosimetry Pasmam

Test device for checking spatial resolution, contrast resolution, signal to noise ratio, dynamic range, image limitation towards the chest wall, AEC performance



▶ Highlights

- Modular construction
- Different test inserts
- -Basic plate with Al step wedge
- -Structural plate with turnable spatial resolution test
- Additional attenuation plates

▶ Konica Minolta CP-1M

Technology Resolution Size CsBr phosphor 43.75 μm

18 x 24 cm / 24 x 30 cm



- ► Highlights

 Crystal column technology
- -Superior DQE
- Superb sharpness, granularity and stability
- Excellent image quality with reduced exposure dose
- Compatible with Regius 110HQ, Regius 190 or Regius 210

▶ PTW Normi PAS

Modular test object for quality control of digital mammography X-ray units



▶ Highlights

- Checks all imaging quality parameters (high contrast, low contrast, spatial resolution, signal-to-noise ratio, dynamic range, artefacts, thoracic wall side limitation, etc.)
- Test element based on the Mammographic Accreditation Phantom of the ACR included
- Acrylic absorbers for AEC testing included

▶ PTW Diados E Dosemeter

High sensitive dosemeter for absolute dosimetry, acceptance testing and quality control



▶ Highlights

- Measures dose, dose rate, dose/pulse, pulses, dose length product, irradiation time
- Wide dynamic measuring range
- New mammography qualities like
 Mo/Mo, Mo/Rh, Rh/Rh, W/Rh, W/Al, etc. available
 Data evaluation by means of the DiaControl *expert* software

▶ Siemens syngo MammoCAD*

CPU | Dual processor Intel Xeon 3.60 HT/800, 1 MB, EM64T

▶ Highlights

- Advanced image processing capabilities with state-of-the-art pattern recognition technology
- Up to 4 DICOM input connections
- Up to 10 DICOM output connections
- CAD processing of a four-image case within less than 90 seconds
- Designed for MAMMOMAT Novation and MAMMOMAT Inspiration*

*not available in the USA





Please see us at ECR, first level, booth 638



Enterprise PACS Solutions



JiveX Radiology

- ▶ Advanced Reading Protocols
- Presentation Workflow
- ▶ Integrated Teleradiology
- ▶ Flexible Archiving

JiveX / PACS and beyond

Your VISUS Sales Team Phone +49 234-936 93-400 E-Mail: sales@visus.com





Shimadzu - Pioneers of Vision

or more than 100 years, Shimadzu has been counted among the most innovative companies in medical technology. Already 10 months after the discovery of X-rays, Shimadzu recorded the first X-ray image in Japan. To this day, the company is one of the pioneers in diagnostic imaging.

Shimadzu focused extensively on DR systems and introduced additional novel diagnostic and interventional radiology products. An essential element of Shimadzu's current product range is the safire flat-panel detector (FPD), which is based on the principle of direct conversion, and this new technology has been further developed for radiography, fluoroscopy, angiography and cardiology systems.

Safire – groundbreaking flat-panel detector technology

The dynamic safire FPD – suitable for Radiography, Fluoroscopy as well as Angiography – offers excellent image quality with an outstanding detailed resolution.

Unlike the safire FPD indirect-conversion FPDs or conventional image intensifiers need several intermediary steps in order to generate an electrical output signal that can be digitally processed. Each of these steps deteriorates the signal-to-noise ratio and therefore the image contrast. This actual loss of information can only be compensated by a higher radiation dose.

Reduced dose and optimized image quality

The RSM-DSA technology

Using the RSM-DSA (Real-time Smoothed Mask DSA) technology patented by Shimadzu, DSA-typical motion artifacts (breathing, heartbeat or patient-induced movement) as well as the noise component in fluoroscopy images are significantly reduced. Multiple-bolus dosing, essential in conventional systems, is superfluous using RSM-DSA. High resolution fluoroscopy images can be generated within shorter patient examination times

– all this at lower radiation doses and lower contrast agent use.

The SUREengine

SUREengine is a real-time image processing software package for the safire flat-panel detector, specifically designed for image improvement and dose reduction. For digital data processing the frequency spectrum is first decompounded into wave bands. The noise components are then subtracted from the individual frequency bands and the image contrast is selectively increased. In order to broadly reduce blooming and to display evenly illuminated image backgrounds, the entire dynamic range is compressed. Moreover, Moiré patterns resulting from scattered radiation grids are filtered out of the original images without any loss in image quality. In this way, for instance, guide wires and catheters become much more clearly distinguishable.

State-of-the-art diagnostic imaging technology

With the introduction of the safire technology, Shimadzu has opened the gateway to new applications such as digital tomosynthesis, slot radiography, dualenergy subtraction and CT-like imaging.

Digital tomosynthesis

Digital tomosynthesis is the next generation development of imaging tomography. Unlike other sectional imaging techniques, tomosynthesis can also be applied to patients in an upright position. This is particularly helpful for functional and orthopedic diagnostics. The advantages of DR tomosynthesis in comparison with computer tomography are considerable. Based on the lower patient radiation dose required, this procedure is suitable especially for lung cancer screening.

Slot radiography

During slot radiography a sequence of 5 cm wide, adjoining bands (slots) is ac-

quired using almost closed collimator geometry and pulsed radiation, while X-ray tubes and flat-panel detector move linearly in one direction. After examination, the individual image bands are digitally reconstructed into a so-called long image (long leg/spine). Therefore, this procedure is used in scoliosis and orthopaedics examinations in order to completely map the spine or the lower extremities.

The component of the elicited image noise due to stray radiation is significantly reduced in slot radiography. The use of a grid is no longer necessary.

Dual-Energy Subtraction

Dual-energy imaging enables fast image acquisition using alternating high and low value kV-pulses. During post processing, the resulting images are subtracted from each other and a soft-tissue image as well as bone image can be calculated after previous standardization with respect to bone or lung values. The main advantage of this procedure is the clear representation of calcifications, extremely helpful in the detection of lung nodules.

CT-like imaging

CT-like imaging closes the gap between angiography and computer tomography systems in interventional radiology.

The digital data sets acquired during C-arm rotation can be reconstructed – similarly to computer tomography – to represent delicate three-dimensional vascular structures. In comparison with a 128-line CT, CT-like examinations using a flat-panel based C-arm can lead to resolution improvements of up to 200%, while the patient dose can be simultaneously reduced up to 90%.

Shimadzu Europa GmbH


Albert-Hahn-Str. 6-10 47269 Duisburg · Germany Phone: +49 (0) 203-7687-0 Fax: +49 (0) 203-7687-680 E-mail: medical@shimadzu.eu www.shimadzu.eu











▶ Philips BuckyDiagnost High-performance room

Design Ceiling-mounted tube carrier for standard bed exposures Table Fixed or hight adjustable (optional) with various table top sizes 30 - 85 kW with several options



▶ Highlights

Optional generator functionalities are AEC, APR, automatic collimation, tracking, and tomography

- Flexible, interchangeable components with a large range of tables, stands, tube carriers, tubes and generators

Same handling and options for floor-mounted system and ceiling-mounted tube carrier

Future-proof with digital upgrade possibilities via PCR Eleva or digital detector
- Ergonomic design enabling easy handling and near

patient control

▶ Philips BuckyDiagnost Standard room

Design Floor-mounted system (floor, floor/wall or floor/ceiling) Table Fixed or hight adjustable (optional) with various table top sizes Power 30 - 85 kW with several options



▶ Highlights

Optional generator functionalities are AEC,


APR, automatic collimation, tracking and tomography

– Flexible, interchangeable components with a large range of

tables, stands, tube carriers, tubes and generators

Same handling and options for floor-mounted system and ceiling-mounted tube carrier


- Digital upgrade possibilities via PCR Eleva or digital detector Ergonomic design enabling easy handling and near patient control














- Flat-panel detector upgradability

centred all the time

▶ Siemens Multix Swing

Design Floor-mounted

Table Weight capacity of 450 kg 30 or 55 kW Power



▶ Highlights

- Cost-efficient, all-in-one X-ray room solution
- Generator is integrated into the table for minimal space requirements
- Floating tabletop with weight capacity of up to 450 kg
- Synchronized tube and bucky tray movements
- Tube can be rotated for cross-table exposures



- Conventional radiography solution with integrated X-ray tube
- Robust and easy-to-use, accommodating a wide range of
- clinical imaging
 Wide range of tabletop movements with table height adjustments
- Automatic exposure control

▶ Siemens Multix Top

Design 3D-ceiling-mounted Height adjustable Table Power

30, 55, 65, 80 kW



▶ Highlights

- Ceiling-mounted, conventional radiography solution
- -Robust and a real workhorse, suitable for high-throughput radiography rooms
- Easy, colour-coded ceiling tube handling
 TOP alignment of x-ray field for dose reduction during chest examinations
- Upgradeable to flat detector technology

▶ Toshiba Radrex-I

Design Ceiling mounted Ceiling mounted Dual Panel System Single Panel System 1 portable FPD 2x fixed FPD

Power 80 kW 80 kW



▶ Highlights

- Easy Operation for high throughput

- In Room Image Review
 Big size OTC Interface for easy Operation
 Automated FollowME Concept e.g. Auto Tracking, -processing, -image distribution, -image stitching, etc.
- Full Dicom functionalities included
- Optional wired portable FPD

R/F SYSTEMS MOBILE

▶ GE Healthcare AMX 4+

12.5 kW Power kV Range 50 - 125 mAs Range 0.4 - 320



▶ Highlights

- High usable battery power storage
- Wide mAs range for variety of applications
- Ûnique column rotation
- AEC for consistent image quality Excellent maneuverability with
- motor drive

▶ GE Healthcare TMX+ / TMX R+ TMX R+ TMX+ Power 30~kW30 kW kV Range 40 - 12540 - 125 mAs Range 0.2 - 2200.2 - 220▶ Highlights - Powerful system for variety of applications Anatomic programs - Dose level selection – Dual focal spot: 0.8 and 1.3 mm -TMX R+: column rotation for easy positioning

▶ Philips Practix Convenio

30 kW Power kV Range 40 - 125

mAs Range 0.63 – 320 (large focal spot) 0.1 - 200 (small focal spot)

▶ Highlights

 Robust electrical and technical concept with swiveling column and telescopic tube arm


- Powerful motor drive for outstanding maneuverability on the spot, on ramps and over obstacles

 Long-life batteries due to intelligent single charging management

- Ergonomic design and intuitive user guidance – winner of the iF and the I.D. design award 2006

 Digital with Philips Computed
 Radiography and UNIQUE image processing










▶ Shimadzu MobileArt eco Power 12.5 kW kV Range 40 - 1250.32 - 100 (200) mAs Range





▶ Highlights

- Telescopic arm Easy positioning Wide coverage









- Lightest system in its class, offers

-Supports a wide range of applications

outstanding maneuverability



▶ Highlights

AutoEx – fully

preferences

flexibility


parameter optimization


Innovative user interface

Seamless digital workflow

Efficient dose management - High resolution image chain








Offers ultimate customization by adjusting the system to

180 degree isocentric c-arm rotation increases projection

– Proven scan principle with c-arm moving around the patients





- ▶ Philips DuoDiagnost Design Free arm system II-format 23/31/38 cm Image system ▶ Highlights
- -2-in-1 system for digital radiography and fluoroscopy which enables you to perform a wide range of digital X-ray applications

 Patient capacity up to 284 kg in all movements combined with
 a large table top size allows for bariatric examinations

 Exams on the table with Source Image Distance (SID) of 180 cm
- Minimum table top height of 62 cm facilitates easy transfer and positioning of patients
 -Scan range of 203 cm without table top movement enables
- examinations without repositioning for increased patient comfort







system to preferences













▶ Highlights

 Complete patient coverage with 8-way tabletop travel and large receptor movements

Single-handed cassette handling:

automatic loading, centering, format sensing and collimation

– FLUOROSPOT Compact high-resolution digital imaging
system with intuitive user interface and DICOM 3.0 interfaces

Seamless integration into DICOM network

- Dose-saving fluoroscopy with SUPERVISION (option)

- Bucky wall stand (option)

▶ Siemens Axiom Iconos R200

Design Technology Resolution II-format Image system Size

Universal digital remote-controlled R/F system 1 k x 1 k matrix

33 cm or 40 cm



Highlights

- Modular design for universal applications including diagnostic angiography (option)

- Complete patient coverage with 8-way tabletop travel and large receptor movements

Single-handed cassette handling: automatic loading, centering,

- Single-nanded cassette handling, date like the format sensing and collimation
- FLUOROSPOT Compact high-resolution digital imaging system with intuitive user interface and DICOM 3.0 interfaces

Seamless integration into DICOM network

- Comprehensive CARE dose reduction package

▶ Siemens Luminos RF Classic Design Remote-controlled R/F system Technology 1 k x 1 k matrix Resolution II-format 23 cm or 33 cm Image system Size

▶ Highlights

 Complete patient coverage with 8-way tabletop travel and large receptor movements

- Single-handed cassette handling: automatic loading, centering, format sensing and collimation

Intuitive and fast operation with innovative control console

– Dose-saving fluoroscopy with SUPERVISION (Option)

- Bucky wall stand (Option)

- Excellent price-performance ratio

▶ Siemens Axiom Luminos TF

Design Technology Resolution II-format Image system

Size

Digital tableside controlled R/F system 1 k x 1 k matrix

33 cm or 40 cm



▶ Highlights

One digital system for

fluoroscopy and radiography from pediatrics to bariatrics - Open design to accommodate bariatric patients up to 272 kg
- Ergonomic Single-handed system operation with OPTI Grip handle

- Ergonomic Single-Handed system operation with 17 City Special System with intuitive user interface and DICOM 3.0 interfaces Comprehensive CARE dose reduction package

Mobile flat detector (option) for fully digital radiography workflow

▶ Siemens Axiom Artis U

Universal, floor-mounted Design II-format 23 cm or 33 cm Image system 1024 x 1024 pixel, 12 bit-CCD



▶ Highlights

High-power output for excellent image quality
High heat capacity x-ray tube virtually eliminates overheating issues

- Compact and »room-mobile« design

- Broad application spectrum

Excellent price-performance ratio

▶ Siemens Axiom Luminos dRF

Design Technology Resolution II-format Image system Remote-controlled 2-in-1 system with dynamic flat detector Amorphous-Silicon with Cesium Iodide scintillator

Up to 3.4 lp/mm

43 cm x 43 cm



- Fully digital 2-in-1 solution for dynamic and static


high-resolution imaging including DSA procedures (option)

 Easy patient transfer at 48 cm lowest table height
 Dynamic Density Optimization (DDO) and DiamondView Plus for excellent detail contrast

-FLUOROSPOT Compact high-resolution digital imaging system with intuitive user interface and DICOM 3.0 interfaces

 Comprehensive CARE dose reduction package
 Limitless projection flexibility with optional ceiling-suspended tube and wireless detector wi-D





▶ Toshiba Zexira True 2 in 1 System for fluoroscopy and radiography Design Detector ▶ Highlights Full clinical flexibility: G.I; Venography; Urology; ERCP; Angiography; Radiography, Oblique Imaging, etc. – Easy Userinterface

R/F SYSTEMS ACCESSORIES



▶ Highlights

- Accurate and precise bone density measurements of the spine, hip, whole body and forearm
- Provides the industry's highest resolution vertebral imaging which dramatically improves the detection of vertebral fractures
- Vertebral imaging (IVA-HD) capability can also be used for the assessment of aortic calcification, a significant indicator of heart disease

▶ IBA Dosimetry Primus

Test device for checking image quality parameters at fluoroscopic units

-Full range of Dose optimization techniques



Highlights

– Modular construction:

- Full patient coverage

- structural plate and separated attenuator
 Check of spatial and contrast resolution, size of the radiation field, artefacts; kV test area
- Compact Al pre-attenuator or PMMA and Cu plates
- Available in two different sizes

▶ IBA Dosimetry DIGI-13

Test device for checking image quality parameters at digital radiographic units



▶ Highlights

- -Compact device with separated Al pre-attenuator
- With integrated copper plate
- Check of homogeneity, spatial and contrast resolution, size of the radiation field, artefacts
- Easy-to-use

▶ IBA Multimeter MagicMax

Simultaneous measurement of dose, dose rate, exposure time. kV, dose/pulse, pulse rate, HVL and total filtration



▶ Highlights

- Small device with separate multifunction detector
- Connected via USB to PC or Notebook
 Intuitive use via PC interface
 Time resolution: 100 µs

- Optimized solutions for all applications

▶ Provotec ProVario Screen

Design Table Power Power unit floor-mounted Control console placed on a desk

50 kW



► Highlights

- High frequency generator for x-ray diagnostic

- Easy operation by monitor or touchscreen
 Digital control of nearly unlimited organ programs
 Safety device against undue radiation for each organ with AEC-technique
- X-ray book for storing patient name with generator exposure data
- Upgradeable for using CR- and DR-systems

▶ PTW Diamentor CM

Miniature dose area product (DAP) meter for patient dosimetry and quality control



Highlights

- Compact solution ideal for integration in mobile units
 Built-in test function for fast calibration and constancy checks
- Easy connection to a RIS or PACS

▶ Radcal ACCU-PRO[™]

X-Ray Analyzer

Simultaneous dose, rate, time, kVp. HVL, filtration, mA/mAs, and more



▶ Highlights

- Use for manufacturing, installation, QA, and service R/F, mammography, CT, dental, leakage Ion chamber based dosimetry, no corrections required Correctly measure AEC fluoro and filtered beams

- Remote control, waveforms, and archiving with XLPRO software
- Compact, easy to use

▶ Radcal RAPIDOSE

PC X-ray Analyzer



Highlights

- Plug into a laptop USB port for an inexpensive X-ray analyzer
 Simultaneous dose, rate, kVp, time, HVL, waveforms, and
- Revolutionary inherent remote measurement operation
- Easy use, genuine time saver
- -Data archiving and analysis using your Excel

▶ Radcal PDC-DAP/KAP verification meter



Highlights

- Newly patented Patient Dose Calibrator
- Use to calibrate DAP/KAP and rate
- Also measures dose and rate Optical and radiographic alignment markers
- Simple to use with optional computer control

▶ RTI Electronics Piranha

The Piranha is designed as a truly self-contained, all-in-one, X-ray multi-function meter that assures accurate results in one shot. kV, time, dose, dose rate, HVL and total filtration



Highlights

- Self-Contained, All-in-One
- Auto-Compensation

- Auto-Compensation
 R&F, Mammo, Dental and CT
 Quick and Simple Set-up
 Enhanced Graphical Display
 Built-In Bluetooth for PC and PDA
- -mA, mAs, and Light Probes
- Fits in the Palm of Your Hand

▶ RTI Electronics Barracuada

The Barracuda X-ray multimeter has a cabinet that can house up to six different application modules, and can measure on all modalities; R/F, mammography, flouroscopy, pulsed flouroscopy, dental, panoramic dental and CT systems



Highlights

- All in One, All at Once
- Auto-Compensation
 Enhanced Graphical PDA Display
- R&F, Mammo, Dental and CT
- Ionization Chambers
- -Built-In Bluetooth for PC and PDA
- -mAs, and Light Probes
- -Fits in the Palm of Your Hand

▶ ulrich medical – CO₂ Insufflator for virtual colonoscopy

Pressure | Insufflation rate

0-30 mmHG, infinitely variable, preselectable 1-4 l/min, arbitrary

Setting

supported by voice confirmation system



▶ Highlights

- Automatic insufflation of CO₂ into the colon for virtual colonoscopy examinations in CT

- Significant improvement of diagnostic results compared to manual room air insufflation

 Increase of patient comfort due to automatic adjustment of over pressure and faster resorption

– Easy setting of gas volume and pressure

– Display of gas consumption

- Four adjustable flow rates



RAD·BOOK 2010



























0:R Technology

C









PHILIPS











World's First Purpose-built Direct Digital RSA X-Ray System

> opened by Susan Hampshire **OBE**

pioneering orthopaedic research direct Digital Radiography system (dDR), the 'Adora RSA', has been officially opened by actress, author and osteoporosis campaigner Susan Hampshire OBE at Aberdeen University.

Xograph Healthcare, the leading independent supplier of diagnostic imaging equipment in the UK and Republic of Ireland helped facilitate the design and construction of this state-of the-art system through collaboration with Aberdeen University, Danish manufacturer Nordisk Røntgen Teknik and market-leading direct Digital Radiography pioneer Canon Inc.

During her visit, Susan Hampshire spoke to members of staff at Aberdeen University and listened intently as they discussed how 'Adora RSA' will greatly assist researchers performing radiostereometric analysis (RSA). This technique is a precise method for determining the migration and wear of orthopaedic implants such as hip and other joint replacements. Two simultaneous X-rays exposures are performed with images being captured on flat panel Canon dDR detectors. The resulting images are analysed using software, developed by specialists at the University, which determines the three-dimensional movements of parts of the body, such as bone fracture components, in relation to one another, or a surgical implant in relation to the host anatomy. The information gathered from RSA demonstrates precisely how much positional change has occurred for example since implantation of a prosthesis. This may help predict how long implants will last, RSA will also benefit future patients by providing important data to help researchers improve joint replacement technology and surgery.

'Adora RSA', tailored to Aberdeen university researchers exacting requirements, is causing quite a stir in the field of radiostereometric analysis and other pioneering authorities in this area of research have expressed an interest in visiting this revolutionary new system for themselves. Whilst other specialist RSA centres are using more traditional forms of technology, Aberdeen University is the world's first that uses a purpose-built direct digital RSA imaging technology gen-

erating high quality digital images that can be viewed just three seconds after exposure. Now, leading researchers at Aberdeen will be able to gain faster results leading to earlier conclusions that will enable more accurate assessment in the healing of broken bones and joint replacement surgery. This in turn will improve the effectiveness of other research and clinical investigations.

contact our representatives in UK: Xograph Healthcare Ltd Hampton Street · Tetbury Gloucestershire · GL8 8LD · UK Tel: +44 (0)1666 501501 Fax: +44 (0)1666 501502 canon@xograph.com

For more information please

▶ Agfa HealthCare DX-S

Slots Capacity Resolution

Up to 130 plates/h 50 μm

▶ Highlights

- High throughput decentralized digitizer, resulting from two major technology shifts (DirectriX and Scanhead)

-Bridges the gap between CR and DR, offering the best of both worlds

- Precision image quality

Dose reduction

-Broad range of applications: general radiography, pediatrics and emergency



▶ Agfa HealthCare CR 85-X

Slots Capacity Resolution 1 – 10 cassettes: drop and go buffer

112 plates/h

50 - 100 - 150 μm



▶ Highlights

- Multi-user digitizer for centralized use

 Unique drop and go buffer system

Broad range of applications: mammography (outside US), orthopedics, dentistry, pediatrics, general radiography

▶ Agfa HealthCare CR 35-X

Slots Capacity

Up to 71 plates/h 50 – 100 – 150 μm Resolution



- Multi application digitizer

- Ideal for decentralized environments

-Broad range of applications: radiotherapy, mammography diagnostic and screening (outside US), orthopedics, dentistry, pediatrics, general radiography
-Three different resolution

modes



▶ Agfa HealthCare CR 30-X

Slots Capacity Resolution

Up to 82 plates/h 10 pixel/mm



Highlights

- Tabletop digitizers

- Broad range of applications: general radiology, orthopedics, chiropractic, dentistry

Low total cost of ownership

Mobile use

▶ Agfa HealthCare DX-G

Slots Capacity 1-5 cassettes: drop and go buffer

Approx. 83 plates per hour (35x43cm cassette) Resolution

6.7-10 pixels/mm

▶ Highlights

- Allows use of both standard phosphor plates and needle-based detectors

- Has the potential for dose reduction for all studies – particularly neonatology and pediatrics

High throughput, delivered by a unique five cassette drop-and-go buffer, and a very fast preview.

- Housed in a compact system - Provided with infuitive NX workstation with the gold standard image software, MUSICA2



Slots Capacity 1-5 cassettes: drop and go buffer

Approx. 83 plates per hour (35x43cm cassette)

Resolution 6.7-20 pixels/mm

▶ Highlights

- Allows use of both standard phosphor plates and needle-based detectors Able to read NIP for

Mammography requirements as well as NIP for general

radiography

- High throughput, delivered by a unique five cassette drop-and-go buffer and a very fast preview

Dedicated mammography

needle-based detector and

compatible with existing breast imaging modalities Provided with intuitive NX workstation with the gold standard image software MUSICA2





▶ Fujifilm FCR Profect CS

Slots Capacity Resolution

165 Imaging plates (IPs)/h 5 – 20 pixel/mm

► Highlights - EUREF & PAS 1054

compliant - First mammography CR system approved by FDA - Fastest mammography system

available

- 50 μm

- Needs 30% less dosage for pediatric exams - Worldwide more than 5.000

FCR Profect installed

▶ Fujifilm FCR XG5000

Slots Capacity

Resolution

165 Imaging plates (IPs)/h

5 – 10 pixel/mm



▶ Highlights

Worldwide more than 70.000 Fujifilm CR systems installed

Universal applicable
IHE certified

- Wide dynamic range

Optimized workflow

▶ Fujifilm FCR Profect One

Slots Capacity Resolution

85 Imaging plates (IPs)/h

5 – 20 pixel/mm



▶ Fujifilm FCR Capsula X / XL II

Slots

Capacity 72/94 Imaging plates (IPs)/h Resolution

5 – 10 pixel/mm



► Highlights - EUREF & PAS 1054 compliant

First mammography CR system approved by FDA
Needs 30% less dosage for

pediatric exams

Compact system

▶ Highlights

- Extremely compact system, mobile model available

- IHE certified - Universal applicable, wide

dynamic range

- Ideal for medium radiologists (e.g. orthopedic doctors)

Optimized workflow

▶ Fujifilm FCR Prima

Slots Capacity Resolution

up to 29 (IPs)/h 10 pixel/mm



▶ Highlights

Seamless image and patient data workflow

-Efficient use of space in examination room

Compact and quiet, Footprint 0.24 m^2

All-in-one workstation: Console, Viewer, Archive

Value for money



▶ Konica Minolta Regius 190 HPSMAMMO

Slots Capacity

188 plates/h

(14"x 14", 175 μm)

Resolution 175 μm/87,5 μm/43.75 μm

▶ Highlights

 Innovative technology and modern design

- Simultaneous reading

of 2 cassettes

- High throughput

performance

– High resolution mammo mode 43.75 µm for mammo cassettes

Normal & high quality mode for all standard cassettes $(175 \mu m/87.5 \mu m)$



▶ Konica Minolta Regius 190^{MAMMO}

Slots Capacity Resolution

94 plates/h (14"x 14", 175 μm) 175 μm/87.5 μm/43.75 μm



▶ Konica Minolta Regius 190

Slots Capacity Resolution

▶ Highlights

design

integration

cassettes

- Innovative technology and modern

-Flexible workflow

Normal & high quality mode for all

 $(175 \mu m/87.5 \mu m)$

94 plates/h (14"x 14", 175 μm) 175 μm/87.5 μm/43.75 μm



▶ Highlights

- Innovative technology and modern design
- High resolution mammo mode 43.75 μm for mammo cassettes -Flexible workflow
- integration
- Normal & high quality mode for all standard cassettes

 $(175 \mu m/87.5 \mu m)$

▶ Konica Minolta Regius 110

Slots Capacity

▶ Highlights - Very compact and flexible design

cassettes

solution

Normal & high

quality mode for all

(175 µm/87.5 µm) Integration in

Regius 190 network

- Cost efficient CR

80 screens/h (14 x 14 cm) Resolution

175 μm/87.5 μm/ 6 – 11 pixel/mm



▶ Konica Minolta Regius 110HQ

Slots

Capacity Resolution 80 plates/h (14" x 14", 175 μm) 175 μm/87.5 μm/43.75 μm

Highlights

- Powerful, compact reader with linear motor technology
- High quality Mammography read function
- Easy operability and maintenance
- Use with standard cassettes and/or CP-



▶ Konica Minolta Regius 210

Slots

Capacity Resolution

100 plates/h (14" x 14", 175 μm) 175 μm/87.5 μm/43.75 μm



- High performance dual bay reader
- Outstanding image quality in both general X-ray and Mammography
- -Low dose imaging for Pediatric use
- Use with standard cassettes and crystal column cassettes (CP-1M, CP-1S)



▶ Philips PCR Eleva CosimaX

Slots Capacity

165 cassettes/h (18 x 24 cm, standard readout) 80 cassettes/h (18 x 24 cm, dual-side readout)

5 – 10 pixel/mm;

Resolution

20 pixel/mm with dual-side reading

for HR-BD and

ST-BD cassettes

▶ Highlights

-Customizable Eleva User Interface combined with superb image quality by

UNIQUE image processing

– Features simultaneous dual-side reading for 18x24 and 24x30 imaging plates

- 40% increase of DQE and enabling high-resolution imaging like mammography (HR-BD cassettes/plates)

 Low-dose imaging for pediatrics (ST-BD cassettes/plates)
- Dedicated for mammography or pediatrics environments
- Orthopedic automatic image stitching



▶ Philips PCR Eleva Corado

Slots Capacity

165 cassettes/h (18 x 24 cm) 143 cassettes/h

(35 x 35 cm, in high-speed mode)

10 pixel/mm,

Resolution

5 pixel/mm in high-speed mode



▶ Highlights

Customizable Eleva User Interface combined with superb image quality by UNIQUE image processing

High-throughput, multi-slot system, for environments using a central reader set-up

- For general radiographic applications including orthopedics

Orthopedic automatic image stitching

▶ Philips PCR Eleva S Hi-res

Slots Capacity

Resolution


94 plates/h (18 x 24 cm, standard readout)

52 plates/h (18 x 24 cm, dual-side readout

5 – 10 pixel/mm;

20 pixel/mm with dual-side reading for HR-BD and

ST-BD cassettes



Highlights

Customizable Eleva User Interface combined with superbimage quality by UNIQUE image processing

Features simultaneous dual-side reading for 18x24 and 24x30 imaging plates

- 40% increase of DQE and enabling high-resolution imaging like mammography (HR-BD cassettes/plates)

- Low-dose imaging for pediatrics (ST-BD cassettes/plates)

- Orthopedic automatic image stitching

▶ Philips PCR Eleva S Plus

Slots Capacity

97 plates/h (18 x 24 cm);

Resolution

94 plates/h (35 x 35 cm, in high-speed mode) 10 pixel/mm, 5 pixel/mm in high-speed mode 20 pixel/mm with HR cassettes (option)



▶ Highlights

- Customizable Eleva User Interface combined with superb image quality by UNIQUE image processing

- For environments with high throughput requirements

- For general applications, including orthopedic and dental applications

Faster read-out than S reader with high speed mode

- Smaller footprint and dimensions

- 50 micron scanning option

Orthopedic automatic image stitching

▶ Philips PCR Eleva S

Slots Capacity

78 plates/h

Resolution

 $(18 \times 24 \text{ cm})$ 10 pixel/mm



▶ Highlights

Customizable Eleva User Interface combined with superb image quality by UNIQUE image processing

Suitable in environments with moderate performance requirements and/or decentral reader set-up

For general applications, including orthopedic and dental applications

Smaller footprint and dimensions

Orthopedic automatic image stitching

▶ Protec Proscan 35E CR-System

Slots Capacity

52 - 110 screens/h Resolution up to 20 pixel/mm



▶ Highlights

- 16 bit grayscale resolution

- Smallest physical pixel size is 12.5 µm

Can read IPs in odd formats, e.g.

18 x 35 cm for lumbar spine images and faster readout

CONAXX image acquisition software included in standard delivery

- Fully DICOM compatible

Independent modality or easy integration to PACS

▶ Protec Proscan 43 CR-System

Slots Capacity

55 - 80 screens/h Resolution Up to 20 pixel/mm



- 16 bit grayscale resolution

- Smallest physical pixel size is 12.5 µm

- Touch-free IP transport

 Installation table top or on cabinet (option)

- Extremely small footprint
- CONAXX image acquisition
software included in standard delivery

Fully DICOM compatible

Network-compatible

RFID chip for cassette identification and workflow optimization

- Independent modality or easy integration to PACS

PROSCAN₂₃

▶ Agfa HealthCare DX-Si

Technology DirectriX and ScanHead

Resolution 50 µm

35 x 43 cm, 24 x 30 cm, 18 x 24 cm, 15 x 30 cm Size

▶ Highlights

- Integrated digital radiography solution

— In room integration/

workflow for all general radiography applications

 Flexibility to create multiple tailor made configurations

-High throughput system: up to 130 plates/h

Dose reduction



▶ Agfa DX-D 300 HC

64 kW Power

kV-Range From 40 to 150 KVp in 1 kVp step mAs-Range From 0.1 to 500 mAs in 38 step

▶ Highlights

Universal direct radiography system, designed for maximum flexibility.

Small footprint enables imaging solution in a single room

Provided with integrated NX workstation for connectivity with RIS/HIS/PACS and imagers

Superior contrast detail provided by MUSICA2 processing producing exam independent consistent image quality
Cesium Iodide DR detector technology

providing significant patient dose reduction potential

▶ Agfa HealthCare DX-D 500

Power 80 kW

kV-Range 40 kV to 150 kV from 0.5 mAs to 800 mAs mAs-Range

Highlights

- Can handle a very comprehensive range of X-ray exams, from general radiography to emergency and pediatrics

- Two detector high productivity

General Radiography DR System

Provided with integrated NX
workstation for connectivity with RIS/HIS/PACS and imagers

Superior contrast detail provided by MUSICA2 processing producing exam independent consistent image quality
- Cesium Iodide DR detector

technology providing significant patient dose reduction potential



▶ Apelem Da Vinci Solo — universal single detector solution

Technology CsI Resolution 143 µm Size 43 x 43 cm



- For gen rad, orthopedics and trauma applications

Maximum space saving

- Fast and easy installation

- Fast and easy patient positioning
Dose reduction and

excellent DQE

Fully DICOM 3 compliant and RIS integration workflow optimization



▶ Apelem Da Vinci Verso – multi-purpose single detector

Technology CsIResolution 143 µm Size 43 x 43 cm

Highlights

- Elevating four way mobile table with auto tracking capability
- Auto positioning
- Dose reduction and excellent DQE
- Single touchscreen console controls generator and ceiling suspension
 Full HIS/RIS and PACS integration

▶ Apelem LEO detector - portable flat panel detector

Technology CsIResolution 143 µm Size 35 x 43 cm



- Wireless/Wi-Fi connected flat panel detector
- Mobile and easy to handle: <5 kg
- Very high quality image: 144 µm pixel size
- Real-time imaging: 3 sec in preview
- Dose reduction thanks to the CsI scintillator




Light.

Light for true portability.

Canon's new range of Digital Radiography Imaging Systems.
Enhance workflow efficiency with our interchangeable flat panel detectors: the thinnest, lightest and most versatile in the industry.
See our latest innovation, the 3-in-1 detector for superior static, serial and fluoro imaging at www.canon-europe.com/medical

We Speak Image



Canon

▶ Canon CXDI-50G/50C

Technology 50 C: Cesium lodide Scintillator;

50 G: Gadolinium OxiSulfide scintillator

Resolution 160 µm Size 35 x 43 cm



▶ Highlights

- For table, upright and portable applications
- 35 x 43 cm imaging area
- 1 tile robust construction
- Powerful image processing software
- Click on grid

▶ Canon CXDI-40EG/40EC

Technology 40 EC: Cesium Iodide Scintillator;

40 EG Gadolinium OxiSulfide scintillator

Resolution 160 µm Size 43 x 43 cm



Highlights

- For table and upright applications
- -43 x 43 cm imaging area -1 tile robust construction
- Powerful image processing software
- No cooling required

▶ Canon CXDI-60G/C

Technology 60C: Cesium Iodide Scintillator

60G: a-silicon Gadolinium OxiSulfide Scintillator

Resolution 160 µm Size

23cm x 28cm

▶ Highlights

- Compact Portable Flat Panel Detector
- -2,5 kg weight
- X-ray image in just a few seconds
- User friendly touch screen operation
- Click on grid
- Powerfull image processing software

▶ Canon CXDI-55G/C

Technology 55C: Cesium Iodide Scintillator

55G: a-silicon Gadolinium OxiSulfide Scintillator

Resolution 160 µm Size 35cm x 43cm

Highlights

- For table, upright and portable applications 3.4 kg weight
- 1 tile robust construction
- -X-ray image in just a few seconds
- User friendly touch screen operation - Powerfull image processing software

▶ Canon CXDI-50RF

Technology Resolution Size Cesium Iodide Scintillator

160 µm

35cm x 43cm



▶ Highlights

- Portable Dynamic Flat Panel Detector
- Static, serial imaging and fluoroscopy
- Acquisition mode 5-30 fps
- Excellent signal to noise ratio
- Powerfull image processing software Available only in integrated applications

▶ Canon CXDI-40G Compact

Technology Resolution Size a-silicon Gadolinium OxiSulfide Scintillator


160 µm


43 x 43 cm




Highlights

- Retrofit system for existing bucky
- 1 tile robust construction
- Economical solution
- No cooling required
 User friendly touch screen operation
 Grid line elimination SW
- Powerfull image processing software









-240 shots/h

Outreaching image quality







▶ GE Healthcare Definium 6000

Technology Resolution Size a-Silicon

2022 x 2022 pixel, 14 bit

41 x 41 cm



▶ Highlights

Fast and proven detector technology
 More flexibility with mobile »flying« detector

- Fully motorized wall stand

OTŠ with vertical auto-tracking

Optional Advanced Applications
 Seamless digital workflow

- Flexible configurations, including cost-effective 1-detector shared solution

▶ GE Healthcare Discovery XR650

Technology Resolution Size a-Silicon

2022 x 2022 pixel, 14 bit



Highlights

High patient throughputFully motorized tube suspension

with Autopositioning

- Advanced Applications:
Autopasting, VolumeRAD, Dual Energy
- Flexible configurations, mobile or fixed detector

Seamless digital workflow

▶ GMM Opera Swing

Technology Resolution Amorphous silicon photodiodes array

148 um Size 43x43 cm



- Highly integrated all-in-one system ensuring enhanced exams

in digital RAD and Fluoro procedures Extraordinary user-friendliness combined with operational efficiency in any application (Emergency, DA, IR, digital Tomosynthesis, etc.)

Efficient execution of any exam either on the completely overhanging tabletop or in direct contact with the detector Easy and precise execution of lateral projections and oblique incidences also on stretchers

Intelligent user interface integrating all the controls of the system components in a unique exclusive Touch Screen combined with a series of joysticks

▶ GMM Opera RT20 "Guitar" and "Harp'

Technology Resolution Size Amorphous silicon photodiodes array

139 µm 43 x 43 cm



Highlights

Radiographic unit with double digital detector for application versatility and full operational efficiency

X-ray tube remarkable movement for the quickest and easiest execution of any exams and oblique incidences also stretchers

Accurate full-length examination of the patient with no need for repositioning

Total comfort for the patient and enhanced diagnostic results in any exam of the spine, thorax, legs, etc.

Ease of installation in any diagnostics room thanks to the extremely compact structure and extraordinary suppleness

▶ GMM Chorus

Technology

Amorphous silicon photodiodes array 139 µm

Resolution Size

43 x 43 cm



▶ Highlights

Integrated DR multifunctional system featured by an extremely flexible configuration and ultimate user-friendliness

Efficient, rapid execution of routine and specialized exams Advanced automatic alignment of the X-ray source to the detector movement on both examination table and wall stand

Examination table ensuring the utmost manoeuvrability and safe positioning

Rapid and accurate execution of any oblique incidence on stretchered patients thanks to the tilting and rotating wall stand

▶ Konica Minolta FlexDR C30

Technology Resolution

3.6 lp per mm Size 43 x 43 cm




Highlights

- User-friendly design provides maximum operability for both patient and technologist

High throughput for efficient workflow


Touchscreen preview display
 Acquisition station controls both DR and CR imaging


Connectible to different generators





- -Full diagnostic web based PACS + CR as an All in One
- Automatic image post processing gives superb quality for both soft tissue and bones
- Scans plates directly into the PACS
- Robust table-top digitizer
 Complete digital x-ray acquisition solution











▶ Mecall Eidos RF 439

Technology Resolution Size Amorphous silicon photodiodes array

148 µm 43 x 43 cm



Highlights

State-of-the-art system with single removable grid with exclusive auto-focusing device

Innovative completely overhanging carbon-fiber tabletop allowing examination of patient from any side

Adjustable height-tabletop ensuring an extraordinary minimum distance from the floor of only 50 cm


Ample variable Focal Distance

Full-length patient examination in both vertical and horizontal position with Stitching function possibility



User-friendly touch screen interface

DICOM services: print, store, query/retrieve, MPPS, WL





- Universal X-ray stand for all general radiology examinations - Swivel arm rotation: -45° - + 135°

Excellent price/performance ratio

User-friendly touch screen interface

DICOM services: print, store, query/retrieve, MPPS, WL





Professional image processing and integrated multimedia



radiographic positioning guide

▶ OR Technology Leonardo DR 4336R

Technology Resolution Size Gadox or CsI scintillator



▶ Highlights

Portable system for mobile X-ray imaging

- Provides all necessary components in one suitcase, including flat panel, notebook, *dicomPACS®DX-R* X-ray acquisition software and cables

▶ Philips DigitalDiagnost High-performance room - Dual detector



a-Si, CsI-Scintillator 3k x 3k image matrix, 143 µm pixel size 43 x 43 cm

▶ Highlights

-50 kW, 65 kW or 80 kW

- Easy switch from table to chest exams with two detectors

- Automated functions such as auto collimation and move-to-position

- Optional wireless portable detector and PCR integration

- Easy orthopedic imaging with automatic

image acquisition and stitching

Vertical stand (moveable or fixed) with integrated detector, digital bucky table with integrated detector and ceiling-based tube carrier

▶ Philips DigitalDiagnost Multi-purpose standard room - Single detector

Technology Resolution

a-Si, CsI-Scintillator 3k x 3k image matrix, 143 µm pixel size

Size

43 x 43 cm

▶ Highlights

–50 kW, 65 kW or 80 kW

- Versatile single detector room for medium to high patient load

- Extended move-to-position functionality
- Easy orthopedic imaging with automatic image acquisition and stitching

Moveable multi-purpose stand with swiveling c-arm and integrated detector, ceiling-based tube carrier and single side suspended table

Optional wireless portable detector and PCR integration

Optional swiveling table for better accessibility

See life more clearly

DR 200 Ma

Cost-Efficient Mobile DR System

With long-time imaging expertise and innovative research in medical engineering, Landwind is launching the most advanced digital mobile radiography system.



Features

- Optimized clinical workflow
- Outstanding operability
- Efficient APR function
- Fast bedside imaging
- Automatic image processing
- Fully Dicom compatible
- Superior DR image quality

▶ Philips DigitalDiagnost Compact room - Single detector



Resolution Size

a-Si, CsI-Scintillator 3k x 3k image matrix, 143 µm pixel size 43 x 43 cm

-50 kW, 65 kW or 80 kW

For multi-purpose use and medium workflow requirements
 Often used as chest room which also serves as back-up

general DR room

Tracking and move-to-position

- Fixed multi-purpose stand with swiveling c-arm and integrated detector, ceiling-based tube carrier plus height adjustable trolley Optional wireless portable detector and PCR integration

 Easy orthopedic imaging with automatic image acquisition and stitching

▶ Philips DigitalDiagnost with wireless portable detector

Technology Resolution Size

a-Si, CsI-Scintillator

3k x 2.4k image matrix, 144 μm pixel size

35 x 43 cm



Highlights

The wireless portable detector is available as an additional detector for all DigitalDiagnost single and dual detector configurations

– More flexibility: The wireless portable detector carries out even the most difficult projections projections at table, patient bed, wheelchair or trolley

More efficiency: Smooth digital workflow with instant results at the Eleva workspot

-More freedom: Convenient handling and high hygienic standards thanks to the wireless detector's cable-free design

Optional wireless portable detector and PCR integration



► **Highlights**- 50 kW, 65 kW or 80 kW

- Highly automated workflow with workstation controlled collimation, asymmetric beam alignment, automatic tracking

-Extended application range for skeletal examinations with tiltable vertical stand

-Optional wireless portable detector and PCR integration

▶ Philips Essenta DR

Technology Resolution

a-Si, CsJ-Scintillator

3k x 3k image matrix, 143 µm pixel size





Highlights

–50 kW, 65 kW or 80 kW

- Cost-effective flat detector technology for state-of-the-art direct digital imaging

- Easy handling through motorized movements

- Tiltable detector and rotatable tube for unlimited patient positioning incl. free cassette exposures

Floor-mounted, u-arm based

▶ Philips Essenta DR Compact

Technology Resolution

Quantmaster, High Stability Scintillator (GoS)

Matrix: 1920 x 2367



Highlights

– 50 kW, 65 kW or 80 kW

- affordable price with low cost of ownership

- Affordable handling through counterbalanced movements - Fits into very small rooms (height 2,50 m)

Floor-mounted, u-arm based

▶ Philips DRF room solution

Technology Resolution a-Si, CsI-Scintillator

3k x 3k image matrix, 143 µm pixel size



▶ Highlights

-Maximized room utilization with high-quality digital radiography and fluoroscopy applications in just one room

- Filmless workflow with DR technology for high throughput
- Excellent image quality with UNIQUE image processing and

DoseWise concept
One room for all types of patients from infants to obese adults

▶ Protec Rapixx 43 WiFi DR-System

Technology CsIResolution 144 µm Size 36 x 43 cm



- 16 bit dynamic range
- Wireless system connection (WiFi)
- Portable and easy to handle: 4,8 kg

- Images in 3 sec. result in high productivity
 Versatile and robust design for long lifespan
 Simple integration and upgrade into existing conventional X-ray units
- Outstanding flexibility: close at hand, close at patients, just
- one panel required for bucky table and wall integration

 Docking station, interface box, power supply and CONAXX image acquisition software included in standard delivery
- Fully DICOM compatible for integration to PACS

▶ Protec Rapixx 4336M DR-System

Gd₂O₂ or CsI Technology Resolution 139 µm Size 36 x 43 cm



Highlights

- 16 bit dynamic range
- Cable connection
- -Lightweight: 3,6 kg
- Minimal cycle time: 8 sec.
- Predestined for simplest retrofitting of existing X-ray units
- due to dimensions equal to conventional X-ray cassette -High shock tolerance and water resistant portable flatpanel
- detector - Interface box, power supply and CONAXX image acquisition software included in standard delivery
- Fully DICOM compatible for integration to PACS

▶ Protec Rapixx 4343M DR-System

Technology Gd₂O₂ or CsI Resolution 139 µm Size 43 x 43 cm



RAPIXX

▶ Highlights

- 16 bit dynamic range
- Cable connection
- -Weight: 7,5 kg
- Minimal cycle time: 6 sec.
- -For integration and upgrade into existing conventional X-ray
- units / intended for constant mounting in a X-ray unit

 Interface box, power supply and CONAXX image acquisition software included in standard delivery
- Fully DICOM compatible for integration to PACS

▶ Provotec PEDS 600

Design Floor-wall mounted Table Prognost XP series, optional 50 kW Power



▶ Highlights

- DR-System with dig. flat panel detector Variable SID 120 180 cm
- Rotatable u-arm 360°
- Rotatable DR-detector

▶ Shimadzu RADSpeed Safire

Technology

Direct-conversion flat-panel detector (a-Se)

Resolution 3.3 lp/mm

Size

17" x 17" (43 x 43 cm)



▶ Highlights

- Superior image quality
- 2880 x 2880 pixel; 16 bit dynamic range
- Parameter setting next to the patient Up to 400 application
- programs Auto-positioning


function

- Automatic tracking functions

▶ Siemens Axiom Aristos FX Plus

Technology Amorphous-Silicon with Cesium Iodide scintillator Detector 143 µm, 3k x 3k, 14 bit

Size 43 x 43 cm



▶ Highlights

- Universal digital flat detector solution
- Fully-automated system positioning via organ programs
- Tube and detector independently mounted on the ceiling
 Auto tracking of X-ray tube and detector
- - in x-, y- and z-direction
- Automated ortho acquisition of entire spine and long legs
- Excellent detail contrast with DiamondView

Siemens Axiom Aristos MX

Technology Detector Size

Amorphous-Silicon with Cesium Iodide scintillator 143 µm, 3k x 3k, 14 bit

43 x 43 cm



▶ Highlights

- Universal digital flat detector solution

 Auto tracking of X-ray tube and detector in all directions

- Motorized grid removal via organ programs

- Pull-out detector for ease of patient positioning

Excellent detail contrast with DiamondView

▶ Siemens Axiom Aristos VX and VX Plus

Technology Detector Size Amorphous-Silicon with Cesium Iodide scintillator 143 µm, 3k x 3k, 14 bit

43 x 43 cm



► **Highlights** – Digital flat detector solutions for chest and skeletal applications

Auto tracking in the vertical direction
 Comprehensive control of parameters via organ programs
 TOP alignment of X-ray field for dose reduction

- Automated ortho acquisition of entire spine and long legs

▶ Siemens Axiom Aristos TX

Technology Detector Size Amorphous-Silicon with Cesium Iodide scintillator 143 μm, 3k x 3k, 14 bit 43 x 43 cm



Highlights

- Dedicated flat detector solution for chest imaging

Automated workflow processes for high patient throughput

- Comprehensive control of parameters via organ programs

- TOP alignment of X-ray field for dose reduction

Excellent detail contrast with DiamondView

▶ Siemens Axiom Multix M

Technology

Amorphous-Silicon with Gadolinium Oxysulfide scintillator

Detector 160 µm,

2688 x 2208 pixel

Size 35 x 43 cm



Highlights

- Universal digital radiography solution with mobile flat detector

-Flexible and easy handling – positions just like a cassette

Electronic tomography possible (option)
 Ceiling-mounted and floor-mounted solutions available

▶ Siemens Axiom Vertix MD Trauma

Technology Detector

Size

Amorphous-Silicon with Cesium Iodide scintillator 160 µm,

2688 x 2208 pixel 35 x 43 cm



Highlights

Digital radiography solution with mobile flat detector

- Ceiling-mounted u-arm for maximal flexibility

X-ray tube is constantly centered to flat detector in all planes
 All exposures with one detector, in or out of the holder

Fast image preview available within 5 seconds

▶ Siemens Ysio

Technology Detector Amorphous-Silicon with Cesium Iodide scintillator Fixed detector (148 µm), 14 bit

Wireless detector, wi-D with (144 µm), 16 bit

Size 43 x 43 cm

35 x 43 cm



▶ Highlights

Flat detector, digital radiography (DR) solutions
 Choice between fully automated or fully synchronised systems

- Digital flat detectors with newest detector technology

 Ceiling-mounted tube with MaxTouch – a color touchscreen for enhanced workflow

Automated system positioning and synchronised tracking of X-ray tube and detector in different planes

-Excellent detail contrast with DiamondView Plus

▶ Siemens Axiom Luminos dRF

Design Technology Resolution II-format Image system Size Remote-controlled 2-in-1 system with dynamic flat detector Amorphous-Silicon with Cesium Iodide scintillator

Up to 3.4 lp/mm

43 x 43 cm



▶ Highlights

- Fully digital 2-in-1 solution for dynamic and static high-resolution imaging including DSA procedures (option)

– Easy patient transfer at 48 cm lowest table height

- Dynamic Density Optimization (DDO) and DiamondView Plus

- FLUOROSPOT Compact high-resolution digital imaging system with intuitive user interface and DICOM 3.0 interfaces

- Comprehensive CARE dose reduction package

- Limitless projection flexibility with optional ceiling-suspended tube and wireless detector wi-D

▶ Siemens Multix Swing mFD

Technology Detector Size Amorphous-Silicon with Oxysulfide scintillator 160 µm 2688 x 2208 pixel

35 x 43 cm



▶ Highlights

- Cost-efficient, all-in-one DR solution with mobile Flat Detector

Flexible positioning of mobile detector in table,

wall stand and for free exposures

Generator is integrated into the table for minimal space requirements
Accomodates wide range of exams for cost-conscious

digital imaging

Synchronized tube and bucky tray movements

▶ Tetenal Vidix U (Universal Type)

Technology Resolution Size a-Se (amorphous Selenium) 168 µm, 14 bit, 2.560 x 2.560 Pixel 43 cm x 43 cm



▶ Highlights

- User-friendly Easy Reliable Efficient
- Low radiation
- Full auto positioning
- Interchangeable grid
- Anti-collision system

▶ Tetenal Vidix II (Single Type)

Technology Resolution Size a-Se (amorphous Selenium) 168 µm, 14 bit, 2.560 x 2.560 pixels 43 cm x 43 cm (17" x 17")



▶ Highlights

- Low radiation
- Preview function
- Automatic tracking function
- Immediately ready without preliminary lead time
 Entry area for bucky table is very low

▶ Tetenal Vidix II (Dual Type)

Technology Resolution Size a-Se (amorphous Selenium) 168 μm, 14 bit, 2.560 x 2.560 pixels 43 cm x 43 cm (17" x 17")



- Improved workflow and working time (fast image acquisition)
- High image quality (excellent dynamic range)
- Decrease of x-ray dose
- Cost-efficient Low radiation

▶ Highlights

▶ Valmex X Store DR

Technology

Selenium type direct conversion detector or Cesium Iodide-Scintillator detector From 129 µm up to 168 µm pixel size

Resolution

20 x 25 cm, 35 x 43 cm, 43 x 43 cm



▶ Highlights

– Plug&Ray® Bucky

System for existing x-ray units

Slim form factor – portable or stationary type
High frame rate with up to 225 exposures per hour

- Dynamic bandwidth - one exposure for all purposes

Moving pendulum grid to exclude artificial grid suppression for better diagnosis

X Store[®] Image Viewer acquisition Software X Store[®] Generator Software

X Store® PACS Software

CR/DR MOBILE



▶ Highlights

16bit dynamic range; Up to 3.7lp/mm

- Fixed multi-purpose stand with swiveling c-arm and Ceiling suspension

- Best image qualityComplete DICOM Data
- -50KŴ,65KW or 80KW

▶ Fujifilm Go


15 kW Power kV Range 40 - 130mAs Range 0.5 - 200 mAs




▶ Highlights

- Flexible & high-performance portable digital X-ray unit

- Dual motor drive allows free and smooth steering
 Easy positioning with telescopic arm
 Smart and easy to use with additional controls at the tube head
- -For use in every department, even NICU and OR










▶ Siemens Mobilett XP Digital

Detector Resolution Power kV Range mAs Range 2688 x 2208 pixel

30 kW 40 – 133



▶ IBA Dosimetry Primus

Test device for checking image quality parameters at fluoroscopic units



▶ Highlights

 Modular construction: structural plate and separated attenuator

- Check of spatial and contrast resolution, size of the radiation field, artefacts; kV test area
- Compact Al pre-attenuator or PMMA and Cu plates
- Available in two different sizes

- Direct organ selection program
- Fully DICOM compatible

- WLAN connectivity for improved workflow

▶ IBA Dosimetry DIGI-13

Test device for checking image quality parameters at digital radiographic units



▶ Highlights

- Compact device with separated Al pre-attenuator
- With integrated copper plate
- Check of homogeneity, spatial and contrast resolution, size of the radiation field, artefacts
- Easy-to-use

▶ IBA Multimeter MagicMax

Simultaneous measurement of dose, dose rate, exposure time, kV, dose/pulse, pulse rate, HVL and total filtration



- ► Highlights

 Small device with separate multifunction detector

 Connected via USB to PC or Notebook

 Intuitive use via PC interface

- Time resolution: 100 µs
- Optimized solutions for all applications

▶ medigration CD-Imager



▶ Provotec Droc

Design Table Power

Placed on a desk Desktop Line connected





- DROC Digital Radiographic Operator Console For DR-detectors and CR-units
- Automatic image optimization and instantaneous acquisition time
- Advanced noise reduction and optimized image calibration
- technology
 Full DICOM 3.0 compatibility
- Support DR-detector and generator
- Remote online system diagnosis

▶ PTW Normi 13

Test object for quality control of digital radiography X-ray units



- Checks all imaging quality parameters (dynamic range, spatial resolution, low contrast, artefacts, radiation field, etc.)
- Convenient use at bucky units
- Patient equivalent absorber (Al or PMMA) included

▶ Radcal Accu-Pro™

X-Ray Analyzer

Simultaneous dose, rate, time, kVp. HVL, filtration, mA/mAs,



▶ Highlights

- Use for manufacturing, installation, QA, and service
 R/F, mammography, CT, dental, leakage
 Ion chamber based dosimetry, no corrections required
- Correctly measure AEC fluoro and filtered beams Remote control, waveforms, and archiving with XLPRO software
- Compact, easy to use

▶ Radcal Rapidose

PC X-ray Analyzer



Highlights

- Plug into a laptop USB port for an inexpensive X-ray analyzer
 Simultaneous dose, rate, kVp, time, HVL, waveforms, and
- Revolutionary inherent remote measurement operation
- Easy use, genuine time saver
- Data archiving and analysis using your Excel

▶ Radcal PDC-DAP/KAP verification meter



Highlights

- Newly patented Patient Dose Calibrator
 Use to calibrate DAP/KAP and rate
- Also measures dose and rate
- Optical and radiographic alignment markers
- Simple to use with optional computer control

▶ RTI Electronics Piranha

The Piranha is designed as a truly self-contained, all-in-one, X-ray multi-function meter that assures accurate results in one shot. kV, time, dose, dose rate, HVL and total filtration.



▶ RTI Electronics Barracuada

The Barracuda X-ray multimeter has a cabinet that can house up to six different application modules, and can measure on all modalities; R/F, mammography, flouroscopy, pulsed flouroscopy, dental, panoramic dental and CT systems.



- ► Highlights

 All in One, All at Once
- Auto-Compensation
 Enhanced Graphical PDA Display
 R&F, Mammo, Dental and CT
- Ionization Chambers
- Built-In Bluetooth for PC and PDA
- -mAs, and Light Probes
- Fits in the Palm of Your Hand

RAD-BOOK 2010

-mA, mAs, and Light Probes

- Fits in the Palm of Your Hand



114 MOLECULAR IMAGING

GE Healthcare Infinia 3.8 mm intrinsic FWHM Resolution Sensitivity 270 cpm/μCi (LEGP) Field of View UFOV: 540 x 400 mm

▶ Highlights

- High performance variable dual head system
- Unmatched productivity time saving up to 15 %
 Excellent clinical versatility & unlimited flexibility
- Advanced image quality
- High reliability and excellent serviceability

▶ GE Healthcare Infinia Hawkeye 4

Resolution Sensitivity Field of View

3.8 mm intrinsic FWHM 270 cpm/μCi (LEGP) UFOV: 540 x 400 mm



▶ Highlights


- True integrated hybrid imaging system
- Four slice axial/helical CT scanning
 Superior image quality & flexibility
 Ultra low dose CT technology

- Leading economic value

▶ GE Healthcare Millennium MG

Resolution Sensitivity Field of View

3.9 mm intrinsic FWHM 262 cpm/µCi (LEGP) UFOV: 510 x 360 mm



▶ Highlights

- -Variable multi-purpose dual head system
- Unmatched scanning efficiency
- Excellent clinical versatility all applications
- Extreme small footprint & gantry size
- High reliability and excellent serviceability

▶ GE Healthcare Millennium MPR

Resolution Sensitivity Field of View

3.9 mm intrinsic FWHM 262 cpm/µCi (LEGP) UFOV: 510 x 360 mm



▶ Highlights

- Multi-purpose single head system
- Operational ease and flexibility
- Open gantry extra large detector
 Excellent clinical versatility all applications
- High reliability and excellent serviceability

▶ GE Healthcare Ventri

Resolution Sensitivity

3.7 mm intrinsic FWHM 325 cpm/μCi (LEGP)



▶ Highlights

- Fixed angle dual head cardio system
 Designed for all patient sizes
 Uncompromized patient comfort

- Superior image qualityExtreme small footprint & gantry size

▶ GE Healthcare Discovery VCT XT

Resolution @ 1cm 5.0 mm NECR @ 6 kBq/cc 64 kcps Transaxial PET Field of View 70 cm



- Volume CT inside
- Uncompromized between high sensitivity and resolution
- -5-beat cardiac CT angio -Snapshot Pulse 70% dose reduction for CT angio
- -VUE point HD 3D iterative reconstruction

▶ GE Healthcare Discovery ST

Resolution @ 1cm 6.1 mm NECR @ 6 kBq/cc 57 kcps Transaxial PET Field of View 70 cm



▶ Highlights

- Optimized for 3D and 2D imaging
 Dynamic and gated PET, CT and PET/CT acquisition
 MotionFree technology
 VUE point HD 3D iterative reconstruction

- Designed for maximum flexibility

▶ GE Healthcare Discovery STE

Resolution @ 1cm $5.0 \; \mathrm{mm}$ NECR @ 6 kBq/cc 64 kcps Transaxial PET Field of View 70 cm



▶ Highlights

- Uncompromized between high sensitivity and resolution
- -Highest NECR in clinical range*
- Optimized for 3D and 2D imaging VUE point HD 3D iterative reconstruction
- MotionFree technology (* based on NEMA 2001)

▶ GE Healthcare Discovery 600

Resolution Peak NECR Transaxial PET Field of View

5.0 mm

75 kcps @ 15 kBq/ml

70 cm




- Highlights
- Motion Free PET/CT
- Optimized for gated and dynamic PET- and CT-acquisitions
 VUE Point HD 3D iterative reconstruction
- IBM Blade Center for fastest reconstruction
- Highest NECR in clinical range

▶ GE Healthcare Discovery 690

Resolution 4.9 mm

Peak NECR 110 kcps @ 20 kBq/ml

Transaxial PET Field of View 70 cm



- -VUE Point HD 3D iterative reconstruction
- with time of flight information
- IBM Blade Center for fastest reconstruction

▶ GE Healthcare Discovery NM/CT 670

Resolution Sensitivity Field of View

3.7 mm intrinsic FWHM ** 270 cpm/uCi (LEGP) *** UFOV 540 mm x 400 mm



- ▶ Highlights
- Combination of SPECT technology power with proven multi-slice CT
- advanced fusion imaging
 State-of-the-art Nuclear Medicine slim detectors
- Ultra-fast gantry robotics
 BrightSpeed™ Elite 16 slice CT


*FDA 510k cleared; EC declaration of conformity pending. ** preliminary data

► GE Healthcare Discovery NM 530c with Alcyone™ Technology

SPECT Acquisition Time 2 min SPECT Resolution 3.1 mm tangential



- Alcyone Technology
 four leading-edge technologies in nuclear imaging: CZT detectors, focused pinhole collimation, 3D reconstruction, and stationary data acquisition.
- -increased diagnostic confidence and efficiency
- improved patient care with low dose and fast scans
 upgradeable
- no motion



-Compact design fits easily into a 2.4 x 3 m room Fixed-90 dual head design and cardiac workflow are

optimized for high throughput

- Comprehensive suite of cardiac image applications

- Vantage Pro clinically validated non-uniform attenuation correction

▶ Philips Forte with JETStream acquisition and AZ detectors 3.3 mm, FWHM intrinsic Resolution

Sensitivity 265 cpm/µm Ci (LEGP) 38.1 x 50.8 cm



Highlights

- High throughput capability for all nuclear medicine procedures

- VantagePro clinically validated non-uniform attenuation

Automatic acquisition workflow

- Concurrent imaging allows a single acquisition step to be saved simultaneously in up to 15 unique data sets

Energy independent detectors provide high performance with reduced quality control time

Philips Gemini TF PET/CT - GXL PET/CT

TF PET/CT Peak NECR 210 kcps **PET Spatial Resolution** 4.3 mm 16-slice or 64-channel CT Configuration

► Highlights TF - World's first

commercially available time-of-flight PET/CT

Fast scans (10 min) with low dose

- Premium Brilliance CT image quality & applications 190 cm PET/CT scan

length

- Exclusive OpenView gantry design

PET Spatial Resolution

Peak NECR

CT Configuration

GXL PET/CT

70 kcps

4.5 mm

16-slice

► Highlights GXL

- Fast scans (15 min) with low dose

- Fully 3D LOR PET reconstruction

- Premium Brilliance CT image

quality & applications
- 190 cm PET/CT scan length

- Exclusive OpenView gantry design

▶ Philips TF Big Bore PET/CT



Radiation Oncology

- Brilliance CT Big Bore
subsystem and exclusive PET TruFlight technology

- 85cm bore diameter to accommodate all positioning devices

for radiotherapy planning.

– Compliant with the AAPM TG-66 standards

for positional accuracy

-State of the art diagnostic image quality

▶ Philips BrightView

Resolution 3.3 mm, FWHM intrinsic Sensitivity 277 cpm/µm Ci (LEGP) Field of view 40.6 x 54 cm



Highlights

-Patient focus for an open experience with all patients and sizes

Maximized image quality with CloseUp technologies

- Improved workflow efficiency, BodyGuard automatic contouring

Rich in capability yet compact in design

- Scalable to match the capabilities with practice

▶ Philips SKYLight

3.3 mm, Resolution

FWHM intrinsic

Sensitivity 265 cpm/µm Ci (LEGP) Field of view 38.1 x 50.8 cm



▶ Highlights

Gantry-free design for clinical flexibility and openness

 DualPlanar capability to acquire two patients simultaneously

Automatic acquisition workflow

- Concurrent imaging allows a single acquisition step to be saved simultaneously in up to 15 unique data sets

- Fully automatic collimator exchange

Philips Precedence SPECT/CT

3.3 mm, FWHM intrinsic Resolution Sensitivity 265 cpm/umCi (LEGP) Field of view 38.1 x 50.8 cm



▶ Highlights Acquisition efficiency

Reconstruction leadership

- Ease-of-use

System efficiency

Superior diagnostic CT image quality

▶ Siemens c.cam

Resolution ≤3.7 mm FWHM in UFOV Sensitivity

290 cpm/μCi (LEAP at 10 cm at 140 keV) Field of view



▶ Highlights

A whole New Angle in Cardiology

Award winning design
Open and patient-friendly

- Reclining position improves image quality

-Easy to install, learn and use

- Comprehensive quantification and viewing software

- Field upgradeable with c.clear attenuation correction

▶ Siemens Symbia E

Resolution ≤3.8 mm FWHM in CFOV Sensitivity 202 cpm/μCi (LEHR 3/8"at 10 cm)

Field of view 53.3 x 38.7 cm



Highlights

- Simply Evolutionary

- Work with confidence through Siemens HD detector, c.clear attenuation correction and Siemens Remote Service

Accelerate your workflow through fast acquistion, autocontour, syngo workflows and physician worklist

Experience versatility through detector flexibility, clinical engines, specialized pallets and investment protection

▶ Siemens Symbia S

Resolution Sensitivity Field of view

≤3.8 mm FWHM in CFOV

202 cpm/μCi (LEHR 3/8" at 10 cm)



▶ Highlights

- The first SPECT that's ready for everything
- Work with confidence through Symbia HD detector, superb image quality with Flash and Siemens Remote Service

Accelerate your workflow through fast acquisition, autocontour, syngo workflows, Integrated/Automated Collimator Changer and Automated Quality Control
- Experience versatility through detector flexibility, clinical

engines, investment protection and bariatric imaging

▶ Siemens Symbia TruePoint SPECT • CT

Resolution Sensitivity ≤3.8 mm FWHM in CFOV

Field of view



Highlights

- A Merger of True Equals

- Work with confidence through Symbia HD detector, superb image quality with Flash, Siemens Remote Service and


SPECT•CT University
- Accelerate your workflow through fast acquisition, autocontour, CT attenuation correction, Integrated/Automated Collimator Changer and Automated Quality Control
- Experience versatility through detector flexibility, clinical

engines and diagnostic CT (2, 6, 16 slice)

▶ Siemens Biograph TruePoint PET•CT

Resolution Sensitivity Field of view

2.0 mm average FWHM at 1 cm (HD•PET) 7.6 cps/kBq at 435 keV (TrueV) 605 mm transaxial, 216 mm axial (TrueV)




▶ Highlights

- -The World's First HD•PET Platform
- Uniform resolution throughout the field of view with HD•PET
- Largest PET field of view in the industry, increasing count rate performance by >70%
- Exceptional lesion detectability with the best NEMA spatial resolution in the industry

Siemens Biograph mCT

Resolution Sensitivity Field of view 2.0 mm average FWHM at 1 cm (HD•PET) 9.5 cps/kBq at 435 keV (TrueV) 700 mm transaxial, 216 mm axial (TrueV)



▶ Highlights

- The World's First Molecular CT 5 minute whole-body PET•CT scans
- unprecedented throughput for true shared service capability in CT and PET
 Exceptional image quality
 More confortable examinations

- Low patient doses

MI ACCESSORIES

▶ Alliance Medical – flexible diagnostic imaging services



▶ Highlights

- Static diagnostic imaging centers MRI, CT, PET, PET/CT
 Interim services for bridging downtimes
- Regular "routing" services

▶ IBA Dosimetry Dosimax plus A HV

Dosimeter for measuring simultaneously dose, dose rate, exposure time and dose length product



Highlights

- Designed according to IEC 61674
- -For use with solid state detectors or ionization chambers
- For CDTI determination in combination with head and body phantom

▶ PTW CT Dosimetry

Quality control equipment for CT dose measurements



Highlights

- Combined head & body phantom available
 CT chamber for precise dose length product measurements and CTDI determination

▶ RTI Electronics CT Slice Probe



The CT Slice Probe is designed to make CTDI measurement more exact and has also the ability to further analyze the result. Following parameters are achived from a single exposure: CTDI100, CTDIvol, DCTIw, CT dose profile, DLP, Point Dose, Performance of the AEC, FWHM and Scatter Index.




- All in One Shot
- Quick and Simple Set up
- Accurate and Sensitive
- No limitations due to the beam width

MOLECULAR IMAGING

▶ TeraRecon Aquarius

Product Technology Resolution Size

Aquarius workstation 3D diagnostic workstation AquariusNET server Client-server 3D architecture



▶ Highlights

- Thin client-server solution enterprise-wide
- VolumePro: handling many 3D sessions at once
 Rendering on central server, results streamed to PC
- Fast and efficient in the reading workflow
- Aquarius NET runs on almost any standard PC

▶ Tomovation - Modular building solutions



- Engineering, rental, sale of modular buildings MRI, CT, PET, PET/CT including or excluding diagnostic equipment

▶ Ziosoft Ziostation

Full Suite of Clinical Applications

3D, 4D, MIP, MPR, CPR, Multi-Modality Viewer

▶ Highlights

- Multi-volume fusion
- Vessel analysis
- Coronary analysis
- Colon analysis
- PET/CT viewer
- CT Cardiac Function Analysis
- CT Brain Perfusion
- Calcium Scoring
- Reporting functions
 Customizable button palettes and keyboard/mouse shortcuts to configure the system individually
- Automatic segmentation algorithms employ anatomical analysis, targeting specific objects to maximize accuracy

 – Full Web-Based Capabilities
- Action macros perform multiple processes with one or zero button clicks

Alliance Medical

Come visit us at ECR, Expo E, 532



Flexible diagnostic imaging services

Static diagnostic imaging centers without financial risk: MRI, CT, PET/CT



Interim services for bridging downtimes or reducing waitlists: MRI, CT, Cardiac CathLab, PET, PET/CT



Regular services make new technologies available for your patients: MRI, CT, PET/CT, Marnmography

vertrieb@alliancemedical.eu www.alliancemedical.de

We scan for you - We scan for you


TOMOVATION®



Innovative products & services MRI, CT, PET, radiation therapy

tomoinject MRI & CT injectors





tomomodular rental and sale of modular buildings MRI, CT, PET, PET/CT


info@tomovation.com Alliance Medical www.tomovation.com

We scan for you - We scan for you



RAD-BOOK 2010





▶ Barco Coronis Fusion 6 MP DL Wide-screen diagnostic color display system -Name Coronis Fusion 6 MP DL

Technology Resolution Size

Color LCD 6 Megapixel (3280 x 2048)


30 inch



▶ Highlights

Bezel-free 30 inch desktop for multi-modality PACS imaging

- Unmatched viewing characteristics High-performance image processing
- Automated Quality Assurance
- -5-year warranty



- 5-year warranty





► Highlights


- Providing consistent DICOM images anywhere, anytime

- Professional LCD Quality

- Approved for medical use

- Backlight output stabilization

- User-friendly Quality Assurance





▶ Image Displays cQ-2MP BASIC

20.1 inch Screensize MegaPixel 2 MP Resolution

1.600 x 1.200 Pixels



▶ Highlights

- Optimized for orthopedic surgeons and imaging centers
- Multiple video interfaces
- Picture-in-picture
- Perfect combination of high performance and affordable
- Available with color, grayscale and extra-bright grayscale panels

▶ Image Displays cQ-3MP PRO

Screensize 21.3 inch MegaPixel 3 MP

Resolution 2.048 x 1.536 Pixels



▶ Highlights

- High-res medical diagnostic displays for professional users
 First medical displays worldwide including true 12-bit image grayscale transfer in the processing line
- Easy and fast installation; IMAGE DISPLAYS are calibrated
- twice before delivery

 Intelligent self-diagnosis function for precise luminance and DICOM calibration
- Excellent image quality even in light rooms

▶ Image Displays cQ-5MP PRO



Screensize 20.1 inch MegaPixel 5 MP Resolution

2.048 x 2.560 Pixels

▶ Image Displays cQ-2in1

Screensize 30 inch MegaPixel 4 MP

1.600 x 2.560 Pixels Resolution



▶ Highlights

- High-res medical diagnostic display for professional users
- Designed for Mammography readings
- Excellent image quality even in light rooms
 First medical displays worldwide including true 12-bit image grayscale transfer in the processing line
- Intelligent self-diagnosis function for precise luminance and DICOM calibration

▶ Highlights

Revolutionary solution for the simultaneous assessment of two X-ray images on one 30" panel

Outstanding magnified views with full quality

- No frameworks in-between when doing dual readings or viewing large images
- Homogenous color and grayscale across the whole display
- Highest performance in processing radiological and nuclear medical images

▶ NDS surgical imaging E-Series

1 MP 2 MP GS E2 3 MP GS E3

4 MP GS E5 5 MP



▶ Highlights

- -FDA 510(k) cleared
- Automatic DICOM calibration
- -3061 unique shades of grey
- Fanless display
- High operating brightness

▶ NDS surgical imaging Ec-Series

1 MP 2 MP E2colour

3 MP E3colour 4 MP E4colour

5 MP



- -FDA 510(k) cleared
- Automatic DICOM calibration
- Fanless display
- High operating brightness
- Flexibility of colour and DICOM calibrated greyscale imaging

▶ NDS surgical imaging Dome Dashboard



▶ Highlights

- Simplifies network display management in healthcare environment
- Monitor and maintenance console for multiple workstations
- Supporting Dome medical imaging displays and the Dome CXtra software
- Operates within
 Windows framework
- Provides intuitive features and functions

▶ NDS surgical imaging Dome GX 2MP

Resolution MP Color Size 20,1"



▶ Highlights

- Superior visual quality with .255mm pixel pitch
- 16 ms response time Outstanding 1000:1 contrast ratio
- 300 cd/m²luminance
- Dome CXtraTM for DICOM Calibration

► NEC MD GS-Series



MD21GS-2MP

21.3" (2.0 MP image resolution: 1200 x 1600 pixels in portrait orientation)

- Unique XLight® Backlight System for long lasting stable white point
- Individual factory calibrated DICOM gamma correction curve
- Precise DICOM GSDF calibration using internal 10-bit look-up table

21.3" (3.0 MP image resolution: 1536 x 2048 pixels in portrait orientation)

▶ Highlights


- Unique Frontsensor for long lasting brightness stability
 Individual factory calibrated DICOM gamma correction curve
- Precise DICOM GSDF calibration using internal 12-bit look-up table

20.1" SA-SuperFine Grayscale TFT glass

▶ Highlights

- -5.0 MP image resolution: 2048 x 2560 pixels in portrait orientation
- Up to 1024 simultaneous shades of gray out of a pallette of 3061 (10-bit image reproduction)
- High adjustable stand (Range: 120 mm)

NEC MD Color-Series



MD212MC

21.3" (2.0 MP image resolution: 1200 x 1600 pixels in portrait orientation)

- Unique X-Light® 3 Backlight System for long lasting stable white point
 Individual factory calibrated DICOM gamma correction curve
- Pre cise DICOM GSDF calibration using internal 12-bit look-up table

MD213MC

21.3" (3.0 MP image resolution: 2048 x 1560 pixels in portrait orientation)

▶ Highlights

- Unique XLight® 3 Backlight System for long lasting stable white point
- Individual factory calibrated DICOM gamma correction curve
- Precise DICOM GSDF calibration using internal 12-bit look-up table

MD304MC

29.8" (4.0 MP image resolution: 2560 x 1600 pixels)

▶ Highlights

- High performance H-IPS TFT display
- Stand-alone calibration and matching to ensure efficient onsite colour temperature matching, with 12-bit LUT and precise DICOM GSDF grayscale
- Medical certification conforming to 93/42/EC (MDD) and FDA510k pending

MD21M

21.3" (2.0 MP image resolution: 1600 x 1200 pixels)

- Factory pre-set DICOM gamma curve
- Digital Uniformity Control (DUC) for enhanced spatial uniformity
- Complies with DIN V6868-578 Cat.B and AAPM TG18 Secondary Class

▶ NEC MDview Series Colour Displays

MDview193

19.0" (1.3 MP image resolution: 1280 x 1024 pixels)

▶ Highlights

- Real Hardware-DICOM Calibration including Ambient Light compensation
- Quick and easy setup utilizing the factory pre-set DICOM curve
- Long distances to modalities and workstations without compromise in display quality font size is to small compared to the others

MDview 202

20.1" (2.0 MP image resolution: 1600 x 1200 pixels)

▶ Highlights

- Factory pre-set DICOM gamma curve
- Digital Uniformity Control (DUC) for enhanced spatial uniformity
- Precise DICOM GSDF calibration using 3 internal 12-bit look-up tables

21.3" (2.0 MP image resolution: 1600 x 1200 pixels)

▶ Highlights

- Factory pre-set DICOM gamma curve
- Digital Uniformity Control (DUC) for enhanced spatial uniformity
- Precise DICOM GSDF calibration using 3 internal 12-bit look-up tables



MDview 242-2

24.1" (21.3 MP image resolution: 1920 x 1200 pixels)

Highlights

- Factory pre-set DICOM gamma curve
- Digital Uniformity Control (DUC) for enhanced spatial uniformity
- Precise DICOM GSDF calibration using 3 internal 12-bit look-up tables

MDview 262-2

26" (2.3 MP image resolution: 1920 x 1200 pixels)

▶ Highlights

- XLight® Pro: Backlight system for stabilised white luminance and
- Factory pre-set DICOM gamma curve
- Digital Uniformity Control (DUC) for enhanced spatial uniformity

▶ Provotec ProVario Screen

Design Power unit floor-mounted Table Control console placed on a desk Power 50 kW



▶ Highlights

- High frequency generator for x-ray diagnostic
 Easy operation by monitor or touchscreen
 Digital control of nearly unlimited organ programs
- Safety device against undue radiation for each organ with AEC-technique
- -X-ray book for storing patient name with generator exposure data
- Upgradeable for using CR- and DR-systems

▶ Totoku Mammography Display

Grayscale 15 MP MS51i2



▶ Highlights

- ISD Technologie to

support »Super high resolution« of 15 megapixel Supports 5 megapixel videoboards thanks to ISD technologie

- Luminance uniformity correction
- Remote calibration and management
- New developed AR coating for an improved resolution
 Increased lifetime with backlight dimming system

▶ Totoku Diagnostic Displays

Gravscale Color 5 MP ME551i2 3 MP ME351i2 CCL354i2, CCL352i2 ME253i2, ME201L CCL240, CCL254i2, CCL252i2, 2 MP CCL208 1 MP ME181L CCL182



Highlights

- Up to 5-years warrantyHigh brightness color displays up to 800 cd/m²
- Luminance uniformity correction
- Remote calibration and management
- Up to true 11 bit grayscale
- Increased lifetime with backlight dimming system

▶ Totoku ISD Displays

Grayscale 3/9 MP MS31i2 2/6 MP MS21i2



- ISD Technologie
- Luminance uniformity correction
- Remote calibration and management
- New developed AR coating for an improved resolution
- Increased lifetime with backlight dimming system

▶ Agfa HealthCare Drystar 5503

Capacity Resolution System

100 films/h (14 x 17") 508 dpi/50 µm pixelsize Direct Digital Imaging



film sorter -Ideal for centralized workflow, can easily be connected

to the network

- Integrated A#Sharp technology for optimized image quality

- 3 multi-format trays, each supporting different film sizes and types

Suitable for CT, MRI, DSA, digital R&F, CR, DR and optional mammography applications

▶ Agfa HealthCare Drystar 5302

75 films/h (14 x 17") Capacity Resolution 320 dpi

System Direct Digital Imaging



▶ Highlights

- Suitable for all applications and ideal for CR/DR
- A#Sharp technology for optimized image quality
 Convenient imaging with two media sizes on-line (multi-format)
- Very short access time ensures fast printing of small print jobs

▶ Agfa HealthCare Drystar 5300

Capacity Resolution System 70 films/h (14 x 17")

320 dpi

Direct Digital Imaging



▶ Highlights

- Tabletop, next-to-application Direct Digital Imager
- Suitable for all applications and ideal for CT/MR
- Reliable, low maintenance printer
- A#Sharp image enhancement for excellent quality
- Very short access time ensures fast printing of small print jobs

▶ Agfa HealthCare Drystar AXYS

Capacity Resolution System

75 films/h (14 x 17") 508 dpi/50 µm pixelsize



▶ Highlights

- -Flexible, tabletop imager delivering mammography-quality images
- Multi-application hardcopy solution
 Integrated A#Sharp technology for optimized image quality
- 2 multi-format trays, each supporting different film sizes and types
- Very short access time for extremely fast delivery of first four prints

▶ CPS - Codonics Horizon Ci / GS / SF

▶ Highlights

- Read in room light with convenient grayscale or color paper (Ci, GS and SF)
- Outstanding image quality, prints on all standard sizes of clear/blue dry film
 Compact design with the fastest print speed in its class
 Economical surrise express swap
- warranty provides a replacement imager
- Network printer with up to 24 DICOM connections, no additional hardware required
- Validated with all major OEMs, modalities, PACS and workstations



Horizon Ci

100 sheets/h 300 copies

Capacity Resolution System

320 dpi – 126 cm spatial resolution, 4096 pixel contrast resolution, shades gray Direct thermographic (grayscale) and dye-diffusion (color)

Horizon GS

Speed Capacity Resolution

Speed

100 sheets/h 300 copies

320 dpi - 126 cm spatial resolution, 4096 pixel contrast resolution, shades gray Direct thermographic (grayscale) System

Horizon SF

Speed Capacity Resolution

100 sheets/h 300 copies

320 dpi – 126 cm spatial resolution, 4096 pixel contrast resolution, shades gray System Direct thermographic (grayscale) and dye-diffusion (color)

▶ CPS DICOM PaperPrint Server

Color Laser Printer Copier Scanner Paper sizes

Print resolution

Printout

yes

DocuColor 242 WorkCentre 7665 yes A3/A4 A3/A4 and many other and many other 2400 x 2400 dpi 2400 x 2400 dpi

color, black & white color, black & white

WorkCentre 7425 yes

yes A3/A4 and many other

1200 x 1200, 1200 x 2400dpi color, black & white

Phaser 5550

► Highlights - Full DICOM 3.0 Basic Grayscale and Color Print Management

- Single licence for unlimited number of modalities and printers

- Individual header and footer including text and graphic logos

Individual LUT in color and black & white for every modality, user and printer

Print-Presentation-LUT

Full immediate onlinesupport for printer, server and software by CPS
-Low average cost per

page, about 1 cent in black&white and 8 cent in color

Color Laser Printer Copier Scanner Paper sizes

Print resolution

Printout

A3/A4 and many other 1200 x 1200 dpi

Phaser 7760

A3/A4 and many other 1200 x 1200 dpi color, black & white color, black & white

Phaser 7500

A3/A4 and many other 1200 x 1200 dpi black & white

► CPS - Codonics Horizon XL

Speed Capacity Resolution

System

100 sheets/h 300 copies

320 dpi - 126 cm spatial resolution, 4096 pixel contrast resolution, shades gray Direct thermographic (grayscale) and dye-diffusion

(color)

▶ Highlights

– 14" x 36" and 14" x 51"

long dry film

Perfect for scoliosis, long bone studies

-»True-size« imaging up to 51" in length

- Saves space and eliminates wet film processing

 Also prints on standard sizes of film

- Low-cost grayscale paper and color paper



Fujifilm DryPix 4000

Capacity Resolution System

160 films/h 50 μm/14 bit Laser



▶ Highlights

5 available formats: from 8"x10" up to 14"x17" from (20x25cm) up to (35x43cm)

-1 or 2 film trays

- On-top sorter optional available

- High resolution for mammography

-High density for mammography (Ď_{max}=4.0)

▶ Fujifilm DryPix 7000

Capacity Resolution System 240 films/h 50 μm/14 bit Laser



▶ Highlights

-5 available formats: from 8"x10" up to 14"x17" from (20x25cm) up to (35x43cm)

-2 or 3 film trays

– 10 bin sorter optional available

- High resolution for mammography

− High density for mammography (Ď_{max}=4.0)

▶ Fujifilm DryPix Prima

Capacity Resolution System up to 70 films /h 100 μm / 14 bit Laser



Highlights

- Tabletop laser printer - 0,38 m² Footprint - Support for 5 different film formats - Fully DICOM compatible



Capacity Resolution System

90 films/h $78,6~\mu m/12~bit$ Laser



▶ Highlights

- Smallest laser imager Fastest time for first film print out (50 s) Ready for up to 2 film trays Support of 5 different film sizes

▶ Konica Minolta DryPro 793

Capacity Resolution System

120 films/h 43.75 µm/14 bit Laser



- ► Highlights

 Central print solution
 for multi-modality environment

 Ready for up to 3 film trays

 Support 5 different film sizes

 Contained sorter available

- Optional sorter available
- Supports mammography

▶ Konica Minolta DryPro 873

Capacity Resolution System

180 films/h 43,75 µm Laser



▶ Highlights

Paper Sizes

(max.) Resolution

System

- Supports all DICOM 3.0 modalities (e.g. CT, MRT, CR, DR, US, NUK etc)
- Supports one or more PostScript printers within the network
- General licence package (no restrictions on how many DICOM modalities are connected)
- Image header and footer customizable incl. physician logo
- Separate LUT (Look Up Table) for each printing system

GŜDF calibration according IHE

for optimal performance - High density printing for mammography - Dmax 4.0 - Fully DICOM compatible

▶ Highlights

- Ready for up to 3 film trays
- Optional sorter available

► Mitsubishi Printer P93E

Capacity Resolution System

Approx. 260 print roll, max. 923 sheets/h 325 dpi/1280 x 600 (PAL/normal) B&w video printer, direct thermal



▶ Highlights

- User-friendly settings with control switches
- Extremely compact dimensions and lightweight design
- Fast print speed due to BAS/FBAS, BNC connector -7 different picture formats
- Picture memory for 10 individually selectable frames

▶ Mitsubishi Printer CP30W

▶ medigration DICOM PaperPrint

DIN A3, 11 x 17

Laser, scan, fax

1.200 x 2.400 dpi (print),

600 x 600 dpi (copy)

Capacity Resolution System 80 prints set, max. 225 sheets/h 423 dpi/1600 x 2100

Color video printer, dye sublimation



- Front-loading system on rails for paper and ink sheet cassette
 Compact and ergonomic design

- Integrated control panel
 Integrated paper tray with illuminated exit slot
 PAL & NTSC compatible, all standard interfaces

► Mitsubishi Printer P93DW

Capacity Resolution System

Approx. 260 print roll, max. 973 sheets/h 325 dpi/1280 x 5760 (panorama) B&w digital printer, direct thermal



- ► Highlights

 USB Version 2.0 guarantees print speed from 3.7 seconds

 Compact size of only 154 x 90 x 256 mm

 Quick and easy adjustment on the front panel

- Panorama-print up to 100 x 150 mm
- Extensive adjustment possibilities using the printer driver

▶ Mitsubishi Printer CP31W

Capacity Resolution System 80 prints set, max. 225 sheets/h

423 dpi/1600 x 2100

Color video printer, dye sublimation



▶ Highlights

- Front-loading system on rails for paper and ink sheet cassette
- Compact and ergonomic design
- Integrated control panel
- Integrated paper tray with illuminated exit slot PAL & NTSC compatible with common S-Video and Composite Video

Mitsubishi Printer P95DE

Capacity Resolution System Approx. 245 print roll, max. 1,895 sheets/h 325 dpi with 16-bit data processing B&w digital printer, direct thermal



▶ Highlights

- Print time per image 1.9 sec. (100 x 75 mm) Compact size of only 154 x 84,5 x 239 mm
- -Exceptional long print format possible (100 x 450 mm)
- New high-gloss thermal paper
- Optimised heating process thanks to Print Control Engine

Mitsubishi Printer CP30DW

Capacity Resolution System 80 prints set, max. 225 sheets/h

423 dpi/1600 x 2100

Color digital printer, dye sublimation



- ► Highlights

 Front-loading system on rails for paper and ink sheet cassette
- Compact and ergonomic design High-speed USB interface (Version 2.0)
- Large integrated paper tray
- Illuminated paper exit slot

Sony UP-DF 750

Capacity Resolution System 75 sheets/h 604 dpi Direct thermal



Highlights

- Multiple Film Sizes for Various Radiology Modalities
 New Blue Thermal Film for Mammography

- New Blac Filmar Film for Mailing Lapiny
 Two Film Trays Compatible with All Film Sizes
 Large Effective Print Area and Edge-to-edge-like Printing
 Easy Network Parameter Settings
 DICOM Connectivity

Sony UP-DR80MD

Capacity Resolution 50 sheets/h 301 dpi


System Dye Sublimation Thermal Printing



- Superb Print QualityStylish Design and Compact Size
- Useful Front Operation
- Easy Color Adjustment Color Look Up Table
- Gray Balance Adjustment



Printing Technology HP 2800 (A3+, A3, A4, A5, A6)



▶ Tetenal spectra™ jet Film

Printing Technology HP 1000 (A4, A5) Resolution b/w 4800_1200 dpi



▶ Highlights

- High resolution quality low costs
- High density
- Fast-drying
- Water-resistant
- Transparent + solid printouts

DISPLAYS/PRINTERS ACCESSORIES

▶ CPS - Codonics Virtua Medical Disc Publisher

▶ Highlights

- -Burns up to 60 CDs or 30 DVDs an hour
- Auto records patient studies and reports without tying up workstation or employee resources
- Touchscreen interface for optimized workflow
- Full-color disc labels creator
- DICOM compliant network appliance
 Burn speeds based on a typical clinical study with full color label. Not all features available on all models. Specifications subject to change

Capacity Drives/Recordable Format Printer

Virtua Medical Disc Publisher

Two 50-disc input bins; 100 disc total capacity DVD±R / CD-R dual-layer drives, CD-R, DVD-R, DVD+R Inkjet 4800 dpi

Capacity Drives/Recordable Format Printer

Virtua XR Medical Disc Publisher

Two 50-disc input bins; 100 disc total capacity Two dual-layer, DVD±R / CD-R drives, CD-R, DVD-R, DVD+R Inkjet 4800 dpi

▶ CPS – Codonics Integrity Medical Image Importer

User Interface Remote web browser access

DICOM 3.0, IHE PDI, ACR NEMA, older DICOM image files Import Formats **Processor** Intel® Core™ 2 Duo

Search Rules Configurable 6.49" (16.5 cm) W, 6.49" (16.5 cm) D, 1.96" (5 cm) H Dimensions Weight

2.41 lbs. (1.46 kg)

▶ Highlights

- -Compact, stand alone solution to read, reconcile and store medical studies from CD/DVD
- Improves workflow by bringing the reconciliation process to the user
- Automatically scans for viruses to protect your data

 Reconciles patient data
- with facility's own modality worklist
- Displays the original imported data as well as the matching MWL or PACS data



▶ IBA Dosimetry LXcan

Spot luminance meter for quality tests at displays



- Luminance and illuminance measurementsDisplay: 1.2" TFT
- Targeting device: integrated camera
- Ultrasonic distance sensor; alignment sensor
- USB interface

RAD-BOOK 2010



► Esaote MyLab25 XVGold

Modes

B-mode, M-mode, color-, high sensitivity power doppler, PW-, CW doppler, TEI, CMM, TVM, TP-View, VPan 3D/4D

Scan format Transducer inputs

2+1 probe connectors. Probes: LA, CA, PA, microconvex, pencil

- ► Highlights

 15"LCD monitor

 high frequency imaging up to 18 MHz

- brilliant images (XView)
 Compound Imaging (MView)
 CnTITM (contrast enhanced ultrasound)
- PC-workstation MyLabDesk™
- Li-Ion battery (up to 1h) mobile system, DICOM



Modes

B-mode, colorized 2D, CFM- mode, TEI, Bmode steering LA, CA and microconvex

Transducer probes

▶ Highlights

Arm-held

- wearable unit
- big 12" LCD real touch screen monitor
- Tutorial
- programmable keys on the LA-probe,
- for activating functions on the scanner remotely
- brilliant images (XView), user defined protocol editor
- $PC\text{-}work station \ MyLabDesk^{\tiny{TM}}$
- 2 USB ports
- Li-Ion battery (up to 3h)
- DICOM

► Esaote MyLab Twice



Modes

B-mode, M-mode, color-, high sensitivity power doppler, PW-, CWdoppler, TEI, TVM, CMM, TP-View, VPan

Scan format


3D/4D linear & convex

Transducer inputs 4+1 probe connectors LA, CA, PA, microconvex, pencil

▶ Highlights

- 19" LCD monitor
- -HF imaging 18 MHz
- XView, Compound Imaging MView, CnTI™, Elastographie, QIMT, Fusion Imaging, MyLabDesk™
- DÍCOM, 4 USB ports, 500 GB, Satellite system

GE Healthcare Logiq E9



Scan format

Modes

Transducer inputs

B-Mode, M-Mode, CFM-Mode, Doppler, Amplitude modulated contrast mode, Realtime4D Linear, convex, microconvex, sector phased array, trapezoid

▶ Highlights

- Matrix array transducer technology
- 3D/4D Volume scan
- Depth independed contrast imaging thanks to new amplitude modulation technology
- True spatial image fusion of CT/MRT Images and realtime ultrasound
- Volume navigation

Modes

- Agile, adaptive beamformer

GE Healthcare Voluson E 8



Modes

Scan format

B-Mode, M-Mode, CFM-Mode, Doppler, HD-Flow, Realtime4D Linear, convex. microconvex, sector phased array

inputs

▶ Highlights

3

- Realtime 4D up to 40 volumes/sec.
- Automatic volumetric analysis
- -STIČ (Realtime 4D view
- of the fetal heart)
 -CRI (Compound Resolution Imaging)
- HD-Flow (high sensitive Power Doppler)

GE Healthcare Logiq 9



Scan format

B-Mode, M-Mode, CFM-Mode, Doppler, B-Flow, coded contrast Harmonic, Realtime4D Linear, convex, microconvex, sector phased array, trapezoid

Transducer inputs

Highlights

Volume ultrasound (3D and contrast harmonic imaging, VOCAL II, 16 Mhz volume probe)

- TUI (Tomographic Ultrasound Imaging)
- CrossBeam realtime compound (up to 9 angles)
- Matrix array transducer technology LOGIQView (panoramic imaging)
- Ergonomic design with swiveling
- keyboard, LCD monitor, VoiceScan

GE Healthcare Logiq 7



Modes

Scan format

Transducer inputs

B-Mode, M-Mode, CFM-Mode, Doppler, B-Flow color, Coded Contrast Harmonic, StressEcho, anatomical M-Mode Linear, convex, microconvex, sector phased array, trapezoid

▶ Highlights

High end shared service system

B-Flow color (digitally subtraction technique)
Matrix array transducer support

CrossBeam realtime compound and speckle reduction imaging

Contrast harmonic imaging with DualView and TIC analysis

► GE Healthcare Logiq S6



B-Mode, M-Mode, CFM-Modes Mode, Doppler, B-Flow

colour, Coded Contrast Harmonic, StressEcho, Anatomical M-Mode Linear, convex, microconvex, sector phased array, trapezoid

inputs

▶ Highlights

- Compact shared service system - B-Flow color (digitally

subtraction technique) Matrix array transducer support

CrossBeam realtime compound and speckle reduction imaging

Digitally archive with RawData support

GE Healthcare Logiq P6



Modes

B-Mode, M-Mode, CFM-Mode, Doppler, B-Flow color, Coded Contrast Harmonic, StressEcho, Anatomical M-Mode Linear, convex, microconvex, sector phased array, trapezoid

Transducer inputs

Scan format

Highlights

Compact shared service system B-Flow color (digitally substraction technique)

CrossBeam realtime compound and speckle reduction imaging

LOGIQView (panoramic imaging) Auto optimize (For B-Mode, color, Doppler)

Digitally archive with RawData support

GE Healthcare Logiq A5 / P5 Premium



Modes

Modular configurable from b/w system up to color triplex system (B-Mode, M-Mode, CFM-Mode, Doppler, B-Flow,

cardiology) Linear, convex. microconvex, sector phased array, trapezoid

Transducer inputs

Highlights

Scan format

- Compact lightweight and modern design with 15" LCD monitor

CrossBeam and speckle reduction

Imaging
- LOGIQView (panoramic imaging)
- Auto optimize (for B-Mode, color,

Doppler) Digitally archive with RawData support

▶ GE Healthcare Logiq e

Modes Scan format B-Mode, M-Mode, CFM-Mode, Doppler Linear, convex, microconvex, sector phased array, trapezoid

Transducer inputs



capabilities High frequency modul

(for vascular and SmallParts diagnostic) Sector phased array modul (for cardiology)

– CrossBeam

▶ Highlights

LOGIQView (panoramic imaging)

▶ GE Healthcare Venue™ 40



Scan format

Modes

Transducer inputs B mode, Color flow Imaging, Power Doppler Linear, Convex, Sector (Phased array)

Highlights

No buttons. No knobs. No keyboard. Easy to use at the point of care.

Concurrent acquisition technology provides fast, high-resolution imaging to easily visualize anatomy and needle placement.
Depth-synchronized optimization

with adjustable gain. CrossXBeamTM and Speckle Reduction Imaging (SRI).

Single-surface screen – no seams, no monitor frame.

ULTRASOUND

Hitachi Medical Systems HI Vision 900



Modes

B & M Mode; Omni Directional M Mode; PW & CW Doppler; Dual Gate Doppler; Colour & power Doppler; Fine Flow Mode; Triplex Mode; TDI; Tissue elastography; Contrast harmonic imaging; 4D; Real-time Virtual Sonography; Real-time Bi-plane

Sector, linear, and convex array, 360° scanning, dual imaging, wideview panoramic, HI definition Zoom, Pan Zoom; Picture in Picture

Transducer

inputs 3 active

Highlights

- 3 types tissue harmonic imaging (6 choice of frequencies)
 Tissue adaptive filtering, HI Rez+ (4 levels) for speckle and noise reduction
 Compound imaging, HI Com (from multiple directions and different fre-
- quencies)
- Coded Imaging for increased penetration from high frequency transducers
- Option of waterproof remote control operation and voice-activated Bluetooth operation

Hitachi Medical Systems HI Vision Preirus



Modes

B & M Mode; Omni Directional M Mode; PW & CW Doppler; Dual Gate Doppler; Colour & power Doppler; Fine Flow Mode; Triplex Mode; TDI; Tissue elastography; Contrast harmonic imaging; 4D; Real-time Virtual Sonography; Real-time Bi-plane Sector, linear and convex array, 360° scanning, trapezoid, B-steer, dual imaging, wideview panoramic, HI definition Zoom, Pan Zoom; Picture in Picture


Transducer inputs

3 active

▶ Highlights

- Unique ergonomic design gives increased system flexibility
- 3 types tissue harmonic imaging (6 choice of frequencies)
- Tissue adaptive filtering, HI Rez+ (4 levels) for speckle & noise reduction
 Compound imaging, HI Com (from multiple directions and different
- Graphical User Interface incorporating Smart Tab menus, Image Thumbnails and Touch Screen panel for image optimisation

litachi Medical Systems EUB-7500 HV



B & M Mode; Omni Directional M Mode; PW & CW Doppler; Dual Gate Doppler; Colour & power Doppler; Fine Flow Mode; Triplex Mode; TDI; Tissue elastography; Contrast harmonic imaging; 4D; Real-time Virtual Sonography; Real-time Bi-plane Sector, linear and convex array, 360° scanning, dual imaging, wideview panoramic, HI definition Zoom, Pan Zoom; Picture in Picture

3 active ports

Highlights

- 3 types tissue harmonic imaging (6 choice of frequencies)
- Tissue adaptive filtering, HI Rez+ (4 levels) for speckle and noise reduction
 Compound imaging, HI Com (from multiple directions and different
- frequencies)
- Option of waterproof remote control operation and voice-activated Bluetooth operation
- Compatibility with wide range of transducers including endoscopic/ bronchoscopic options

SonoScape SSI-8000 The Cutting-edge Technology for Accuracy and Clarity An Ever-higher Standard of Echo HCU Winter of the "2009 Product Quality Leadership Awar from Frest & Sollivan. **A6**




Ultimate Performance. Like Nothing Else.







Hitachi Medical Systems HI Vision Avius



B & M Mode; Omni Directional M Mode; PW & CW Doppler; Colour & power Doppler; Fine Flow Mode; Triplex Mode; TDI; Tissue elastography; Contrast harmonic imaging; 4D Sector (phased), linear, and convex array, 360° scanning, trapezoid, B-steer, dual imaging, Wideview panoramic, HI definition Zoom, Pan Zoom; Picture in Picture

Transducer inputs

3 active ports

▶ Highlights

- 3 types tissue harmonic imaging (6 choice of frequencies)
 Tissue adaptive filtering, HI Rez+ (4 levels) for speckle & noise reduction
 Compound imaging, HI Com (from multiple directions and different frequencies)
- Graphical User Interface incorporating Smart Tab menus, Image Thumbnails for image optimisation
- PSS, patient specific scanning selector

▶ Hitachi Medical Systems EUB-7000 HV



Modes

Scan format

Transducer

B & M Mode; Omni Directional M Mode; PW and CW Doppler; Colour and power Doppler; TDI; Tissue elastography; 4D

Sector (phased), linear, and convex array, 360 degrees radial scanning, dual imaging, wide view panoramic, HI definition Zoom, Pan Zoom; Picture in Picture


inputs 3 active ports


▶ Highlights


- 3 types tissue harmonic imaging (6 choice of frequencies)
 Tissue adaptive filtering, HI Rez (4 levels) for speckle and noise reduction
- Option of waterproof remote control operation and voice-activated Bluetooth operation
- Compatibility with wide range of transducers including endoscopic/bronchoscopic, surgery and urology options
- Real-time tissue elastography for breast, prostate, pancreas, thyroid, musculoskeletal, and many more

▶ Landwind Mirror2Plus Modes B-mode, M-mode, color-, high sensitivity power doppler, PW Doppler£" CW Doppler Scan format Phased, Linear, convex, micro-convex Transducer 3 active ports inputs ▶ Highlights - Tissue Doppler imaging - Color Doppler panoramic imaging - High resolution write zoom - MCI (Multi-angle Compound Imaging) - Phase inversion tissue harmonic imaging - Color M mode





















Medison Sonoace X4



Modes

Scan format Transducer

inputs

M-mode and PW Spectral Doppler Linear, convex and micro convex

B-mode, freehand 3D.

1 (opt. 2)

Highlights

- Premium B/W and Doppler ultrasound system
- Digital multi-beam-forming channels
- Full Spectrum™ imaging
- -Free hand 3D imaging

▶ Medison Sonoace X6



B-mode, Tissue- and Pulse inversion Harmonic Imaging, DynamicMR™, M-mode, Anatomical M-mode, CFM-mode, Power Doppler, Spectral Doppler (PWD/CWD), Pulsed Wave Tissue Doppler Imaging Linear, trapezoidal, compound linear, convex, micro convex, 3D convex, phased array sector and pencil

▶ Highlights

B-mode, Tissue- and Pulse

DynamicMR™, M-mode,

inversion Harmonic Imaging,

Tissue Doppler Imaging, Live

Linear, trapezoidal, compound

linear, 3D linear, convex, micro

convex, 3D Convex, phased

array sector and pencil

3D™/4D and 3D XI™

Anatomical M-mode, CFM-mode, Power Doppler, Spectral Doppler (PWD/CWD), Pulsed Wave

- Economical multi-speciality digital color ultrasound system
 Full Spectrum™ and Tissue Harmonic Imaging
- Color and Power Doppler
- Free hand 3D imaging
- High resolution (LCD monitor)

Medison Sonoace X8

Modes



Scan

format

Transducer

B-mode, Tissue- and Pulse inversion Harmonic Imaging, DynamicMR[™], M-mode, Anatomical M-mode, CFM-mode, Power Doppler, Spectral Doppler (PWD/CWD), Pulsed Wave Tissue Doppler Imaging, Live 3D™/4D and 3D XI™

Linear, trapezoidal, compound line ar, convex, micro convex, 3D convex, phased array sector and pencil

3 + 1

inputs

▶ Highlights

- Popular class multi-speciality live 3D™/4D ultrasound system High resolution (1280 x 1024) 17" LCD monitor

- Live 3D™ with extreme volume rates 3D XI™, Multi Slice View™, Oblique view™ and VolumeCT™
- Highly sensitive directional Power Doppler Auto IMT, CEUS, Elastoscan

Medison Accuvix V20 Prestige

Modes



▶ Highlights

Scan format

Transducer

Multi-application inputs 3D/4D ultrasound system

High resolution 19" LCD monitor

- Auto IMT, CEUS, Elastoscan

4 + 1

Auto IMT, CELOS, Elastoscali
Panoramic View 2.0, Spatial Compound Imaging, Speckle Reduction
Filter, DynamicMR™ PLUS
XI STIC™, BW STIC, Color STIC™, STIC + MSV™, 3D MXI™
,Multi Volume Slice™, Mirror View™, Multi-OVIX™
3D DMRTM, VCETM (Volume Contrast Enhancement)
Highly conciding directional Payor Departs

Highly sensitive directional Power Doppler

Scan format

Medison Accuvix V10

B-mode, Tissue- and Pulse inversion Harmonic Imaging, Dynamic MR^{TM} , M-mode, Anatomical M-mode, CFM-mode, Power Doppler, Spectral Doppler (PWD/CWD), Pulsed Wave Tissue Doppler Imaging, Live $3D^{\text{TM}}/4D$ and $3D XI^{\text{TM}}$ Linear, trapezoidal, compound linear, convex, micro convex, 3D convex, phased array sector and pencil

Transducer inputs ▶ Highlights

Multi Application live 3D/4D ultrasound system
 High resolution 19"LCD monitor
 Live 3D™ with extreme volume rates
 3D XI ™ Multi Slice View™, Oblique view™ and VolumeCT™
 3D XI STIC™, VOCAL and XI VOCAL™

- Highly sensitive directional Power Doppler

- CEUS, Elastoscan

Medison MySono U5



Transducer inputs B-mode. Tissue- and Pulseinversion Harmonic Imaging, Dynamic MRTM, Mmode, Anatomical M-mode, CFM-mode, Power Doppler, Spectral Doppler(PWD/CWD). Pulsed WaveTissue Doppler Imaging, Live3DTM/4D and

Linear, trapezoidal, compound linear, convex, micro convex. 3D convex, phased array sector

1 port

▶ Highlights

Multi Application portable live 3D/4D ultrasound system
 High resolution 15"LCD monitor

- High resolution 13 Lee Informer

- Live 3D[™] with extreme volume rates

- 3D XITMMulti Slice View[™], Oblique view[™] and
VolumeCT[™],-3D XI STIC[™], VOCAL, XI VOCAL[™], Auto IMT

- Highly sensitive directional Doppler

► Medison Sonoace X1



Modes Scan format Transducer inputs B-mode, BB mode, 4B, BM, M mode Linear, convex

2 port

▶ Highlights

- Multi Application portable ultrasound system– 10" CRT monitor
- Large image storage
- USB direct storage
- High resolution zoom

► Mindray DC-3

Modes

B,2B,4B, B/M, M, Color Doppler Flowing Imaging, HPRF, Power, Dirpower

Scan format

Linear Array, Convex Array, Micro Convex, Endocavity, Intrarectal, Phased array, T-type linear, Biplanar

Transducer inputs



▶ Highlights

- Compact system with full ergonomic design
- iClear™: Speckle reduction technology
- THI(Tissue Harmonic Imaging), Smart 3D™, Trapezoid imaging, iScape™ view and Free Xros™
- imaging iTouchTM: Intelligent one-touch image optimization

 iStation™: Intelligent patient
- management platform DVD-R/W, USB, DICOM 3.0,
- and ECG module
- Height adjustable and rotatable control panel

► Mindray DC-7

Technology

iClear, iBeam, Frequency Compounding Imaging,

Adaptive Frame Averiging

Modes

B, C, M, PW, CW, Power (DirPower),

TDI, CM (Color M), 4D

Scan format

Convex, Linear, phased array, convex

volume, endovaginal volume

Transducer inputs 2-15 Mhz



Highlights

- High-resolution LCD screen, color touch screen
- Four active probe connectors, wide bandwith and high density integrated probes
- Simplified workflow with fast system response
- Value-added performance in versatile fields of applications
- Rich use-define packages and specific reports

► Mindray M7

Microsoft Windows XP Professional Technology

with ASIC Design Platform B / M / Anatomical M /

Color Doppler Velocity / Power (DirPower) / PW / CW /Smart 3D /Static 3D / 4D/

iScape (panoramic imaging) / TDI (Tissue Doppler Imaging) Color M (CM)

Linear, Convex, Phase,

Endocavity, Volume

Scan format

Modes

Transducer inputs

2-15 Mhz



Phase Inversion THI

iBeam™ Spatial compound imaging

iClear™ Adaptive Speckle Suppression Imaging iTouch™ Intelligent Image Optimization

iZoom™ Automatically expand the image to full screen

iRoam™, 802.11b/g wireless data transfer solution

► Mindray M5

Technology Modes

Microsoft Windows XP Embedded B / M / Color Doppler Velocity / Power (DirPower)

/ PW / CW / Smart 3D / Static 3D / iScape (Panoramic imaging)

Scan format

Linear, Convex, Phase, Endocavity

Transducer inputs 2-15 Mhz



▶ Highlights

- -THI

- iBeam™ Spatial compound imaging
 iClear™ Adaptive Speckle Suppression Imaging
 iTouch™ Intelligent Image Optimization
 iZoom™ Automatically expand the image to full screen
- iStation Powefull patient data management platform

▶ Philips HD7

Modes

Microfine 2D focusing, Color Doppler, Color Power Angio imaging, Pulsed Wave Doppler, 3D, Stress Echo, M-mode and Anatomical M-mode, pulse inversion, Tissue Harmonic Imaging' DICOM SR and contrast Curved, linear and sector arrays

Scan format Transducer

inputs up to 4

- User centric functions and design
- Exceptional HD performance with Broadband beamformer and broadband transducers
- Optimized workflow with simple 2D, color and spectral Doppler optimization
- All type of clinical applications:
- Data Management with USB access through user interface
- Transducer compatibility with other Philips systems



▶ Philips iU22

Modes

2D, Tissue Harmonic, M-Mode, advanced volumetric modes including Live xPlane (simultaneous display of 2 live images), PW-Doppler, CW-Doppler, panoramic, contrast, Color Power Angio

Scan format Transducer inputs Curved, linear, sector (PureWave) and xMATRIX arrays



- Exceptional image quality as a result of leading-edge technologies

Advanced volumetric imaging with freehand, automated and electronic (xMATRIX) acquisition Unique system design provides

unparalleled ergonomics

Built-in automation - optimal image quality with the push of a button

Full range of applications including abdominal, ob/gyn, TCD, musculoskeletal, adult cardiac (including stress echo), vascular and small parts

▶ Philips HD11 XE

Modes

2D, Tissue Harmonic, M-Mode, PW-Doppler, Color Power Angio, CW-Doppler, panoramic, trapezoidal, contrast

Scan format Curved and linear arrays

Transducer inputs

Highlights - Fully equipped to cover a full range of applications: abdominal, small parts/musculokeletal, ob/gyn, vasclar, cardiac and TCD

Advanced imaging capabilities with volumetric imaging and manipulation tools: iSlice, STIC, Colour invert

Cardiac imaging supported with stress echo, anatomical M-mode and tissue Doppler imaging
Superb ergnomics and mobility – the

system goes where you need it



▶ Philips EnVisor HD

Modes

2D, Tissue Harmonic, PW-Doppler, CW-Doppler, M-Mode, trapezoid, panoramic, 3D, contrast

Scan format Transducer inputs

▶ Highlights


- Optimized workflow for the demanding practice

- Full range of transducers

support a wide variety of exams, including abdominal, cardiac (including TEE), ob/gyn, vascular and musculoskeletal/small parts

Workstation quality data management at the point of

- Affordably priced





Modes

Scan format Transducer inputs B-Mode, M-Mode, Color Doppler, VeloPower, PW-Doppler, HPRF, CW-Doppler, Tissue Harmonic Imaging, Color Power Angio Imaging, Grayscale 3D Curved and linear arrays

3 + 1



▶ Highlights

- Powerful capabilities in a compact package

– HD performance in a small,

easy-to-use, mobile system

Easy to maneuver and ergonomically adjustable

Extremely affordable

Siemens Acuson S2000



Modes

B-mode, Color Doppler, Power Doppler, PW Doppler (Duplex, Triplex), Doppler Tissue Imaging (Color and PW), CW spectral Doppler, M-mode and Color Doppler Mmode

Curved array, phased array, linear, endocavity, 3D/4D imaging, pencil

Transducer inputs

Scan format

▶ Highlights

Advanced transducer technology including micro-pinless connectors, Hanafy lens and matrix arrays, and silicon-ready

Advanced breast imaging application with eSieTouch™ elasticity imaging and Fatty Tissue Imaging technologies including option to add ABVS Automated Breast (see Mammography)

 Advanced SieClear™ spatial compounding with dynamic TCE™ technology with speckle reduction in 3D

Advanced fourSight™ technology

Automatic measurement of lesions with syngo® e-Sie Calcs native tracing software

▶ Siemens Acuson Antares



High-end ultrasound system

application with eSieTouch™ elasticity imag-

ing and fatty tissue imaging

Advanced breast imaging

▶ Highlights

3D/4D imaging

technologies

Scan for mat

Modes

Transducer inputs B-mode, Color Doppler, Power Doppler, PW Doppler (Duplex, Triplex), Doppler Tissue Imaging (Color and PW), CW spectral Doppler, M-mode and Color Doppler M-mode Curved array, phased array, linear,

endocavity, 3D/4D imaging, pencil

- Cadence CPS Contrast enhanced imaging

– Hanafy lens transducer technology

- MultiHertzTM multiple frequency

 MultiHertz¹ infutiple frequency imaging technology
 Advanced SieClearTM spatial compounding with dynamic TCETM technology
 Advanced fourSightTM technology
 TEQTM ultrasound technology: ClarifyTM vascular enhancement technology, syngo® auto OB measurements

ULTRASOUND

Siemens Acuson X300 Premium Edition

Modes

B-mode, Color M-mode, Mmode, Color Doppler velocity mode, Power Doppler mode, Pulsed Wave spectral Doppler mode (PW), Continuous Wave spectral Doppler mode (CW), Duplex mode, Triplex mode Curved array, phased array, linear, endocavity, 3D/4D imaging

Scan format

Transducer inputs

▶ Highlights

- Excellent imaging performance through excellent detail and contrast resolution high temporal resolution in 2D
- TGOTM tissue grayscale optimization technology
- for more consistent image quality
- High quality 4D imaging through Advanced fourSight™ technologies
 Exceptional clinical performance across a variety of
- applications and patient body types
- Easy-to-use ErgoDynamic imaging system design

► Siemens Acuson X300

Modes



B-mode, Color M-mode, M-mode, Color Doppler Velocity mode, Power Doppler mode, Pulsed Wave (PW) spectral Doppler mode, CW Continuous Wave spectral Doppler mode Phased array, curved array, endocavity. linear array

Transducer inputs

Scan format

▶ Highlights

- Hanafy lens transducer technology
- Tissue harmonic imaging
- DTI™ Doppler tissue imaging capability
- Multi-beam formation technology
- Streamlined clinical workflow with integrated DIMAQ-IP workstation, a user customizable control panel, and TGOTM tissue grayscale optimiza-
- tion technology
 ErgoDynamic™ imaging system design with flat
 panel display and articulating arm

Siemens Acuson X150

Modes



B-mode, M-mode, Color Doppler Velocity mode, Power Doppler mode, Pulsed Wave (PW) spectral Doppler mode, Duplex mode, Triplex mode

Scan format

Phased array, curved array, endocavity, linear arrav 2 + 1 optional

Transducer inputs

▶ Highlights

- Top diagnostic performance and scalability
- Superior 2D-mode imaging
- Color imaging option
- Cardiac screening option and phased array transducers fully integrate 3-Scape™ real-time 3D imaging application to easily acquire realtime 3D images during freehand acquisition
- Intuitively simple, yet powerful user interface with highly functional ergonomics

RAD-BOOK 2010

see it all o MEDISON

3D/4D Beyond Imagination The President Associated Imagination States of Property of Proper

technological advancement on 30/40 ultrascond imaging at MEDISON, Inferioring a tradition of excellence in 30/40 technology and diagnostic application from the ACCLMIX line, the flagship (pestige) sets the standard in 30/40 altranound imaging with breakthrough volume acquisition technologies, faster and innocts operation, and significantly enhanced imaging capabilities.

Discover the new 'Prestige' and experience 3D/4D imaging beyond your imagination ACCUVIX V20

The legacy continues....

With Introducing the MySono US-sitrassund system, MEDISON has ach New Dimensions, MEDISON's well known 2D image quality and advanced 30/4D technologies are now available in portable, small and sheek designed. MySono US. Its completeness, small size and light weight makes the MySono US an ultrascond machine that is unique in its class. Its origineness lies also in the intuitive user interface, fast response and easy operation, MySono US. proves that size and design are not compromised to achieve excellent im-



3D XITM(3D+XXxxxid

Imaging) - Multi-SliceTM Oblique ViewTM

Volume CTTM

DHR SWTH(D) ic HR lite?

Full package

Wireless LAN DICOM 3.0

Pointer Function TM -USB2.0





- Offline image review software

Siemens Acuson P10 B-Mode, harmonic Modes modes Scan format Phased array Single handheld unit Transducer with integrated inputs transducer ▶ Highlights - Excellent image quality - Instant power-up - Removable, rechargeable battery Simple, intuitive user interface TGO™ tissue grayscale optimization technology Application presets - SD memory card and USB port

► SonoScape A6 Portable B/W

Modes Scan format Transducer inputs

4B, Tissue Harmonic Imaging, B, B/M, M, B/B Linear-, Convex-, Micro Convex Array



▶ Highlights

 Portability with Position adjustable large 12" LCD monitor

Weight less than 6kg with a convertible cart-based system design

Super high physical processing channel
 HD vary broadband 5-frequency-selection probes

- Full patient database solutions: DVD RW, PDF report, Dual USB 2.0 and DICOM 3.0, etc.



► SonoScape SSI-8000 Color Doppler

Technology Modes

Scan format

Micro Scan Processing B-mode, Steer M-mode, TDI, CFM, PDI, PWD, CWD, 3D/4D, Color M-mode, Tissue Doppler Imaging

Linear-, Convex-, Micro Convex-, Phased Array, TEE, Volume Convex

Transducer inputs

▶ Highlights

- Premium multi-specialty 4D ultrasound

- A full range of cardiovascular transducers and Cardiovascular analysis

- Integrated with state-of-the-art imaging technologies, like microscan, multiple-beam processing, etc.

200 degree transvaginal image

- Full patient database solutions: DICOM3.0, USB2.0, HDD, etc.



► SonoScape A8 B/W

Modes Scan format Transducer inputs

patient database solutions

4B, B, B/M, M, B/B, Tissue Harmonic Imaging Linear-, Convex-, Micro Convex Array

▶ Highlights

- Compact and agile trolley design with standard three transducer
- Intuitive operation with M-Tuning one-key image optimization
- 180 degree HD endocavity probe with temperature-detection technol-
- Intuitive operation with M-Tuning one-key image optimization.
- Full patient database solutions: DICOM3.0, USB2.0, HDD, PDF report, etc



▶ SonoScape S6BW Portable B/W

Modes

B-mode, Steer M-mode, TDI, CFM, PDI, PWD, CWD, 3D/4D, Color M-mode, Tissue Doppler

Imaging

Linear Array, Convex array, Micro Convex Array,

Phased Array, TEE

Transducer inputs

Scan format

▶ Highlights

- Premium B/W HCU with color option
- HCU with two standard transducer sockets and convertible cart-based system design
- Micro scan real-time compound imaging

Full cardiovascular transducers:
TEE, Pediatric TEE, high frequency phased array, linear, etc.

- HD transducers (frequency range: 1.9 to 15MHz) and full patient database solutions

► SonoScape SSI-6000 Color Doppler

Technology

Modes

Micro Scan Processing, 200 degree transvaginal image with transducer temperature detection technology B, Steer M, CFM, PDI, PWD, CWD, 3D, Color M,

4D(optional)

Scan format

Trapezoid Linear, Electronic Linear, Convex-,

Micro Convex-, Phased Array, Volume Convex(optional)

Transducer inputs

▶ Highlights

-Premium 2D/3D images stem from micro scan real-time compound imaging

- A wide choice of transducers: Convex. Sector, Linear, transvaginal, etc. Maximum 15MHz frequency linear probe

- Integrated with state-of-the-art imaging technologies, like micro-scan, multiple-beam processing, etc.

200 degree transvaginal image and temperature-detection technology for endocavity transducers

▶ Sonosite MicroMaxx

Modes

2D, Tissue Harmonic Imaging, M-mode, Velocity Colour Dopppler, Color Power Doppler, PW, PW Tissue Doppler and CW Linear array, curved array, phased array, multiplane TEE and micro-convex

Transducer inputs

Scan format



Highlights

- Image quality
- Portable
- Ease-of-use
- Durable
- Reliable

► Sonosite M-Turbo

Modes

2 D, Tissue Harmonic Imaging, M-Mode, Velocity Colour Doppler, Colour Power Doppler, PW, PW Tissue Doppler and CW

Scan format

Linear array, curved array, phased array, micro convex, multidriven multiplane TEE, pencil probe, intraoperative transducer and vaginal probe

Transducer inputs 1 and 3

▶ Highlights **Highlights** – Portability – weight 3.8 kg; lithium-ion battery; full activity at the point of care

 Ease-of-use – booth up in a few seconds; clearly arranged user guidance Reliability – drop tested to withstand the real world (1.50 m)

Durability – unique 5-year warranty SonoHD™ – high resoluted image quality; 8 GB storage capability

(pictures and videos up to 60 seconds)
SonoMB™ – multibeam, real-time compound imaging

 SonoAdaptTM - tissue optimization
 SonoCalc IMT - intima media thickness measurement tool, automatic edge detection with mean and maximum thickness reporting, internal and external

Auto-gain for 2D imaging

▶ SuperSonic Imagine Aixplorer

Technology Modes

MultiWave™ Ultrasound B-mode, Color Flow, Color Power, PW Doppler, ShearWave™ Elastography, UltraFast™ Imaging, Tissue Harmonic Imaging, SuperCompound™, SuperRes™, TissueTuner™,Trapezoidal Imaging Linear array, curved array

Scan format Transducer inputs

▶ Highlights

SonicSoftware™, an all software based

architecture, providing impeccable image quality SonicTouch™ and UltraFast™ Imaging enabling

ShearWave™ Elastography to measure and display, in real-time, local tissue elasticity in kilopascal

ShearWave™ Elastography / user-skill independent, real time, quantitative, reproducible

Últra fast acquisition rate, up to 20.000 Hz

- Q-Box™ Quantitative Tools Elasticity Measurement

- Integrated BI-RADS® and new Thy-RADS™ Clinical Reporting

- Easy to operate with Interactive Touch Screen

▶ Toshiba Aplio XG

Modes Scan format 2D, 3D, 4D, M modes; PW/CW Doppler; high PRF; color flow Doppler Linear, convex, matrix, and

phased arrays; biopsy & 4D volume probes (linear and convex); Motorised-TEE; rectal, vaginal & pencil probe

Transducer inputs

3 + 1 (pencil)

▶ Highlights

- Precision Imaging, MicroPure and Elastography ApliPure Plus: Advanced realtime

compound imaging
Differential THI: better resolution

and depth of penetration Advanced Dynamic Flow:

Broadband colour flow Doppler Contrast imaging: Low MI, VRI,

microflow imaging Whole body 4D imaging with linear and convex transducers; Volume view; Multiview





▶ Highlights

Premium image quality

5 seconds bootup time

Toshiba Viamo

Modes Scan format Transducer inputs

2D, M modes; spectral Doppler; high PRF; color flow Doppler Linear, convex and phased arrays

1 + 1

- Hybrid operation with touch screen and programmable panel

Versatile mounting in desktop, cart and tablet modes

- Instant image optimisation One-click workflow control

▶ Toshiba Xario XG

Modes

Scan format

2D, 3D, 4D, M modes; PW/CW Doppler; HPRF; color flow Doppler Linear, convex and phased arrays; biopsy probe; 4D volume probes (linear and convex); Motorised-TEE; rectal and

Transducer inputs

3 + 1 (pencil)

Highlights

- Precision Imaging and MicroPure
- Aplipure Plus: advanced real-time compound imaging
- Differential THI:
- better resolution and depth of penetration
- Advanced Dynamic Flow: broadband color flow doppler
- Quick Scan: image optimisation with just one click
- Whole body 4D imaging with linear and convex transducers; Volume view; Multiview



▶ Toshiba Xario

Modes

Scan

format

2D, 3D, 4D, M modes; PW/CW Doppler; high PRF; color flow Doppler

Linear, convex, and phased arrays; biopsy probe; 4D volume probe; Motorised-TEE; rectal and vaginal probe; pencil probe

3 + 1 (pencil)

Transducer inputs



- Precision Imaging and 4D ConvexApliPure: Realtime compound
- imaging
- Advanced Dynamic Flow:
- Broadband colour flow Doppler Quick Scan: image optimization
- with just one click
- User defined programming of operating console IASSIST: Remote control via
- handheld Bluetooth controller

▶ Toshiba Nemio XG



Scan

format

2D, 3D, 4D, M modes; PW/CW Doppler; HPRF; color flow Doppler Linear, convex and phased arrays; biopsy probe; 4D volume probe; TEE; rectal and vaginal probe; pencil probe, endoscopic FNA, Laparoscopic





Highlights

- ApliPure:
- Realtime Compound Imaging
- Advanced dynamic flow: broadband color flow Doppler

 – SonoSet: Workflow control
- with just one click
- User defined programming of menus and buttons
- Onboard reporting, DICOM, DVD, USB, and export to PC

▶ Toshiba Famio XG

Modes Scan format 2D, M modes; THI Linear and convex arrays; biopsy probe; rectal and vaginal probe

Transducer inputs



Highlights

- 100% digital signal processing
- Broadband technology
- Image optimization option
- -User defined programming of the operating console
- Image and loop storage on hard-disk and CD, DICOM

Ultrasonix SonixMDP

Modes

B, Dual B, Quad B, Spatial Compound. Trapezoidal. B-Mode, Extended Sector, CFM Colour Doppler, Directional Power Doppler, Colour Steering, PW Doppler, Pulsed Doppler steering, CW Doppler, M, 3-D, Panoramic, Harmonic. Duplex, Triplex, 4D/Live 3D, Elastography, Cardiac Linear, Convex, Microconvex, Phased Array, Endocavity

Scan format Transducer inputs

▶ Highlights

- Intuitive console design, QSonix quick exam button
 Premium Image Quality, Compact Design
- OpenSONIX Platform Increased Expansion Options
- 17" LCD Monitor, DVI Input
 Dicom Compliant, Online Updates, Wireless Connectivity
- Sonix LIVE Live Streaming
- ECG Module

▶ Ultrasonix SonixOF

Modes

B. Dual B. Quad B. Spatial Compound. Trapezoidal. B-Mode, Extended Sector, CFM Colour Doppler, Directional Power Doppler, Colour Steering, PW Doppler, Pulsed Doppler steering, CW Doppler, M, Harmonic. Duplex, Triplex Linear, Convex, Microconvex, Phased

Scan format

Transducer inputs




▶ Highlights

- Intuitive console design, QSonix quick exam button

Array, Endocavity

- Premium Image Quality, Compact Design
 OpenSONIX Platform Increased Expansion Options
- 17" LCD Monitor, DVI Input
- Dicom Compliant, Online Updates



▶ Ultrasonix SonixSP B, Dual B, Quad B, Spatial

Compound. Trapezoidal. B-Mode, Extended Sector, CFM Colour Doppler, Directional Power Doppler, Colour Steering, PW Doppler, Pulsed Doppler steering, CW Doppler, M, 3D, Panoramic,

Modes

Harmonic. Duplex, Triplex, 4D/Live 3D, Elastography Scan format Linear, Convex, Microconvex, Phased Array, Endocavity

Transducer inputs

▶ Highlights

- Intuitive console design, QSonix quick exam button - Premium Image Quality, Compact Design - OpenSONIX Platform Increased Expansion Options

– 17" LCD Monitor, DVI Input

calculation software packages

-Patient data managing system

-4X windows - Pseudo color

-Dicom 3.0

-12 inch LCD

- Full range of probes

- Dicom Compliant, Online Updates, Wireless Connectivity

Ultrasonix SonixTOUCH



B, Dual B, Quad B, Spatial Compound, Trapezoidal, B-mode steering, Extended Sector, CFM Color Doppler, Power Doppler, Directional Power Doppler, Color Steering, PW Doppler, Pulsed Doppler steering, CW Doppler, M, 3-D, Panoramic, Harmonic, Duplex, Triplex, 4D/Live 3D, Elastography Linear, Convex, Microconvex, Phased Array, Endocavity

Highlights

Touch Screen Technology

- Premium Image Quality

- Application Specific Packages (Anesthesia, Critical Care, IVF, and more)


Easy to use

- Easy to Maneuver


- Compact design, small foot print, fold down 17" DVI display


- Battery Powered (option)


Online updates














> z.one *ultra*

Technology Mode

Zone SonographyTM Technology 2D / B- & M-Mode, Tissue Harmonic Doppler Imaging, Compound Harmonics, Color- and Color Power Doppler, PW- & CW-Doppler, Simultaneous Dual Imaging, 3D-Imaging, Elastography, Real-Time Triplex, Contrast Imaging Curved Array (Micro-convex), Linear Array, Phased Array, Virtual Apex Array (trapezoidal), TEE 1 – Z.ONE Scan Engine only (hand-held use) 3 – Z.ONE ULTRA or ULTRA SP (Scan Engine or Scan

Module combined with SmartCart or SmartCart SP

Scan Transducer inputs

Workstation)



- ZST Zone Speed Technology - ZSI Zone Speed Index - AUTO-OPT Automatic Optimization

IQ Scan / Retrospective Imaging (The Virtual Patient)

Utilizing the power of DSP Digital Signal

Processing chip technology Convertible / Hybrid Ultrasound Concept

Battery Pack for SmartCart Workstations

▶ Civco Ultrasound Transducer Covers

CIVCO's CIV-FlexTM Transducer Covers have been recommended for over two decades as the latex-free cover of choice and are a reliable way to provide patient and staff safety and prevent the spread of infection.



Highlights

Offers extended sterile protection when performing puncture and drainage procedures Soft, pliable, distortion-free, latex-free material

- Available telescopically-folded for easy application to transducer and in extended lengths of 36 and 58 inches

Select covers offer a three-dimensional box end

▶ Civco General Purpose Needle Guidance Systems

CIVCO's needle guidance systems offer physicians reduced technique variability, a shorter learning curve and reduced procedure time during ultrasound-guided procedures. They are designed for use

with ultrasound systems from leading equipment manufacturers



Ultra-Pro IITM Needle Guide is designed for ultrasound-guided procedures including catheter placement, core tissue biopsy, drainage procedures and fine needle aspirations

Accust TETM Needle Guide is a transverse guide designed

specifically for guidance in central line placement procedures

InfinitiTM Needle Guide offers a unique open channel with infinite angle capabilities and is ideal for deep or shallow access applications including breast biopsy and regional anesthesia

▶ Civco eTRAX™ Needle Guidance System

eTRAX allows physicians safe and precise placement of instruments in the interventional suite by tracking the tip of the needle using

electromagnetic technology in real-time navigation.



Highlights

- Enables precise percutaneous targeting of lesions without radiation or open surgery

- Ideal for ablations, core tissue biopsies, drainage, fluid aspirations, therapeutic delivery, vascular access and anesthesia

 – Accepts 18 GA or smaller instruments
- Compatible with validated ultrasound equipment and software only

▶ Siemens Acuson S2000 Automated Breast Volume Scanner



Highlights

Ideally suited to image patients with dense breast tissue and/or a history of breast disease
 Acquisition of full-field volumes of the breast automatical-


ly, quickly and comfortably


Efficient and comprehensive analysis of the volume data
 Comprehensive BI-RADS® reporting capabilities


- Patient friendly - minimal compression

No radiation


COMPANIES/SUPPLIERS Mammography Computed Tomography Speech Recognition **C**R Magnetic Resonance DR CT/MR Accessories Interventions Ultrasound Molecular Imaging R/F Display Systems Injectors RIS Printers **PACS** Workstations Agfa HealthCare Septestraat 27 B – 2640 Mortsel AGFA 🐠 HealthCare Ø +32 3 444 94 44 www.agfa.com Alliance Medical Interim Solutions Iceni Centre, Warwick Technology Alliance Park, Warwick, CV34 6DA. Medical Ø +44 (0)192 648 2000 www.allianceinterim.com www.alliancemedical.eu Barco N.V. President Kennedypark 35 8500 Kortrijk, Belgium Ø +32 56 233211 www.barco.com Canon Europa NV PO Box 2262 NL – 1180 EG Amstelveen Canon Ø +31 20 545 8831 www.canon-europe.com/medical **CHILI** GmbH Burgstrasse 61 D - 69121 Heidelberg CHILL, Ø +49 6221 18079-10 www.chili-radiology.com CIVCO Medical Solutions 102 First Street South Kalona, Iowa 52247, USA Ø +1 319-656-4447 www.civco.com info@civco.com **Covidien** Formerly Tyco Healthcare/Mallinckrodt Corso Matteotti 3 I – 20121 Milan COVIDIEN Ø +39 02 76264223 www.covidien.com CPS Color Printer Systems Vertriebs GmbH Lange Zaun 1 D – 57319 Bad Berleburg Ø +49 2751 444629









		Mammography Computed Tomography Magnetic Resonance CT/MR Accessories Interventions R/F Injectors RIS PACS			S	Speech Recognition CR DR Ultrasound Molecular Imaging Display Systems Printers Workstations							
MMS Medicor Medical Supplies GmbH Heinrich-Hertz-Strasse 6 D − 50170 Kerpen Ø +49 2273 9808-0 www.mms-medicor.de	ONI Medical Bystems, Inc.		•										
National Display Systems, Inc. Olarinluoma 9, P.O. Box 46 FIN − 02201 Espoo Ø +35 89 42 001 www.ndssi.com	SURGICAL - IMAGING												
NEC Display Solutions Europe GmbH Landshuter Allee 12-14 D − 80637 Muenchen Ø +49 89 99699-0 www.nec-displays-solutions.com	NEC												
MMS Medicor Medical Supplies GmbH Heinrich-Hertz-Strasse 6 D − 50170 Kerpen Ø +49 2273 9808-0 www.mms-medicor.de	Nemoto Nemoto				•								
OR Technology Waldemarstraße 20 g/h D – 18057 Rostock Ø +49 381 20 36 126 www.or-technology.com	O R Technology					•						•	
Philips Healthcare, Boschdijk 525, 5621 JG, NL – Eindhoven Ø +31 402785109, e.de.wilde@philips.com, www.philips.com/healthcare	PHILIPS	•				•			•			•	
Planmed Oy Asentajankatu 6 00880 Helsinki, Finland Ø +358 20 7795 300 www.planmed.com	Planmed	•											
PROTEC GmbH & Co. KG Lichtenberger Str. 35 D − 71720 Oberstenfeld Ø +49 7062 92550 www.protec-med.com	PROTE					•			•			•	








Mountain View, CA 94043-5209 Ø +1-(650) 230-2800 www.zonare.com

	Client-s		Spec	COM support ech recognition ord Party API	RIS PACS SOLUTIONS RAD-BOOK 2010		
AGFA 4960 HealthCare	• • • • •	IMPAX 6	• • • •	IMPAX 6	Agfa HealthCare ∅ +32 3 444 94 44 www.agfa.com		
CHILL	• • • • •	CHILI/Report	• • •	CHILI/Report	CHILI GmbH Ø +49 6221 18079-10, info@chili-radiology.com, www.chili-radiology.com		
FUJIFILM	• • • • •	Synapse	• • • •	Synapse	FUJIFILM EUROPE GMBH Ø +49 211 5089-246 www.fujifilm.de/medical		
GE Healthcare	• • • • •	Centricity RIS/PACS	• • • •	Centricity RIS/PACS	GE Healthcare Ø +49 7348 9861-0 www.gehealthcare.com		
GEMED Gesellschaft für medizinisches Datenmanagement mbH	• • • • •	GEMED-PACS®	• • • •	GEMED-PACS®	GEMED GmbH Ø +49-7304-9191-60 info@gemed.de, www.gemed.de		
IMAGE Information Systems	• • • • •	iQ-VIEW PRO / 3D – iQ-WEBX	• • • •	iQ-SYSTEM PACS / iQ-RIS	IMAGE Information Systems Europe Ltd. ∅ +49 381 203 38 58 / +49 381 660 80 43, Fax +49 381 203 38 59, www.image-systems.biz		
● INFINITT	• • • • •	INFINITT PACS	• • • •	INFINITT PACS	INFINITT Europe GmbH Ø 069 / 58 3000 200 www.infinitteu.com, E-mail: sales@infinitteu.com		
isoft	• • • • •	iSOFT PACS	• • • •	iSOFT Radiology / iSOFT PACS	iSOFT Health GmbH Ø +49 621 3928-0 www.isofthealth.com		
NONICA MINOLTA	• • • • •	ACIES	• •	ACIES	KONICA MINOLTA Medical Ø +31 20 659 260 www.konicaminolta.eu		
medavis 🏰	• • • • •	medavis PACS	• • • •	medavis RIS / medavis PACS	medavis GmbH ∅ +49 721 92910-360 marketing@medavis.com, www.medavis.com		
medigration	• • • • •	medigration PACS / ImageVision	• • •	medigration PACS / ImageVision	medigration GmbH ∅ +49 9131 69087-40 info@medigration.de, www.medigration.de		
MERGE Healthcare	• • • • •	FUSION PACS GL	• • • •	FUSION PACS GL	Merge Healthcare Ø +31 40 299-0771 EMEA_Sales@merge.com, www.merge.com		
0 R Technology	• • • • •	dicomPACS® Diagnostic Workstation	• • • •	dicomPACS® Diagnostic Workstation	OR Technology ∅ +49 381 20 36 126 info@or-technology.com, www.or-technology.com		
PHILIPS	• • • • •	iSite Radiology	• • • •	XRIS / iSite PACS	Philips Healthcare, ∅ +31 402785109, e.de.wilde@philips.com, www.philips.com/healthcare		
PROTEC	• • •	PROPAXX	•	CONAXX / PROPAXX	PROTEC GmbH & Co. K ∅ +49 7062 92550 protec@protec-med.com, www.protec-med.com		
SECTRA	• • • • •	IDS	• • • •	IDS	Sectra Imtec AB ∅ +46 13235200 info.imtec@sectra.se, www.sectra.se/medical		
SIEMENS	• • • • •	syngo®.via	• • • •	syngo® Workflow / syngo®.via	Siemens AG, Healthcare Sector Ø +49 9131 84-0 www.siemens.com/healthcare		
VISAGE IMAGING®	• • • • •	Visage® 7	• • • •	Visage® 7	Visage Imaging GmbH Ø +49 030 / 700968-0 www.visageimaging.com		
VIBUB	• • • • •	JiveX	• • • •	JiveX	VISUS Technology Transfer GmbH \varnothing +49 234 93693-0 sales@visus.com, www.visus.com		
VITAL	• • • • •	Vitrea® / VitalConnect®	• • •	Vitrea® / VitalConnect®	Vital Images Europe B.V. ∅ +31 704 262 181 info@vitalimages.com, www.vitalimages.com		



iSOFT Radiology. Integrated RIS-PACS solution.

iSOFT Radiology enables your organisation to engage the challenges of the future. Here an example: iSOFT Radiology in combination with iSOFT PACS integrates your whole IT environment in one optimised radiological workflow structure. An efficient referrer management and the integration of imageguided modalities offer optimised processes. Consistent interaction with existing patient management systems simplifies the planning of resources and deletes treatment delays. The result: sustained economical growth transforming your organisation into a real profit center. An easy integration by use of standard technologies strengthens your economic position and protects existing investments.

www.isofthealth.com



SKF your medical partner

A focus on quality and reliability, coupled with supporting and monitoring services, can put a device manufacturer at the forefront of its industry.

Thanks to its knowledge of the medical industry, SKF has developed a unique range of products and services. Our advanced actuation and guiding systems, bearings and seals are manufactured to meet the industry's requirements and to comply with all specifications and regulations of the medical industry, such as EN/UL and IEC 60601 approval.

SKF is integrating mechanical, electronic and information technologies to high level functions and cooperates with customers to design a new generation of smart products. Employing sensor-equipped bearings or sensorized linear actuators and mechatronic slides enable designers to substitute traditional hydraulic systems, simplifying design and minimizing environmental impact.














PoleStar® N30 Surgical MRI System


Experience Excellence

- · Intraoperative MRI and StealthStation® navigation for accuracy throughout surgery
- · Confidence in achievement of surgical goals
- New improved T2/FLAIR imaging for low grade gliomas
- · Auto-registration of images for navigation
- Easy retrofit into existing OR's
- Minimal impact on conventional surgical equipment and technique









Find a full article on the PoleStar® N30 inside this issue of RadBook